

# CR EASYLIFT

*Unique wall-mounted, floor-standing lifting device especially designed for smooth and safe positioning of the CR Full Body Cassette Holder in full body CR imaging*



- Quick and easy lifting assuring improved safety for user and patient
- Ergonomic and user-friendly design
- Increased accuracy and ease of operation for improved workflow in full body CR imaging



Cassette holder positioning made safe and easy

The CR Easylift has been developed to work with the CR Full Body Cassette Holder, allowing smooth and easy lifting at greatly reduced effort.

During installation, the CR Easylift can be modified for left or right hand operation. Its ergonomic design and user-friendly single handle adjustment system ensure fewer handling movements together with quick and easy height adjustment.

The CR Easylift is a robust system that greatly improves safety for both the user and the patient.

Smooth and ergonomic workflow

The CR Easylift brings increased accuracy and ease of operation to Full Leg Full Spine exposures, enabling precise positioning of the CR Full Body Cassette Holder with respect to the X-ray tube, regardless of patient height. This results in a smoother and more ergonomic workflow and enhanced productivity.

## Key part of a unique complete solution

The CR Easylift brings the user all the benefits of the CR Full Body Cassette Holder, the CR Full Leg Full Spine License and the CR Full Leg Full Spine Sets (plates and cassettes), with the added value of easy and safe height manipulation.

It is a key part of a complete range of integrated solutions,

guaranteeing ease of use and uniformly high quality images for all Full Leg Full Spine exposures in full body CR imaging.

# Technical

## SPECIFICATIONS

### DIMENSIONS

- Size (W x D x H): 39 x 13.5 x 230 cm (15 x 5.5 x 90.5 in)

### WEIGHT

- 73 kg (161 lbs)
- With CR Full Body Cassette Holder: approx. 100 kg (220 lbs)

### INSTALLATION REQUIREMENTS

- An Agfa technician will perform the installation of the CR Easylift, based on local guidance and instructions. Please refer to the installation manual for important installation requirements and preconditions.

### OPERATION

- During installation, the CR Easylift can be modified for left or right hand operation

### APPROVALS

- CE

### TRANSPORT DETAILS

- Packaging dimensions (W x D x H): 60 x 39 x 238 cm (24 x 15 x 94 in)
- Shipped in a wooden crate

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringement. The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications, which must be met by Agfa HealthCare. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative at [agfa.com](http://agfa.com) for availability information. Agfa HealthCare diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

Copyright 2009 Agfa HealthCare NV

All rights reserved

Printed in Belgium

Published by Agfa HealthCare NV

B-2640 Mortsel-Belgium

NFWD1 GB 200905