




770-182

Airgo Bariatric Walker

This heavy duty, lightweight walker is designed for bariatric users up to 600lbs. Featuring 1" diameter heavy gauge anodized aluminum tubing, extra-wide geometric design, and angled front legs. The walker comes with a double crossbar and a two-touch folding mechanism making it easy to open and close.

- Two Button Folding Mechanism
- Maximum weight capacity: 600 lb

AADL C007

Item number				Size	Colour	N. W. Kg
Two touch button folding walkers with regular tips						
770-150	23"	16.5"	31" - 38"	Adult	Silver	2.8
770-182	25.5"	16.5"	31" - 41"	Bariatric	Silver	3.7
Two touch button folding walkers with 5" wheels and glide tips						
770-155	23.5"	21"	31" - 38"	Adult	Silver	3.0

