



#### DON'T FORGET...

- Aseptic technique and proper skin preparation of the site are essential.
- Discard BD SafetyGlide™ in the nearest sharps collector. Follow your institution's policy for safe disposal of all medical waste.

## BD SafetyGlide™ Product Offerings

BD SafetyGlide Needles		
CATALOG #	NEEDLE	PACKAGING
<b>1 mL Syringe with Detachable BD SafetyGlide Needle</b>		
305927*	27G 5/8"	50/Box - 400/Case
305903	25G 5/8"	50/Box - 400/Case
<b>3 mL Syringe with Detachable BD SafetyGlide Needle</b>		
305904	25G 5/8"	50/Box - 400/Case
305924*	25G 1"	50/Box - 400/Case
305905	23G 1"	50/Box - 400/Case
305923*	22G 1"	50/Box - 400/Case
305906	22G 1-1/2"	50/Box - 400/Case
305909	21G 1-1/2"	50/Box - 400/Case
<b>5 mL Syringe with Detachable BD SafetyGlide Needle</b>		
305907	22G 1-1/2"	50/Box - 400/Case
<b>10 mL Syringe with Detachable BD SafetyGlide Needle</b>		
305908	22G 1-1/2"	50/Box - 400/Case
<b>BD SafetyGlide Needle</b>		
305921*	27G 5/8"	50/Box - 500/Case
305901	25G 5/8"	50/Box - 500/Case
305916	25G 1"	50/Box - 500/Case
305902	23G 1"	50/Box - 500/Case
305919*	22G 1"	50/Box - 500/Case
305900	22G 1-1/2"	50/Box - 500/Case
305915	21G 1"	50/Box - 500/Case
305917	21G 1-1/2"	50/Box - 500/Case
305918	18G 1-1/2"	50/Box - 500/Case
<b>BD SafetyGlide Needle &amp; Syringe</b>		
CATALOG #	NEEDLE	PACKAGING
<b>Insulin Syringe with Permanently Attached Needle</b>		
305930	1mL 29G 1/2"	100/Box - 400/Case
305932	1/2mL 29G 1/2"	100/Box - 400/Case
305934	1/2mL 30G 5/16"-Thinwall	100/Box - 400/Case
305935	3/10mL 29G 1/2"	100/Box - 400/Case
305937	3/10mL 31G 5/16"-Thinwall	100/Box - 400/Case
<b>1 mL Tuberculin Syringe with Permanently Attached Needle</b>		
305945	27G 1/2"	100/Box - 400/Case
305946	26G 3/8" (IB)	100/Box - 400/Case
<b>1 mL Allergist Tray with Permanently Attached Needle</b>		
305950	27G 1/2"	25/Tray - 1,000/Case
305951	26G 3/8" (IB)	25/Tray - 1,000/Case

\*Available Summer 2004.

**1.888.237.2762**  
**IN CANADA 1.800.268.5430**



BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. ©2004. BD 4/04-0420

BD Medical  
 1 Becton Drive  
 Franklin Lakes, NJ 07417  
 www.bd.com



Helping all people live healthy lives

## Safety Without Compromise

### BD SafetyGlide™ Needle

# Step-by-Step Safety

## 1 Aspiration

- Aspirate medication into syringe as per your usual technique.

## 2 Injection

- Administer injection following established technique.
- For user convenience, the needle "bevel-up" position is oriented to the lever arm. (Bevel up = lever arm up)

## 3 Activation

- After injection, immediately apply a single finger stroke to the Activation-Assist™ lever arm to activate the shielding mechanism.

# BD SafetyGlide™ Needle: where safety meets performance

## EASE OF USE

- Single-handed activation
- Minimal change in technique
- Minimal training

## SAFETY

- Secure and visible lock
- Virtually no splatter upon activation


## PATIENT CARE

- Activation occurs after needle is withdrawn: no negative impact to the patient

## EFFICIENCY

- Detachable needles fit on any size conventional syringe
- Facilitates standardization and SKU reduction

### PRIOR TO ACTIVATION



### AFTER ACTIVATION



Use BD SafetyGlide™ Syringes  
for all your insulin, TB and allergy needs.

## Setting the standard in safety

More than a decade ago, BD pioneered safety injection with the introduction of BD Safety-Lok™, the world's first safety-engineered syringe. BD has since set a new standard in safety with the BD SafetyGlide™ Needle, a single-handed, needle-based safety device and the new smaller-dose BD SafetyGlide™ Syringe.



BD, the industry leader in safety injection, offers the most extensive array of safety-engineered injection products available. BD's record of innovation spans more than a century and our continued investment in safety technology shows our commitment to providing innovative products and solutions for the future.