

ErgoFit® Bariatric sling

The ErgoFit® Bariatric sling is a premium design sling that provides reinforced material and stitching in the most stressed areas of the sling including the shoulder, legs and spine.

The extra padding provides more comfort and less pressure on skin. Re-inforcement in leg pieces prevent leg roping, enhancing comfort for the patient.

# Product	Description	Dimensions	
		cm	in
MAA8000-L	Bariatric Sling	A	138 54
	Capacity	B	105 41
	1000 lb / 455 kg / 71.43 st	C	160 63


General Characteristics

- Guarantee against manufacturing defects: 1 year
- Maximum loading capacity: 1000 lb / 454 kg / 71.43 st
- Positioning handles at the back to help with transfers, and adequately position the patient above a chair, to simply and quickly transfers from a sitting to an inclined position and to help pivoting the user.
- Machine washable; easy to care for and clean.
- Double padding in the thigh section with ultra flexible nylon/polyester, which prevents skin irritations. Easily installed on a person already sitting or lying down.
- Strong nylon & polyester straps allowing the person to be in a sitting, semi-inclined or inclined position.
- User positioning system based on the color of the straps.
- Straps especially designed for the legs, prevents falls and allows the adjustment and positioning of the legs.

Technical Specifications

Conform to ISO 10535 & CE marking

*Conditionally to a normal use of an appropriate BHM Medical patient lift and maintenance in accordance with the labels.



INSTALLATION OPTIONS

Bridge, leg closed with straps interlaced



Independent, legs open



Cradled, legs closed

