




DESIGNED FOR DIGITAL SKIN MARKERS FOR MAMMOGRAPHY

Highlight areas of concern without distraction in Digital Mammography with Beekley's exclusive Designed for Digital® SPOTS® Skin Markers.

Exclusive lower density markers allow complete visibility of tissue detail through the marker without distraction or artifact. Available in Beekley's most popular styles of skin markers for communicating **nipples**, **moles**, **scars**, **palpable masses** and **areas of concern or pain** on mammographic images.

*Designed for Digital® SPOTS® Skin Markers are a product based on the needs of the imaging center (preference product); alternatively, try Beekley's **Light Image™ SPOTS®**.*



2.3mm non-metallic pellet




N-SPOT® Designed for Digital® Soft 'n Stretchy®

Florals™ order code: 790

Contents per box: 94

- 2.3mm non-metallic pellet
- pinch-free™ material with adhesive-free center
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



2.3mm non-metallic pellet




NEW! N-SPOT® Designed for Digital® Soft 'n Stretchy® plus

Florals™ order code: 790plus

Contents per box: 94

- 2.3mm non-metallic pellet
- larger adhesive-free area for increased patient care; sticks securely on the breast under compression
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



2.3mm non-metallic pellet




N-FREE SPOT® Designed for Digital®

Florals™ order code: 732

Contents per box: 110

- 2.3mm non-metallic pellet
- extended side tabs and adhesive-free center
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



2.3mm non-metallic pellet



V-SPOT® Designed for Digital®

Florals™ order code: 703

Contents per box: 130

- 2.3mm non-metallic pellet
- easy pick-up™ tab
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



see-through mole marker



O-SPOT® Designed for Digital®

Florals™ order code: 752

Contents per box: 58

- uniform ultra low-density see-through ring
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



low-density see-through ring




NEW! O-SPOT® Designed for Digital® Soft 'n Stretchy®

Florals™ order code: 791

Contents per box: 55

- 1mm uniform low-density see-through ring
- pinch-free™ material
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



1.0mm diameter see-through scar marker




S-SPOT® Designed for Digital®

Florals™ order code: 750
Contents per box: 327cm cut to measure roll

- 1mm diameter uniform ultra low-density see-through line
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



1.0mm diameter see-through scar marker




S-SPOT® Designed for Digital® Soft 'n Stretchy®

Florals™ order code: 777
Contents per box: 280cm cut to measure roll

- 1mm diameter uniform ultra-sheer see-through line
- pinch-free™ material
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



see-through palpable mass marker



A-SPOT® Designed for Digital®

Florals™ order code: 751
Contents per box: 120

- uniform ultra low-density see-through triangle
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)



see-through square



Altus™ Designed for Digital®

Florals™ order code: 754
Contents per box: 60

- uniform ultra sheer low-density see-through square
- designed for use with the most sensitive digital equipment

[Request sample](#) | [Pricing](#)