

## General Imaging Skin and Cassette Markers

### General Imaging Products:

Beekley SPOTS for use in General Imaging serve a variety of purposes within radiology and the emergency room.

Products within Beekley's General Imaging category can be used to communicate points of pain on film, document right or left on x-rays and identify entry and exit sites.

All Beekley SPOTS® skin markers are lead-free, latex-free, and hypoallergenic.



#### X-SPOTS®

order code: 101

contents per box: 150

- 1.5mm pellet with easy pik-up tab
- perfect for general imaging



#### Y-SPOTS®

order code: 102

contents per box: 132

- 2.0mm pellet with easy pik-up tab
- perfect for general imaging



#### D-SPOTS® arrow

order code: 104

contents per box: 99

- radiopaque arrow with easy pik-up tab
- perfect for injury sites