

# Which ultrasound table is right for you?



## Sound Pro Combination Table

### SUITABLE FOR:

- Abdomen • Renal • Aorta • Pelvis
- Lower Extremity Arterial • Endovaginal • OB
- Hysterosonography • Thyroid • Testicular
- Breast • Musculoskeletal • Carotid
- DVT/Lower Extremity Venous • Venous Reflux
- Upper Extremity Venous
- Upper Extremity Arterial • Renal Artery
- Renal Transplant • Abdominal Doppler
- SMA • Echocardiography • Stress Echo
- Transesophageal Echo

### SPECIFICATIONS:

**Dimensions:** 71" l x 30" w (180.3 x 76.2 cm);

35" w (88.9 cm) with side rails

**Drop-Down Leg Section:** 12.6" l (32 cm)

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Tabletop:** Three primary sections: torso section features sonographer and patient cutouts with Fowler back positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access

**Cardiac Scanning Cutout:** with pivoting cushion that doubles as an armrest

**Adjustable Sonographer Cutout:** drops down or folds up for maximum access to the patient

**Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®

**Upholstery Color:** Graphite

### Motions:

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Trendelenburg:** 0° to ±15°

**Fowler Back:** 0° to 80° infinitely adjustable

### Controls:

**Hand Control:** Activates height, Trendelenburg motions, Fowler back and auto level motions

**Foot Control:** Activates height, Trendelenburg motions and Fowler back

**Wheels:** 5" (12.7 cm) Individual locking swivel casters; central floor locking system

**Patient Restraints:** One body strap

**Finish:** Powder coat

**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.

**Weight:** 350 lb (159 kg)

**Power:** 115 VAC or 230 VAC

**Warranty:** Two year

**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



**058-710** Table, Combination, Sound Pro, 115 VAC

*Includes: hand and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.*

**058-715** Table, Combination, Sound Pro, 230 VAC

*Includes: hand and foot controllers, folding side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.*

### Replacement:

**056-612** Table Paper, Hygienic, 12 rolls  
18" w x 225' l (46 cm x 69 m)



## Echo Pro Echocardiography Table

### SUITABLE FOR:

- Abdomen • Renal • Aorta • Pelvis • OB
- Lower Extremity Arterial • Thyroid • Testicular
- Breast • Musculoskeletal • Carotid
- DVT/Lower Extremity Venous • Venous Reflux
- Upper Extremity Venous
- Upper Extremity Arterial • Renal Artery
- Renal Transplant • Abdominal Doppler
- SMA • Echocardiography • Stress Echo
- Transesophageal Echo

### SPECIFICATIONS:

**Dimensions:** 73" l x 30" w (185.4 x 76.2 cm);

35" w (88.9 cm) with optional side rails

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Tabletop:** Two primary sections: torso section: features sonographer and patient cutouts with Fowler back positioning, and lower body section remains fixed

**Cardiac Scanning Cutout:** with pivoting cushion that doubles as an armrest

**Adjustable Sonographer Cutout:** drops down or folds up for maximum access to the patient

**Mattress:** Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®

**Upholstery Color:** Graphite

### Motions:

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Trendelenburg:** 0° to ±15°

**Fowler Back:** 0° to 80° infinitely adjustable

### Controls:

**Hand Control:** Activates height, Trendelenburg motions, Fowler back and auto level motions

**Foot Control (optional):** Activates height, Trendelenburg motions, and Fowler back

**Wheels:** 5" (12.7 cm) Individual locking swivel casters; central floor locking system

**Patient Restraints:** Two body straps

**Finish:** Powder coat

**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.

**Weight:** 335 lb (152 kg)

**Power:** 115 VAC or 230 VAC

**Warranty:** Two year

**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



**058-700** Table, Echocardiography, Echo Pro, 115 VAC  
*Includes Hand Controller.*

*All options available separately\**

**058-705** Table, Echocardiography, Echo Pro, 230 VAC  
*Includes Hand Controller.*

*All options available separately\**

### Optional:

**058-740** Foot Controller

**058-632** Side Rails, Folding

**058-735** Arm Board, Vascular Scanning

**058-656** IV Pole

**058-610** Dispenser, Paper  
*with Cutter Strap*

**058-640** Foot Support

### Replacement:

**056-612** Table Paper, Hygienic, 12 rolls  
18" w x 225' l (46 cm x 69 m)

*\*Note: Options are installed by the customer.*



## Ultra Pro Ultrasound Table

### SUITABLE FOR:

- Abdomen • Renal • Aorta • Pelvis
- Lower Extremity Arterial • Endovaginal • OB
- Hysterosonography • Thyroid • Testicular
- Breast • Musculoskeletal • Carotid
- DVT/Lower Extremity Venous • Venous Reflux
- Upper Extremity Venous
- Upper Extremity Arterial • Renal Artery
- Renal Transplant • Abdominal Doppler
- SMA • Transesophageal Echo

### SPECIFICATIONS:

**Dimensions:** 70" l x 30" w (185.4 x 76.2 cm);

35" w (88.9 cm) with side rails;

**Drop-Down Leg Section:** 12.6" l (32 cm)

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Tabletop:** Three primary sections: torso section with Fowler back positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access

**Mattress:** Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®

**Upholstery Color:** Graphite

### Motions:

**Height Adjustable:** 23" to 39" (58.4 to 99 cm)

**Trendelenburg:** 0° to ±15°

**Fowler Back:** 0° to 80° infinitely adjustable

### Controls:

**Hand Control:** Activates height, Trendelenburg motions, Fowler back and auto level motions

**Foot Control (optional):** Activates height, Trendelenburg motions, and Fowler back

**Wheels:** 5" (12.7 cm) Individual locking swivel casters; central floor locking system

**Patient Restraints:** Two body straps

**Finish:** Powder coat

**Patient Capacity:** 500 lb (227 kg); weight tested to four times the patient load rating.

**Weight:** 335 lb (152 kg)

**Power:** 115 VAC or 230 VAC

**Warranty:** Two year

**Certifications:** ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



**058-720** Table, Ultrasound, Ultra Pro, 115 VAC  
*Includes Hand Controller.*

*All options available separately\**

**058-725** Table, Ultrasound, Ultra Pro, 230 VAC  
*Includes Hand Controller.*

*All options available separately\**

### Optional:

**058-740** Foot Controller

**058-632** Side Rails, Folding

**058-651** Stirrups, Retractable

**058-735** Arm Board, Vascular Scanning

**058-656** IV Pole

**058-610** Dispenser, Paper  
*with Cutter Strap*

### Replacement:

**056-612** Table Paper, Hygienic, 12 rolls  
18" w x 225' l (46 cm x 69 m)

*\*Note: Options are installed by the customer.*