

# BARIATRIC SERIES



**Bariatric Full Cabinet Treatment Table**

	<i>Length</i>	<i>Height</i>	<i>Width</i>
<b>6160</b>	78"	31"	33"

- Plywood top secured in place with tee nut fastening system
- Steel reinforced stress points
- Cabinet constructed with steel fasteners
- Furniture-grade laminated plywood cabinet frame
- 3 storage areas, each with an adjustable shelf
- 6 doors with concealed, self-closing Euro-style hinges
- 600 lbs. capacity under normal use



**Bariatric Treatment Table**

	<i>Length</i>	<i>Height</i>	<i>Width</i>
<b>6190</b>	78"	31"	30"

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- 600 lbs. capacity under normal use



**Bariatric Treatment Table/Low Height**

	<i>Length</i>	<i>Height</i>	<i>Width</i>
<b>6196</b>	78"	26"	30"

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- Lower height makes patient access easier
- 600 lbs. capacity under normal use



**Bariatric Treatment Table**

	<i>Length</i>	<i>Height</i>	<i>Width</i>
<b>6198</b>	72"	31"	30"

- Plywood top and welded steel frame
- 6-leg design with bolted-through assembly
- Leg levelers
- Laminate full length shelf with H-braces
- 600 lbs. capacity under normal use



### Extra Wide Bariatric Straight Top Power Table

<b>84100-34</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	72"	34"	18"-35"

<b>84100-40</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	72"	40"	18"-35"

- All welded uni-frame constructed steel base and top frame
- 1" plywood top for extra support
- Hands-free foot control
- Optional locking casters (*see other options, page 56*)
- 550 lbs. load capacity under normal use



### Extra Wide Bariatric Power Table with Adjustable Backrest

<b>84200-34</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	72"	34"	18"-35"
Backrest	38"	34"	

<b>84200-40</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	72"	40"	18"-35"
Backrest	38"	40"	

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded uni-frame constructed steel base and top frame
- 1" plywood top for extra support
- Hands-free foot control
- Optional locking casters (*see other options, page 56*)
- 550 lbs. load capacity under normal use



### Extra Wide Bariatric Power Table with Adjustable Backrest & Drop Section

<b>84430-34</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	76"	34"	18"-35"
Backrest	38"	34"	
Drop Section	13"	34"	

<b>84430-40</b>	<i>Length</i>	<i>Width</i>	<i>Height Range</i>
Overall	76"	40"	18"-35"
Backrest	38"	40"	
Drop Section	13"	40"	

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- All welded uni-frame constructed steel base and top frame
- Hands-free foot control
- 1" plywood top for extra support
- Optional locking casters (*see other options, page 56*)
- 550 lbs. load capacity under normal use

# COLOR SELECTOR

## Upholstery Colors



3SB-Slate Blue



3DT-Desert Tan



3BK-Black



3PG-Purplegray



3VP-Viscaya Palm



3AS-Allspice



3SP-Sapphire




3GP-Graphite



3AB-Aubergine



3AL-Alabaster



3NL-Neutral

## Base Laminate Colors



1NT.\* Natural



1DC.\* Dark Cherry




1MP-Maple

## CLINTON VINYL UPHOLSTERY

Clinton uses only premium, woven knit-backed upholstery, that resists stretching and bunching associated with ordinary vinyl upholstery. Plus, Clinton upholstery is topcoated to resist stains, scuffs and abrasions. *Standard flammability specifications meet CAL 117. For additional flammability requirements, call for quote.*

### ORDERING INFORMATION


1. **PICK A TOP COLOR**  
*Upholstery for table or laminate for cabinet*
2. **CHOOSE A LAMINATE FRAME COLOR**  
*Legs, frame and table backs*
3. **SELECT DOOR & DRAWER LAMINATE COLOR**  
**NOTE:** *If one color is selected, it will be used on the entire base*
4. **SELECT DOOR & DRAWER PULLS**  
*(See page 22 for complete selection)*  
**NOTE:** *If no selection is noted on the order, Gray pulls will be used on solid colors, and Silver Wire pulls will be used on wood tones.*



## Color Coordinating Upholstery



3GM-Gunmetal



3DG-Dove Gray



3CM-Cream



3CS-Clamshell



3WW-Wedgewood



3RB-Royal Blue



3EM-Emerald




3MB-Mulberry




3BG-Burgundy


## Color Coordinating Laminates



2SG-Slate Gray



2GR-Gray



2FS-Fossil



2SD-Sand



2AB-Aztec Blue




2DB-Deep Blue



2IG-Imperial Green



2WR-Wild Rose



2SGA-Sangria

*Due to the limitations of the printing process, actual colors may vary. Actual upholstery and laminate samples are available.*

### Quick Ship Log—3 to 5 days

*Quick ship items do have similar item quantity restrictions and also have some option limitations. Call for more information.*