

## Custom Comfort Child Booster Chair for Phlebotomy Chairs

### Child Booster Seat DB4998



- Added Items:
- Child Booster Seat
  - Elevates the Child an Additional 6"
  - Padded Seat and Back
  - Fits Standard Width 1201 Series Chairs
  - 100 lbs. Patient Weight Capacity

Our DB4998 Child Booster Seat is easily placed on any of our standard width phlebotomy chairs with no tools required. The 10" x 12" upholstered seat and back provide an additional 6" of extra seat height. A safety belt aids in securing the pediatric patient into the seat.