

Custom Comfort Inc XL Bariatric Blood Chair

**Blood Draw Chair
1201-SXXL**



- Added Items:
- Fully Upholstered
 - 20" Standard Seat Height
 - Extra Wide 35" Patient Width
 - 1,000 lbs. Patient Weight Capacity

Our Bariatric series chairs are built to safely accommodate your larger patients needs or parents with children. The 1201-LXL series and the 1201-SXXL series have 10" and 15" of additional seating area respectively. The chairs feature a fully upholstered seat and back for comfort and are upholstered in a commercial grade vinyl fabric. The fully adjustable arms on the 1201-LXL series have two flip-up arms for patient safety and easier vein puncture. The 1201-SXXL series feature heavy duty straight arms to help support the largest patients.