

## Custom Comfort Infant Phlebotomy Station

### Infant Station MA3091-FS



- Added Items:
- 36" Working Height (Without Casters)
    - Arm Extensions
    - 18" x 33" Bed Area
  - 150 lbs. Patient Weight Capacity
  - Shown With Optional Locking Casters

Our MA3091 series infant stations have all the same features as the MA3014 with the addition of arm extensions. The arm extensions allow the infant's arm to be supported and secured in a flat position. Adjustable Velcro straps provide additional safety and accuracy for the technician while comfortably restraining the arm of the infant during the phlebotomy procedure.