

Custom Comfort Tall Height Blood Drawing Chair

Tall Height Blood Draw Chair 1201-LU/XT



Added Items:

- Fully Upholstered
- 26" ERGO Seat Height
- 8" Adjustable Arm Height
- Extra Wide Padded "L" Arms
- 325 lbs. Patient Weight Capacity

Our blood draw chairs are ideal for clinics, labs, outpatient centers, hospitals or anywhere you need to draw blood. Optimum vein puncture level and patient comfort can be achieved with the infinite adjustable armrests. The armrests also have the capability of rotating left or right to help accommodate any size patient. The flip-up padded arm helps protect your patient in the event they slump forward or faint. The fully welded steel frame is epoxy coated and has a 5 year warranty. The ERGO height chairs, like the 1201-LUXT, feature a chair seat and frame that is 5" taller than most standard phlebotomy chairs. Patients, especially those with limited mobility, will find the higher seat enables easier entry and exit.