




Dragon X

SPS-HF4.0



All Purpose

Portable X-ray Unit

The Dragon X is designed to carry the X-ray unit from door to door of patient room in clinics and hospitals and therefore, can be reduced patient positioning time and efforts.

Standard applications are: medical radiographic use in nursing homes, private homes, field clinics, correctional facilities, military field operations, intensive care units, emergency units, extremity orthopaedic or anywhere an x-ray machine must be brought to a patient.

- Small in size, big in performance (4 kW).
- High frequency inverter (300 kHz).
- Automatic line voltage compensation (90-264 VAC).
- Automatic output power adaptation to power line conditions (line impedance).
- Shorter exposure times than conventional X-ray units.
- Dual focal spot provides superior image quality.
- Easy to transport - friendly user.

