

# Surgical Accessory Products

## SURGICAL CLIPS AND CLAMPS

### Multi-Purpose Tube Holder


DeRoyal is pleased to offer a Multi-Purposed Tube Holder that frees clinicians from the task of securing lines with clumsy, make-shift tape constructions that are hard to secure and harder to keep in place and remove.

- Designed for quick and easy placement of all sizes of tubing and lines
- Easy-to-apply adhesive base attaches firmly to skin, cloth or drapes
- Gentle to skin and easily removed without painful pulling
- Holds lines securely in place once line is positioned on patient
- Secure placement for positioning to create slack without bladder traction, thus decreasing risk of urinary track infection
- Latex safe

#### For use with:

- Jejunostomy Tubing
- Leg Bag Tubing
- Subclavian Port Lines
- Wound drainage tubing
- Peritoneal dialysis exit site tubing
- A variety of other tubes and catheters

Product No.	Description	Quantity
30-1002	Tube Holder, Medium, Non-Sterile	50 Pks/Cs
30-1003	Tube Holder, Large, Non-Sterile	50 Pks/Cs




#### Competitors' Products

Other tube holders allow the tubing line to slide and dislodge from place.

#### DeRoyal Tube Holder

DeRoyal Tube Holder's adhesive surface holds tubing securely in place without slipping or displacement.



Tab folds to Velcro® center for use with or without adhesive center.

Adhesive backing on central tab prevents tube dislodgement.



DeRoyal Multi-Purpose Tube Holder can be attached to skin, drapes or cloth for use during a variety of procedures to accommodate numerous styles of lines.



## SURGICAL CLIPS AND CLAMPS

### Towel Clamp

- Positions drapes, towels and sponges in sterile field
- Built-in wire holder for cautery pencil or other wires
- Special bridge in curved jaws clamp and occlude a variety of tubes
- Strong and durable, the nylon reinforced plastic will not break
- Latex safe

Product No.	Description	Quantity
30-901	Towel Clamp	1/Pk, 80 Pks/Cs
30-902	Towel Clamp	2/Pk, 120 Pks/Cs
30-900NS	Towel Clamp, Non-Sterile	100/Cs
30-9010NS	Towel Clamp, Stainless Ball & Cup, 5 1/4", Non-Sterile	12/Cs



### Tube Positioner

- Allows stationary placement of tubes, lines, and cords in sterile field
- Extended adhesive peel-off backing for anchored placement
- Easily removed by gloved personnel
- Raypaque element ensures x-ray detectability
- Latex safe

Product No.	Description	Quantity
30-802	Tube/Cord Holder	25 Pks/Bx, 2 Bxs/Cs




# Surgical Accessory Products

## SURGICAL CLIPS AND CLAMPS

### Breathing Circuit Support Clamp

- Holds breathing circuits and suction tubes in place
- Relieves the "tug" on the ET (endotracheal tube) connection
- Reduces risk of kinks in ET tube, more comfortable for ICU patient
- Reduces chance of a disconnection
- Provides suction handle support
- Suction handle is more easily accessible to staff and patient
- Easy to use and holds tubing securely to patient's gown, bedding or drape
- Convenient support and easily positioned by staff
- Latex safe

Product No.	Description	Quantity
JLC-3NS	Breathing Circuit Support Clamp, Non-Sterile	100/Cs



### Anesthesia Drape "Peach" Clip

- Used to attach surgical drape to IV poles to create barrier at patient's head
- Can be used for a variety of other clamping purposes
- Heavy duty holding power that will not slip
- Wide, easy-to-grip surface
- Reusable, cost-effective, washable
- Non-rusting components
- Multiple colors
- Latex safe

Product No.	Description	Quantity
JPC-3NS	Assorted Colors, Non-Sterile	8/Bx



## SURGICAL CLIPS AND CLAMPS

### Tubing & Wire Holder

- Securely positions and holds suction tubing, cautery pencil wires and bipolar wires
- Eliminates need for expensive drapes with ineffective built-in tube holders
- Simple, one-handed operation
- Helps prevent instruments from falling off sterile field
- Latex safe

Product No.	Description	Quantity
MSC-1	Suction Tubing & Cautery Wire Holder	10/Bx, 2 Bxs/Cs
MSC-1NS	Suction Tubing & Cautery Wire Holder, Non-Sterile	100/Bx

