



11130-2

BARIATRIC  
**650 LB**  
WEIGHT CAPACITY



**ENVIRONMENT FRIENDLY PRODUCT**  
Recycled materials were used  
to make this product.

### Deluxe Bariatric Commode

11130-2	Assembled, 2/cs
11129	Replacement Wing Nut, 2/bx
11130B	Replacement Back Rest, 1/ea
11130S	Replacement Seat, 1/ea
11135TIP	Replacement Tip, 1" Grey, 1/ea

#### FEATURES:

- Accommodates individuals up to 650 lbs.
- Back can be removed without tools.
- Large, padded backrest provides patient comfort.
- Durable, heavy duty, grey powder coated steel tubing is sturdy and easy to maintain.
- Large, durable "snap-on" seat.
- Plastic armrests provide comfort and stability.
- Comes complete with 12 qt. commode bucket and cover that removes from behind.
- **Limited Lifetime Warranty.**



11135-1

BARIATRIC  
**1000 LB**  
WEIGHT CAPACITY



**ENVIRONMENT FRIENDLY PRODUCT**  
Recycled materials were used  
to make this product.

### Deluxe Bariatric Drop-Arm Commode

11135-1	Assembled, 1/cs
11135-S	Replacement Seat, 1/ea
11135TIP	Replacement Tip, 1" Grey, 1/ea
11135B	Replacement Back Support, 1/cs

#### FEATURES:

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Accommodates individuals up to 1000 lbs.
- Durable, heavy duty, powder coated steel tubing is sturdy and easy to maintain.
- Large, durable "snap-on" seat.
- Padded armrests and backrest.
- Pail easily slides in from the front of the commode.
- Comes complete with 12 qt. commode bucket and cover that removes from behind.
- Large backrest provides patient comfort.
- **Limited Lifetime Warranty.**

#### SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11130-2	23.5"	26"	23"	18.5"	22"	15"	20.5 lbs.	46 lbs.	650 lbs.
11135-1	26"	28"	23"	18.5"	23"	19"	22 lbs.	26 lbs.	1000 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.



11138-1

## Large Bariatric Drop Arm Commode

11138-1	Assembled, 1/cs
11138B	Replacement Back Rest, 1/ea
11138SEAT	Replacement Seat, 1/ea
10374BK-8	Replacement Tip, 1" Black, (Has replaced 11100B), 8bx/cs

### FEATURES:

- 850 lb. Weight Capacity
- Extra Large opening in toilet seat, 8 1/2" x 12".
- Durable heavy duty, powder coated steel construction tubing is sturdy and easy to maintain.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Extra large, padded, deep seat provides comfort.
- Padded armrests and backrest provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- **Limited Lifetime Warranty.**



11132-1

## Bariatric Extra Large Drop Arm Commode

11132-1	Assembled, 1/cs
11132SEAT	Replacement Seat, 1/ea
10374BK-8	Replacement Tip, 1" Black, (Has replaced 11100B), 8bx/cs

### FEATURES:

- 1000 lb. Weight Capacity
- Extra large opening in toilet seat, 8 1/2" x 12".
- Heavy duty, white powder coated steel construction.
- Extra large, padded deep seat provides comfort.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Padded armrests provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- **Limited Lifetime Warranty.**

### SPECIFICATIONS

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11138-1	25.5"	27.5"	24"	21"	18.5"	18.5"	37 lbs.	41 lbs.	850 lbs.
11132-1	36.5"	40"	36"	21"	18.5"	NA	49 lbs.	53 lbs.	1000 lbs.