

REFER TO PAGE 182-190  
FOR LIFT ACCESSORIES  
AND SLINGS.

RECOMMENDED  
FOR USE WITH 2 &  
4 POINT SLINGS!

## Battery Powered Patient Lift with Six Point Cradle

13240 1/cs

BARIATRIC  
**450 LB**  
WEIGHT CAPACITY

### FEATURES:

- Fast, easy charge from AC outlet.
- **To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.**
- Heavy gauge construction.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Audible warning when battery power falls below 50%.
- The "CHARGE" light will automatically shut off when battery is fully charged.
- **Ideal for Homecare use.**
- For use with 4 point and 2 point slings.
- Wide base provides additional stability.
- Lifting capacity: 450 lbs.
- Limited Lifetime Warranty.



NEW & IMPROVED  
CRADLE

13240  
with Sling  
(Sold Separately)

CAN FIT UNDER  
A LOW BED!

## Battery Powered Patient Lift with Rechargeable, Removable Battery

13242 1/cs

BARIATRIC  
**450 LB**  
WEIGHT CAPACITY

### FEATURES:

- Same features as 13240.
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

*To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.*

### REPLACEMENT PARTS:



**External Charger  
with Power Cord**

13255 1/ea



**Wall Mounted  
Bracket**

13260 1/ea



**Rechargeable  
Battery**

13258 1/ea



NEW & IMPROVED  
CRADLE

13242  
with Sling (Sold Separately)

CAN FIT UNDER A LOW BED!

### SPECIFICATIONS

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13240	76"	27.5"	41"	15"	4.5"	2"	5" Rear	106 lbs.	116 lbs.	450 lbs.
13242	76"	27.5"	41"	15"	4.5"	2"	5" Rear	106 lbs.	116 lbs.	450 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.



**NEW & IMPROVED CRADLE**

13244  
with Sling  
(Sold Separately)

**CAN FIT UNDER A LOW BED!**

**BARIATRIC**  
**600 LB**  
WEIGHT CAPACITY

**REFER TO PAGE 182-190 FOR LIFT ACCESSORIES AND SLINGS.**

**Bariatric Battery Powered Lift with Six Point Cradle**

13244 1/cs

**FEATURES:**

- 600 lb. weight capacity.
- To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Base legs adjust easily and lock securely into open position with padded shifter handle.
- **Ideal for homecare use.**
- Extra wide coated cradle with 360° rotation with 6 point hook up.
- High lift range allows transfer to high surfaces.
- Dual, 3" heavy duty front casters.
- Warranty: Lifetime on frame, 1 year on components and electrical parts.

**NEW & IMPROVED CRADLE**

**BARIATRIC**  
**600 LB**  
WEIGHT CAPACITY



**CAN FIT UNDER A LOW BED!**

13245  
with Sling  
(Sold separately)

**Bariatric Battery Powered Patient Lift with Rechargeable, Removable Battery**

13245 1/cs

**FEATURES:**

- Same features as 13244
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

*To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.*

**REPLACEMENT PARTS:**



**External Charger with Power Cord**

13255 1/ea



**Wall Mounted Bracket**

13260 1/ea



**Rechargeable Battery**

13258 1/ea

**SPECIFICATIONS**

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13244	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.
13245	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.