

BARIATRIC SENTRA EC HEAVY DUTY



**NOW
AVAILABLE
IN UP TO 24",
AND 450 LB.
WEIGHT
CAPACITY.**

**BARIATRIC
450 LB
WEIGHT CAPACITY**

**NOW
COMES WITH
RIGGINGS IN
CARTON!**

Bariatric Sentra EC Heavy Duty

Deluxe Dual Axle (K0007)

FEATURES:

- **Standard dual axle provides easy transition of seat to hemi-level.**
- **"Comes with Dual Cross Bar".**
- Carbon Steel frame with triple coated chrome for an attractive, chip-proof, maintainable finish.
- Double embossed vinyl upholstery is durable, light-weight, attractive and easy-to-clean.
- Urethane tires mounted on composite wheel provides durability and low maintenance.
- Aluminum footplates.
- **Comes with padded arm rests.**
- No maintenance, sealed, precision bearings.
- **8" front casters, adjustable in 3 positions.**
- Chrome handrims standard.
- Comes with push to lock wheel locks.
- Color: black upholstery
- Limited 3 Year Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*

BARIATRIC SENTRA HEAVY DUTY, EXTRA, EXTRA WIDE, (DELUXE)



AVAILABLE
IN 26", 28"
AND 30"
SEAT
WIDTHS.

SHOWN WITH OPTIONAL
HDELR SWING-AWAY
ELEVATING LEGRESTS.

Deluxe Sentra Heavy Duty, Extra, Extra Wide Dual Axle (K0007)

BARIATRIC
700 LB
WEIGHT CAPACITY

FEATURES:

- Steel magnesium wheels are strong and maintenance free.
- 700 lb. Weight Capacity.
- Powder coated silver vein frame.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- **Durable triple reinforced nylon upholstery.**
- Extra large, contoured plastic armrests provide comfort and support.
- Precision sealed magnesium wheel bearings in front and rear ensure long lasting performance and reliability.
- Extra large heavy duty, aluminum tapered footplates provide maximum support, are attractive and crack proof. (Must be ordered separately)
- Footrests are tapered for comfort.
- **Heavy duty, flat free 8" front casters are adjustable in 2 positions.**
- Extra heavy duty front forks.
- Dual axle provides easy transition of seat to hemi-level.
- Reinforced side frame, caster journals and dual cross bars.
- Comes with push to lock wheel locks.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

BARIATRIC SENTRA EXTRA HEAVY DUTY

COMES WITH
DUAL CROSS
BAR!

AVAILABLE
IN 20", 22"
AND 24"
SEAT
WIDTHS.

NOW
COMES WITH
RIGGINGS IN
CARTON!



SHOWN WITH OPTIONAL STDELR-TF
SWING-AWAY ELEVATING LEGRESTS.

Sentra Heavy-Duty

Dual Axle when used with STDS1074 Caster
(Optional Caster Sold Separately) (K0007)

FEATURES:

- Carbon steel frame with TRIPLE coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable, heavy gauge, reinforced naugahyde upholstery.
- Padded armrests provide added patient comfort.
- **Composite Mag-style wheels with chrome handrim are lightweight and maintenance free.**
- Dual position front fork.
- Dual position front frame.
- Comes standard with carry pocket on back rest.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Aluminum footplates are attractive, crack-proof and lightweight.
- **"Hemi" adjustable caster is optional (Part # STDS1074).**
- **Reinforced side frame, caster journals and dual cross bars.**
- Comes with push to lock wheel locks.
- Black upholstery, chrome frame.
- Limited Lifetime Warranty.

BARIATRIC

500 LB

WEIGHT CAPACITY

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

BARIATRIC SENTRA EC HEAVY DUTY, EXTRA, EXTRA WIDE, (STANDARD)



AVAILABLE
IN 26", 28"
AND 30"
SEAT
WIDTHS.

SHOWN WITH OPTIONAL
PH-SF SWING-AWAY
FOOTRESTS.

Sentra EC Heavy Duty, Extra, Extra Wide (K0007)

BARIATRIC

700 LB

WEIGHT CAPACITY

FEATURES:

- Steel, spoke wheels are strong and maintenance free.
- 700 lb. Weight Capacity.
- Carbon steel frame with TRIPLE coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable reinforced nylon upholstery.
- Extra large, padded armrests provide comfort and support.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Extra large heavy duty, aluminum footplates provide maximum support, are attractive and crack proof. (Must be ordered separately)
- Heavy duty, flat free 8" front casters are adjustable in 3 positions..
- Extra heavy duty front forks.
- Reinforced side frame, caster journals and dual cross bars.
- Comes with push to lock wheel locks.
- Color: Black upholstery, chrome frame.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.