

**522KL  
Digital Pediatric Scale**

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb 0.5 oz  
0 - 9 kg / 5g; 9 - 23 kg / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG conversion, LB /KG lock out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Tray Size: 24 1/16" (w) x 14 1/8" (d) x 3" (h)
- Measuring Tape: 0" - 23" / 0 cm - 58 cm
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty

**522KL**



Built in tape measure

Hygienic ,  
easy clean tray

**522KL Display**



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
522KL	1	17.13" x 26.25" x 14.25" / 3.71	23 lb	892076002541

**524KL**  
**Digital Pediatric Seat Scale**

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb 0.5 oz  
 0 - 9 kg / 5g; 9 - 23 kg / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG conversion, LB /KG lock out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Seat Size: 18" (w) x 14-1/8" (d) x 18-1/2" (h)
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty

**524KL**



**524KL Display**



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
524KL	1	17.13" x 27" x 20.5" / 5.49	26 lb	892076002534

**553KL  
Digital Tray Scale**

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Tray Size: 25" (w) x 13-7/8" (d) x 2" (h)
- Measuring tape: 0" - 24" / 0 cm - 58 cm
- 6 AA batteries included, optional power adapter (order # ADPT30)
- Optional carrying case (order # 553CASE)
- Optional Rolling Scale Cart (order # 2210CART)
- 2 Year Limited Warranty



Carrying Case for 553KL  
Item # 553CASE

**553KL**



**2210CART  
Rolling Scale Cart for 553KL**



- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h)
- Lockable cabinet: 23" (w) x 16-3/4" (d) x 17" (h)
- Storage Shelf
- 1 Year Limited Warranty

**553KL Display**



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
553KL	1	27.75" x 19.13" x 5.25" / 1.61	12 lb	892076002329

**549KL  
Digital Tray Scale**

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- 0.6" LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out / Tare, Auto Zero, Auto Off
- Easy to clean plastic tray
- Tray Size: 19-3/8" (w) x 12-1/8" (d) x 2" (h)
- 9 Volt battery included
- Optional carrying case (order # 64769)
- Optional Rolling Scale Cart (order # 2210CART)
- 2 Year Limited Warranty



Carrying Case for 549KL  
Item # 64769

**549KL**



**2210CART  
Rolling Scale Cart for 549KL**



- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h)
- Lockable cabinet: 23" (w) x 16-3/4" (d) x 17" (h)
- Storage Shelf
- 1 Year Limited Warranty

**549KL Display**



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
549KL	1	21.75" x 15.5" x 6" / 1.17	5 lb	892076002336

# Health o meter®

## Health o meter® Infant/Toddler Scale



**Model: HDC100-01**

- 5 year limited warranty
- Includes growth chart to track height and weight against national averages
- Replaceable 9V battery
- Calibration feature allows parent to "zero out" weight of clothing or diaper
- Weighs in increments of 1/2 oz. up to 60 lb. capacity
- Large, easy to read 1.2" LCD front-facing display so parent can easily read weight
- "Grow with me" teddy bear scale: first-of-its kind consumer baby and toddler scale converts from tray to platform as child grows