

AMG Self-Taking Home Blood Pressure Kit

This easy-to-use, self-taking blood pressure kit is ideal for home use. It features a certified no-stop pin, a 300 mm Hg aneroid gauge, a self-adjust D-ring cuff and a built-in stethoscope. Comes complete with vinyl zipper case, instruction booklet, and recording chart.

Related items:

[106-410](#) AMG Self-Taking Home Blood Pressure Kit *Large adult*

Warranty:

Three year limited warranty.



106-967

Physio Logic Automatic Blood Pressure Monitor, with Large Cuff

Track your blood pressure anywhere!

With fuzzy logic technology that “personalizes” your inflation level for a more comfortable reading!

Measurements obtained with this unit are equivalent to those obtained by a trained person using the cuff/stethoscope auscultation method, as per the American National Standard, Electronic or automated sphygmomanometers.

One-Touch Operation

Operates automatically at the touch of a button, for ease of use.

85 Memories

Stores up to 85 readings, because trending is important.

- Date and time indications
- Low Battery Indicator

- Auto Power-Off 150 seconds after last button operation

Includes:

Blood Pressure Monitor, Large Arm cuff with connector, Instruction Manual, Carry Bag, and 4 x AA batteries.

Arm cuff size (measure upper arm circumference): 12.6" to 16.5" (32 cm to 42 cm)

