

# Foam Positioning Aids




Infab Corporation's experience with foam products dates back to 1970 when we were instrumental in designing new, innovative specialty seating fabrications for aircraft and helicopter pilots. Many of our designs revolutionized the industry, providing comfortable and ergonomically correct seating where none had existed previously. We have taken those many years of experience and since 1981, applied it to the wide range of products we now manufacture for the medical community. Not all foam products are alike. There is a difference in quality. We use only the finest raw materials and fabricate every piece in our factory. For the best in patient positioning, insist on Infab!

Material	Features
Super Foam AM	<ul style="list-style-type: none"> <li>• <b>Antimicrobial</b> (see next page)</li> <li>• <b>Artifact Free</b> (see next page)</li> <li>• Our exclusive formula closed cell EVA foam</li> <li>• Will not stain or absorb liquids of any kinds, including blood, barium or betadyne</li> <li>• Impervious to liquids, even when punctured or torn</li> <li>• More firm than Standard Foam</li> <li>• Excellent positioning block (Not recommended for table pads)</li> </ul>
Standard Foam	<ul style="list-style-type: none"> <li>• A first quality, fine cell, charcoal gray ester foam</li> <li>• Soft to the touch</li> <li>• Strong enough to withstand repeated washing, boiling or autoclaving</li> <li>• Named for it's reputation as the industry standard</li> <li>• One of the lowest density foams available</li> <li>• Will not shadow or artifact</li> <li>• Satisfy all patient positioning block needs</li> <li>• Exceptional choice for table pads</li> </ul>
Temper Foam	<ul style="list-style-type: none"> <li>• Originally designed for NASA and used in the Space Shuttle Program</li> <li>• Utilized for an extensive array of products where absolute comfort is a must</li> <li>• "Molds" to the body</li> <li>• Ideally suited for table pads in Cath Labs, the O.R. or any environment where extended procedures are common</li> <li>• The absolute best in patient comfort</li> </ul>

# Super Foam AM

- **Antimicrobial**

Our new "Super Foam Antimicrobial" works continuously, 24/7 to reduce the presence and transference of **MRSA** and other microbes by up to **99.97%**




- **Artifact Free**

For more complete and artifact free images. New Super-foam AM contains no glue joints and therefore is primarily artifact free at normal kVp levels.



- **Exclusive formula closed cell EVA foam**
- **Will not stain or absorb liquids of any kind, including blood, barium or betadyne**
- **Impervious to liquids, even when punctured or torn**
- **More firm, for better support and higher quality image**
- **Excellent positioning block (Not recommended for table pads)**
- **We can make any size block and any custom shape**
- **We can duplicate any manufacturers size or style**
- **See Custom Foam Form**

# Triangle




Model Number	Description	A	B	C	D	Qty.
583005	Standard Foam Triangle	10deg	2"	9"	9"	2pcs.
111005	Super Foam Triangle	10deg	2"	9"	9"	2pcs.
583010	Standard Foam Triangle	12deg	2.25"	9.50"	7.250"	2pcs.
111010	Super Foam Triangle	12deg	2.25"	9.50"	7.250"	2pcs.
583015	Standard Foam Triangle	12.5deg	2.25"	8"	10"	2pcs.
111015	Super Foam Triangle	12.5deg	2.25"	8"	10"	2pcs.
583116	Standard Foam Triangle	45deg	7"	7"	24"	2pcs.
111116	Super Foam Triangle	45deg	7"	7"	24"	2pcs.
583027	Standard Foam Triangle	15deg	7"	10"	2"	2pcs.
111027	Super Foam Triangle	15deg	7"	10"	2"	2pcs.
583025	Standard Foam Triangle	15deg	3"	9"	9"	2pcs.
111025	Super Foam Triangle	15deg	3"	9"	9"	2pcs.
583030	Standard Foam Triangle	16deg	2.50"	7"	10"	2pcs.
111030	Super Foam Triangle	16deg	2.50"	7"	10"	2pcs.
583035	Standard Foam Triangle	18deg	2.50"	8"	7"	2pcs.
111035	Super Foam Triangle	18deg	2.50"	8"	7"	2pcs.
583040	Standard Foam Triangle	18.75deg	3.75"	9"	8"	2pcs.
111040	Super Foam Triangle	18.75deg	3.75"	9"	8"	2pcs.
583045	Standard Foam Triangle	20deg	4"	9"	9"	2pcs.
111045	Super Foam Triangle	20deg	4"	9"	9"	2pcs.
583050	Standard Foam Triangle	21deg	4.625"	10"	10"	2pcs.
111050	Super Foam Triangle	21deg	4.625"	10"	10"	2pcs.

## Positioning Products

Model Number	Description	A	B	C	D	Qty.
583060	Standard Foam Triangle	25deg	5"	9"	9"	2pcs.
583080	Standard Foam Triangle	33deg	6.625"	9"	12"	2pcs.
111080	Super Foam Triangle	33deg	6.625"	9"	12"	2pcs.
583087	Standard Foam Triangle	45deg	7"	21"	7"	2pcs.
111087	Super Foam Triangle	45deg	7"	21"	7"	2pcs.
583100	Standard Foam Triangle	45deg	10"	10"	10"	1pc.
111100	Super Foam Triangle	45deg	10"	10"	10"	1pc.
583105	Standard Foam Triangle	45deg	7"	7"	10"	1pc.
111105	Super Foam Triangle	45deg	7"	7"	10"	1pc.
583110	Standard Foam Triangle	45deg	7"	7"	14"	1pc.
111110	Super Foam Triangle	45deg	7"	7"	14"	1pc.
583115	Standard Foam Triangle	45deg	7"	7"	21"	1pc.
111115	Super Foam Triangle	45deg	7"	7"	21"	1pc.
583117	Standard Foam Triangle	45deg	9"	9"	24"	1pc.
111117	Super Foam Triangle	45deg	9"	9"	24"	1pc.
583135	Standard Foam Triangle	45deg	7"	7"	7"	1pc.
111135	Super Foam Triangle	45deg	7"	7"	7"	1pc.
583140	Standard Foam Triangle	45deg	9"	9"	12"	1pc.
111140	Super Foam Triangle	45deg	9"	9"	12"	1pc.
583145	Standard Foam Triangle	45deg	4"	4"	7"	1pc.
111145	Super Foam Triangle	45deg	4"	4"	7"	1pc.
583070	Standard Foam Triangle	27deg	6"	10"	10"	2pcs.
111070	Super Foam Triangle	27deg	6"	10"	10"	2pcs.
583075	Standard Foam Triangle	30deg	6"	9"	12"	2pcs.
111075	Super Foam Triangle	30deg	6"	9"	12"	2pcs.


# Rectangle



Model Number	Description	A	B	C	Qty
583195	Standard Foam Rectangle	2	10	10	1
111195	Super Foam Rectangle	2	10	10	1
583197	Standard Foam Rectangle	6	6	2	1
111197	Super Foam Rectangle	6	6	2	1
583201	Standard Foam Rectangle	0.5, 1, 1.5 and 2	4	6	4
111201	Super Foam Rectangle	0.5, 1, 1.5 and 2	4	6	4
583205	Standard Foam Rectangle	2	8	8	1
111205	Super Foam Rectangle	2	8	8	1
583210	Standard Foam Rectangle	3	8	10	1
111210	Super Foam Rectangle	3	8	10	1
583225	Standard Foam Rectangle	2	8	12	1
111225	Super Foam Rectangle	2	8	12	1
583235	Standard Foam Rectangle	4	8	12	1
111235	Super Foam Rectangle	4	8	12	1
583240	Standard Foam Rectangle	4	10	10	1
111240	Super Foam Rectangle	4	10	10	1
583243	Standard Foam Rectangle	8	10	2	1
111243	Super Foam Rectangle	8	10	2	1
583250	Standard Foam Rectangle	4	18	24	1
111250	Super Foam Rectangle	4	18	24	1


# Cylinder

Model Number	Description	A	B
583255	Standard Foam Cylinder	4	12
583260	Standard Foam Cylinder	6	12



# Miscellaneous

Part #	Description
A	583425 Standard Foam Concave Skull Disc
B	583895 Standard Foam Tomo Skull Immobilizer
C	583760 Calcaneus strap
D	583420 Standard Foam Dual Lateral Skull Block
E	583791 Auto traction strap
F	583896 Infants Body Pediatric Immobilizer
G	583898 Infants Body Immobilizer Replacement Pads
G	583899 Childs Body Immobilizer Replacement Pads
H	583770 Superior patella strap
I	583750 Toe separator device
F	583897 Childs Body Pediatric Immobilizer
B	111895 Superfoam Foam Tomo Skull Immobilizer



F




G



H

I



E

Part #	Description
A	583270 Standard Foam Circle
A	111270 Superfoam Circle
A	583272 Standard Foam Circle Head Block
A	111272 Superfoam Foam Circle Head Block
B	583280 Standard Foam Torso Block
B	111280 Superfoam Torso Block
C	583330 Standard Foam Multiposition Head Block
C	111330 Superfoam Foam Multiposition Head Block
D	583350 Standard Foam Basal Block
D	111350 Superfoam Foam Basal Block
E	583460 Standard Foam Knee Rest

E	111460 Superfoam Foam Knee Rest
F	583285 Standard Foam Myleography Block (adult)
F	111285 Superfoam Myleography Block (adult)
F	583290 Standard Foam Myleography Block (child)
F	111290 Superfoam Myleography Block (child)
G	583355 Standard Foam Neck & Head Support
G	111355 Superfoam Foam Neck & Head Support
H	583340 Standard Foam Pressure Block
H	111340 Superfoam Foam Pressure Block
I	583371 Endovaginal U/S Wedge w/20mil



A



B



C



D



E



F



G



H



I

## Positioning Products

Part #	Description	Part #	Description
A 583375	Standard Foam AP-PA Head Immobilizer (adult)	D 111440	Superfoam Foam Dual Leg Immobilizer
A 111375	Superfoam Foam AP-PA Head Immobilizer (adult)	E 583405	Superfoam Foam Oblique Finger Block (Dual)
A 583380	Standard Foam AP-PA Head Immobilizer (child)	E 111405	Standard Foam Oblique Finger Block
A 111380	Superfoam Foam AP-PA Head Immobilizer (child)	F 583411	Standard Foam Lateral Finger Block (L/H)
B 583400	Standard Foam Oblique Finger Block (1 hand)	F 111411	Superfoam Foam Lateral Finger Block (L/H)
B 111400	Superfoam Foam Oblique Finger Block (1 hand)	G 583410	Standard Foam Lateral Finger Block (R/H)
C 583450	Standard Foam Single Leg Brace	G 111410	Superfoam Foam Lateral Finger Block (R/H)
C 111450	Superfoam Foam Single Leg Brace	H 583415	Standard Foam Lateral Finger Block (Dual)
D 583440	Standard Foam Dual Leg Immobilizer	H 111415	Superfoam Foam Lateral Finger Block (Dual)




A



B



C



D



E



F




G



H


### General X-Ray Kit



Contains: 2 each 45 deg Triangle 3" x 3" x 7"  
 1 each 45 deg Triangle 7" x 7" x 7"  
 1 each Circle 1.5" x 7"  
 2 each 12 deg Triangle 2.25" x 9.5" x 7"  
 1 each Rectangle 2" x 8" x 10"

Model Number	
INF-1	Standard Foam
INF-1SF	Super Foam


### General X-Ray Kit 2



Contains: 1 each Circle 1.5" x 7"  
 1 each Rectangle 2" x 8" x 10"  
 2 each 34 deg Triangle 3" x 4" x 7"  
 2 each 18 deg Triangle 2.25" x 8" x 7"  
 1 each 45 deg Triangle 7" x 7" x 7"

Model Number	
INF-2	Standard Foam
INF-2SF	Super Foam


### General X-Ray Kit 3



Contains: 1 each Circle 1.5" x 7"  
 1 each Rectangle 2" x 8" x 10"  
 2 each 34 deg Triangle 3" x 4" x 7"  
 2 each 18 deg Triangle 2.25" x 8" x 7"  
 1 each 45 deg Triangle 7" x 7" x 7"  
 1 each Oblique finger block, 1 hand

Model Number	
INF-3	Standard Foam
INF-3SF	Super Foam

### General X-Ray Kit 4



Contains: 1 each Circle 1.5" x 7"  
 1 each Rectangle 2" x 8" x 10"  
 2 each 34 deg Triangle 3" x 4" x 7"  
 1 each 45 deg Triangle 7" x 7" x 7"  
 1 each Oblique finger block, 1 hand  
 2 each 22 deg Triangle 3.5" x 7" x 8"

Model Number	
INF-4	Standard Foam
INF-4SF	Super Foam

## Positioning Products


### Deluxe X-Ray Kit



Contains: 2 each 45 deg Triangle 3" x 7" x 7"  
 1 each 45 deg Triangle 7" x 7" x 7"  
 1 each Circle 1.5" x 7"  
 2 each 12 deg Triangle 2.25" x 9.5" x 7"  
 1 each Rectangle 2" x 8" x 10"  
 1 each Rectangle 4" x 18" x 24"  
 1 each Table Pad 1" x 24" x 72"

Model Number	
INF-5	Standard Foam
INF-5SF	Super Foam


### Trauma Kit



Contains: 1 each Basal Block  
 1 each AP/PA Skull Block  
 1 each Pressure Block  
 2 each 18.75 deg Triangle 3.75" x 9" x 8"  
 1 each 45 deg Triangle 6.75" x 6.75" x 12"  
 1 each Rectangle 3" x 8" x 10"  
 1 each 27 Torso Block

Model Number	
INF-6	Standard Foam
INF-6SF	Super Foam

### Trauma Kit 7



Contains: 1 each Basal Block  
 1 each AP/PA Skull Block (adult)  
 1 each Pressure Block  
 1 each Rectangle 4" x 10" x 10"  
 2 each 12.5 Triangle 2.25" x 8" x 10"  
 2 each 21 deg Triangle 3.75" x 8" x 10"

Model Number	
INF-7	Standard Foam
INF-7SF	Super Foam


### Special Procedure Kit



Contains: 1 each Waters Positioner  
 1 each Tomo Skull Immobilizer  
 1 each AP/PA Head Immobilizer  
 1 each AP/PA Skull Block

Model Number	
INF-8	Standard Foam
INF-8SF	Super Foam


### Triangle Kit



Contains: 2 each 21 deg Triangle 4.5" x 10" x 10"  
 2 each 33 deg Triangle 6.67" x 9" x 12"  
 1 each 45 deg Triangle 10" x 10" x 10"  
 2 each 27 deg Triangle 6" x 10" x 10"

Model Number	
INF-9	Standard Foam
INF-9SF	Super Foam


### Deluxe Kit 10



Contains: 1 each AP/PA Skull Block (adult)  
 1 each Pressure Block  
 2 each 13.50 deg Triangle 3" x 10" x 10"  
 2 each 21 deg Triangle 4.67" x 10" x 10"  
 1 each 45 deg Triangle 10" x 10" x 10"  
 1 each Rectangle 4" x 10" x 10"  
 2 each 12.5 deg Triangle 2.25" x 8" x 10"  
 2 each 21 deg Triangle 3.75" x 8" x 10"  
 2 each 27 deg Triangle 6" x 10" x 10"  
 1 each Basal Block

Model Number	
INF-10	Standard Foam
INF-10SF	Super Foam

### Orthopedic Kit



Contains: 1 each 27 deg Torso Block  
 1 each AP/PA Skull Block (adult)  
 1 each Pressure Block  
 1 each Rectangle 2" x 8" x 10"  
 1 each Oblique Finger Block (adult)

Model Number	
INF-11	Standard Foam
INF-11SF	Super Foam


### Orthopedic Kit 12



Contains: 1 each AP/PA Skull Block (adult)  
 1 each Oblique Finger Block (adult)  
 1 each Cervical Rest  
 1 each 45 deg. Triangle 7" x 7" x 21"  
 1 each Rectangle 3" x 8" x 10"

Model Number	
INF-12	Standard Foam
INF-12SF	Super Foam


**Orthopedic Kit 13**



Contains: 1 each AP/PA Skull Block (child)  
 1 each Neck and Head Support  
 1 each Cervical Rest  
 2 each 12.5 deg Triangle 2.25" x 8" x 10"  
 2 each 21 deg Triangle 4.5" x 10" x 10"  
 1 each 45 deg Triangle 7" x 7" x 10"  
 1 each Rectangle 3" x 8" x 10"

Model Number	
INF-13	Standard Foam
INF-13SF	Super Foam


**Deluxe Mammo Kit**



Contains: 1 each 8 deg Triangle 1.5" x 10" x 10"  
 1 each 10 deg Triangle 1.5" x 8" x 8"  
 1 each 11 deg Triangle 2" x 10" x 10"  
 1 each 14 deg Triangle 2.5" x 10" x 10"  
 1 each 15 deg Triangle 2" x 8" x 8"  
 1 each 16 deg Triangle 3" x 10" x 10"  
 1 each 19 deg Triangle 3.5" x 10" x 10"  
 1 each 20 deg Triangle 3" x 8" x 8"  
 1 each 21 deg Triangle 4" x 10" x 10"  
 1 each 24 deg Triangle 3.5" x 8" x 8"  
 1 each Rectangle 1" x 10" x 10"  
 1 each Rectangle 2" x 10" x 10"  
 1 each Rectangle 3" x 10" x 10"

Model Number	
INF-14	Standard Foam
INF-14SF	Super Foam

**Mammo Triangle Kit**



Contains: 1 each 10 deg Triangle 1.5" x 8" x 8"  
 1 each 15 deg Triangle 2" x 8" x 8"  
 1 each 20 deg Triangle 3" x 8" x 8"  
 1 each 24 deg Triangle 3.5" x 8" x 8"

Model Number	
INF-15	Standard Foam
INF-15SF	Super Foam

**Mammo Rectangle Kit**



Contains: 1 each Rectangle 1" x 10" x 10"  
 1 each Rectangle 2" x 10" x 10"  
 1 each Rectangle 3" x 10" x 10"

Model Number	
INF-16	Standard Foam
INF-16SF	Super Foam


# Custom Foam Form

(Can also be downloaded from [www.infabcorp.com](http://www.infabcorp.com))


Can't find the size or shape you are looking for from our wide selection of stock blocks? No problem. Just fill in the form below, and we'll provide you a high quality positioning block cut exactly to your specifications.



Rectangle



Triangle



Cylinder

Block Type	Dimension A	Dimension B	Dimension C	Dimension D	Standard or Super Foam
Rectangle					
Triangle					
Cylinder					
Other					

**Special instructions:**

---



---



---



---

# Sandbags

Model	Description	Size
583840	Cervical Sandbag Set	11" x 11"
583850	3 Pound Sandbag	7" x 7"
583851	5 Pound Sandbag	7" x 12"
583860	7 Pound Sandbag	6" x 17"
583865	10 Pound Sandbag	7" x 20"
583870	12 Pound Sandbag	7" x 24"
583871	15 Pound Sandbag	11" x 15"
583881	20 Pound Sandbag	8.5" x 22"



## Sandbags Set:

Model Number	Description	Qty	Size	Weight
583890	Pediatric Set Contains:	3 (per set)	6" x 6"	3 lbs

Model Number	Description	Qty	Size	Weight
583882	Sandbag Set Contains:	1 (per set)	7" x 12"	5 lbs
		1 (per set)	7" x 20"	10 lbs
		1 (per set)	11" x 15"	15 lbs
		1 (per set)	8" x 22"	20 lbs

Model Number	Description	Qty.	Size	Weight
583884	Sandbag Set Contains	1 (per set)	7" x 7"	3 lbs
		1 (per set)	7" x 12"	5 lbs
		1 (per set)	6" x 17"	7 lbs
		1 (per set)	7" x 20"	10 lbs
		1 (per set)	7" x 24"	12 lbs
		1 (per set)	12" x 12"	13 lbs
		1 (per set)	11" x 15"	15 lbs
		1 (per set)	8" x 22"	20 lbs