



Invacare® Heavy-Duty Dual-Release Paddle Walker And Bariatric Rollator

Invacare® Heavy-Duty Dual-Release Paddle Walker and Bariatric Rollator

Specifications

(all specifications are approximate)

Model no. 6291-HDA	
Width inside handgrips	20"
Width inside base legs	24"
Depth opened at base	17.5"
Depth folded	4"
Height adjustment (min. -max.)	30.38" - 37.38" (adjusts in 1" increments)
Approx. patient height (min. - max.)	5'0" - 6'3"
Product weight	6 lb. 11 oz.
Weight capacity	500 lb.
Limited warranty	lifetime
Quantity	2 per carton

Features

- 500 lb. weight capacity
- Extended frame for extra stability
- Long-lasting PVC handgrips



Model no. 6291-HDA

The Invacare Heavy-Duty Dual-Release Adult Paddle Walker is strong, stable and simple to use, helping individuals to walk with confidence and ease. This walker also is compatible with Invacare leg accessories.



NEW Model no. 66500

The Invacare Bariatric Rollator is made of a durable and stable steel frame designed for individuals up to 500 lb. who need an extra wide and deep ambulatory device. The flip-up padded seat allows quick access to the storage basket while the 8" wheels provide improved mobility over rougher indoor or outdoor terrains.

Features

- 500 lb. weight capacity
- Wide, deep frame
- Flip-up padded seat
- Removable padded backrest
- Folds for ease of storage and transport
- Large 8" wheels for improved mobility over rougher terrain
- Height adjustable loop handbrakes with foam handgrips
- Durable steel construction
- Basket included

Specifications

(all specifications are approximate)

Depth opened at base	32.875"	Seat width x depth	18" x 13"
Depth folded	19.5"	Wheels	8"
Width inside handgrips	22.50"	Product weight	26 lb.
Width inside base of legs	21.875"	Weight capacity	500 lb.
Width overall	31.125"	Quantity	1 per carton
Height adj. (min. -max)	34.25" - 38.50"	Limited warranty	
Approx. client height (min.-max)	5'8" - 6'5"	frame	3-years
Seat-to-floor height	21.625"	brakes	3-years
		other components	3-years

Additional Product Suggestions

Additional product suggestions are for recommendation only. Please refer to the individual product specification to verify weight capacity.



1. Bariatric Offset Handle Cane #8940-6, page 28
2. Bariatric Shower Chair w/back #9785-2, page 29
3. Bariatric Drop Arm Commode #9785-2, page 32

Up to 500 lb

Invacare Customer Service (800) 333-6900 www.invacare.com Invacare Bariatric Price List Form no. 99-277



Invacare® Bariatric Quad Canes and Offset Cane

NEW Model nos. 7831-2, 7835-2, 8940-6

Invacare Bariatric Canes are built with durable steel construction that provides maximum stability for individuals weighing up to 700lb. All of the canes are height adjustable for a proper fit to the consumer and are designed with a twist locking ring providing additional security during height adjustment.

Features

- 700 lb. weight capacity
- Right- or left-hand use
- Height adjustable
- Comfortable PVC hand grip
- Locking ring for additional stability
- Durable steel construction



Small base quad cane model no. 7831-2



Large base quad cane model no. 7835-2

Specifications

(all specifications are approximate)

Model no.	7831-2	7835-2
Base	Small	Large
Color	black	black
Adjustment (min. -max.)	29.50" - 38.50"	29.50" - 38.50"
Client height (approx.)	4'11" - 6'5"	4'11" - 6'5"
Footprint	8" x 6"	12" x 8.50"
Weight	3 lb.	4 lb.
Weight capacity	700 lb.	700 lb.
Quantity	2 pair per carton	2 pair per carton
Limited warranty	3-years	3-years



Offset base quad cane model no. 8940-6

Features

- 700 lb. weight capacity
- Height adjustable
- Comfortable PVC hand grip
- Locking ring for additional stability
- Durable steel construction

Specifications

(all specifications are approximate)

Model no.	8940-6
Style	Offset
Color	black
Adjustment (min. -max.)	29" - 38"
Client Height (approx.)	4'10" - 6'4"
Weight	3 lb.
Weight capacity	700 lb.
Quantity	6 pair per carton
Limited warranty	3-years

Additional Product Suggestions

Additional product suggestions are for recommendation only. Please refer to the individual product specification to verify weight capacity.



1. Bariatric Walker #6441-A, page 30
2. Bariatric Transfer Bench #9675-2, page 29
3. Bariatric Drop-Arm Commode #6599, page 32

Invacare Customer Service (800) 333-6900 www.invacare.com Invacare Bariatric Price List Form no. 99-277

Up to 700 lb.

Invacare® Bariatric Dual Release Walker and Walker Kit



Invacare® Bariatric Dual-Release Walker

Specifications

(all specifications are approximate)

Dimensions	
Width inside handgrips	20"
Width inside base legs	23.50"
Depth opened inside base	16"
Depth folded	4.25"
Height adjustment (min. - max.)	31.25" - 41.25"
Approx. client height (min. - max.)	5'3" - 6'10"
Product weight	7 lb.
Weight capacity	700 lb.
Quantity	2 per carton
Limited warranty	3-years

Features

- 700 lb. weight capacity
- Wide, deep frame
- Dual-push button mechanisms provide visual and audible "locked" cues
- Large number of height adjustments
- Comfortable aluminum frame
- Compatible with Invacare model no. 6372 Bariatric 5" fixed wheels and leg extensions




NEW

Model no. 6441-A

The Invacare Bariatric Adult Walker provides a wide and deep frame that can support an individual up to 700 lb. The dual-release buttons provide both visual and audible "locked" cues for the user. The durable aluminum frame provides maximum stability for consumers and the legs have a large number of height adjustments.

NEW

Model no. 6372



Features

- Only to be used with Invacare bariatric walker model no. 6441-A
- Dual wheels for added stability
- Large number of height adjustments

Specifications

(all specifications are approximate)


Dimensions	
Height adjustment (min. - max.)	33.50" - 43.50"
Approx. client height (min. - max.)	5' 6" - 6' 10"
Weight capacity	700 lb.
Quantity	1 per carton
Limited warranty	3 years

Additional Product Suggestions

Additional product suggestions are for recommendation only. Please refer to the individual product specification to verify weight capacity.



1.



2.



3.

1. Bariatric Bed (1000 lb.) #BAR1000, page 34
2. Bariatric Small Base Quad Cane #7831-2, page 28
3. Bariatric Shower Chair #9875-2, page 29

Invacare Customer Service (800) 333-6900 www.invacare.com Invacare Bariatric Price List Form no. 99-277

Up to 700 lb.