

The Invacare® WalkLite™ Walker revolutionized the world of walkers by addressing the needs of providers, clinicians and consumers at an affordable price. The line of walkers covers a wide range of consumer heights from Junior to Tall Adult. The WalkLite models continue to feature all the benefits of the original WalkLite – a comfortable flip-up seat, durable front 5" fixed wheels, rear spring-loaded SureGlide™ brakes, a lightweight aluminum frame with flexible backrest and ergonomic dual-paddle folding mechanisms.



Model no. 6300-JRA
(Junior to Adult)



Model no. 6300-ATA
(Adult to Tall)



Folds into a compact unit for storage. Built-in seat handle allows for easy transport.

Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanisms
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Rear spring-loaded SureGlide brakes
- Adjustable seat-to-floor height

Basket and Tray Accessories

- Basket (model no. 65110 (shown) for 6300-ATA; 65109 for 6300-JRA)
 - 10 lb. weight capacity
 - Position basket underneath seat to secure valuables
- Tray (model no. 65408)
 - For food and beverage placement
 - Tray fits on both 6300-JRA and 6300-ATA



Specifications

	6300-JRA	6300-ATA
Depth opened at base	23"	23"
Depth folded	7"	7"
Width inside handgrips	18"	18"
Width inside base of legs	20"	24 ½"
Width overall	27"	27"
Height adjustment min. - max.	29 ¼" - 36 ½"	32" - 39"
Approx. client height* min. - max.	4' 11" - 6' 0"	5' 5" - 6' 5"
Seat-to-floor height min. - max.	16 ⅝" - 23 ⅝"	19" - 25 ¾"
Seat width x depth	17" x 9 ½"	17" x 9 ½"
Wheels	5" fixed	5" fixed
Product weight	11 lb.	11 lb.
Weight capacity	300 lb.	300 lb.
Units/carton	1	1
Limited warranty - frame	lifetime	lifetime
Limited warranty - other components	1 year	1 year

*Client height is approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment.



Flip-up seat allows user to step inside the frame.



*Rear spring-loaded SureGlide brakes provide additional security when seated.***

** The rear SureGlide brakes engagement may vary depending upon the weight distributed on the seat and the positioning of the user on the seat.

Invacare contemporary walkers feature wide, deep frames with a large number of height adjustments and a lower side brace for added stability. They are light and easy to lift and maneuver. Walkers include PVC handgrips that provide comfort and long-lasting wear and anti-rattle silencers for quiet operation.



Invacare® Single-Release Walkers

Model no. 6281-A - Single-Release Paddle Adult
Model no. 6281-JR - Single-Release Paddle Junior



Invacare® Blue-Release™ Walkers

Model no. 6240-A - Dual Blue-Release Adult
Model no. 6240-JR - Dual Blue-Release Junior



Invacare® I•Class™ Paddle Walkers

Model no. 6291-A - Dual-Release Paddle Adult (shown)
Model no. 6291-1 - Dual-Release Paddle Adult
Model no. 6291-JR - Dual-Release Paddle Junior




Invacare® Wheeled Walkers

Model no. 6291-3F - Dual-Release Paddle Adult Walker
Model no. 6291-JR3F - Dual-Release Paddle Junior Walker
Model no. 6291-5F - Dual-Release Paddle Adult Walker
Model no. 6291-JR5F - Paddle-Release Adult Walker - **New**
Model no. 6240-5F - Dual Blue-Release™ Adult Walker (shown)
Model no. 6240-JR5F - Dual Blue-Release Junior Walker

Features

Invacare® Folding Mechanisms

Our folding mechanisms make it easy to open and close the walkers even from a seated position. Walkers fold easily for transport or storage.



Single-Release Button

- One button to press for folding
- Button provides both visual and audible “locked” cues



Dual-Release Paddle

- Ergonomic paddles “pop up” when walker is opened
- Nylon paddles are resistant to cracking or shattering



Dual-Release Button

- Dual Blue-Release™ buttons provide both visual and audible “locked” cues



Invacare® Heavy-Duty Paddle Walker

Model no. 6291-HDA - Heavy-Duty Dual-Release Paddle Adult

- Compatible with most Invacare leg accessories



Invacare® Bariatric Dual Release Walker

Model no. 6441-A

- Wide, deep frame
- Dual-push button mechanisms provide visual and audible "locked" cues
- Large number of height adjustments
- Comfortable aluminum frame
- Compatible with Invacare model no. 6372 bariatric 5" fixed wheels and leg extensions



Invacare® Adult Hemi Walker/Sidestepper

Model no. 6252 - Hemi Walker

- Wider base provides more stability than a quad cane
- 14 ¼" x 17 ¾" width at high height setting
- 12 ¼" x 16 ½" width at low height setting

Specifications

Model no.		Width inside hand-grips	Width inside base legs	Depth opened at base	Depth folded	Height adjustment min. - max.	*Approx. client height min. - max.	Product weight	Weight capacity	Units/ carton	Limited warranty
6291-A (4); 6291-1	Dual-Release Paddle Adult	17"	21 ½"	17 ½"	4"	30 ¾" - 37 ¾"	5' 3" - 6' 4"	5 lb. 5 oz.	300 lb.	4	lifetime
6291-JR	Dual-Release Paddle Junior	17"	20 ½"	16 ½"	4"	25 ¾" - 32 ¾"	4' 4" - 5' 7"	5 lb. 1 oz.	300 lb.	4	lifetime
6291-3F	Wheeled Dual-Release Paddle Adult	17"	21 ¼"	17 ½"	4 ½"	32" - 38"	5' 4" - 6' 4"	5 lb. 5 oz.	300 lb.	4	lifetime
6291-JR3F	Wheeled Dual-Release Paddle Junior	17"	20 ¼"	16 ½"	4 ½"	27" - 33"	4' 6" - 5' 6"	5 lb. 1 oz.	300 lb.	4	lifetime
6291-5F**	Wheeled Dual-Release Paddle Adult	17"	21 ½"	17 ½"	4 ½"	33" - 39"	5' 6" - 6' 6"	5 lb. 6 oz.	300 lb.	4	lifetime
6291-JR5F**	Wheeled Dual-Release Paddle Junior	17"	20 ¼"	16 ½"	4 ½"	28" - 34"	4' 4" - 5' 6"	5 lb. 6 oz.	300 lb.	4	1 year
6291-HDA	Heavy-duty Dual-Release Paddle Adult	20"	24"	17 ½"	4"	30 ¾" - 37 ¾"	5' 3" - 6' 4"	6 lb. 11 oz.	500 lb.	2	lifetime
6441-A	Bariatric Dual-Release	20"	23 ½"	16"	4 ¼"	31 ¼" - 41 ¼"	5' 3" - 6' 10"	7 lb.	700 lb.	2	3 years
6240-A	Dual Blue-Release Adult	17"	21"	16"	4"	32" - 39"	5' 4" - 6' 6"	5 lb. 7 oz.	300 lb.	4	1 year
6240-JR	Dual Blue-Release Junior	17"	20"	16"	4"	25" - 32"	4' 4" - 5' 7"	5 lb. 5 oz.	300 lb.	4	1 year
6240-5F	Wheeled Dual Blue-Release Adult	17"	21"	18 ½"	4"	32" - 39"	5' 4" - 6' 6"	6 lb. 6 oz.	300 lb.	4	1 year
6240-JR5F	Wheeled Dual Blue-Release Junior	17"	20"	18"	4"	25" - 32"	4' 4" - 5' 7"	6 lb. 12 oz.	300 lb.	4	1 year
6281-A	Single-Release Adult	17"	21"	17 ½"	4 ½"	30 ¾" - 37 ¾"	5' 3" - 6' 4"	6 lb. 2 oz.	300 lb.	4	1 year
6281-JR	Single-Release Junior	17"	20 ½"	16 ½"	4 ½"	25 ¾" - 32 ¾"	4' 4" - 5' 7"	5 lb. 14 oz.	300 lb.	4	1 year
6252	Hemi Walker Adult	n/a	n/a	n/a	3"	29 ½" - 34 ½"	5' 1" - 5' 10"	2 lb. 14 oz.	250 lb.	2	3 years

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

"F" designation in Model no. 3" or 5" fixed wheels with glide tips. ** With 13" rear leg extensions.

Court-Side Glides™



Model no. 6290: Sold as a carton; 6 pairs per carton. Each pair complete with aluminum leg extensions, housings, and tennis balls.

Invacare® Court-Side Glides™ give consumers durable walker glide tips that are long-lasting, safe and easy to install, and provide a smooth glide experience. Consumers can install their Court-Sides Glides in seconds on Invacare® fixed-wheeled walkers and WalkLite™ models.* The Court-Side Glides feature a specially designed housing that allows each ball to rotate so it lasts longer. The Court-Side Glides have spring-loaded brakes for additional security and lock-out clips for pure gliding. The tennis balls pop in and out for quick and effortless replacement.

*Court-Side Glides change the approximate client height range of walkers (see table).

Do not use on Model no. 6291HDA

Client-Walker Height Range Table

Walker Model No.*	New Approx. Client Height Range w/Court-Side Glides*
6291, 6240, 6281 w/ 5" fixed wheels	5' 6" - 6' 4"
6291, 6240, 6281 w/ 3" fixed wheels	5' 4" - 6' 4"
6291JR, 6240JR, 6281JR w/ 5" fixed wheels	4' 8" - 5' 6"
6291JR, 6240JR, 6281JR w/ 3" fixed wheels	4' 6" - 5' 6"
6300-JRA Adult to Junior WalkLite	4'11" - 5' 6"
6300-ATA Tall Adult WalkLite	5' 4" - 6' 0"

All client height ranges are approximate. Do not use on 6291HDA

Invacare® Walker Accessories

Invacare® walker accessories enhance walkers for safety, maneuverability and durability. Wheel rims are made of durable, easy-to-clean polypropylene; tires are made of non-marring thermoplastic rubber; wheels and bearings are corrosion protected.

Swivel Wheel Attachments



5" Swivel Wheels with 13" Cambered SureGlide™ Brakes

Model no. 6268



5" Swivel Wheel Attachments with Rear Glide Tips

Model no. 6267



3" Swivel Wheels with 13" Cambered SureGlide™ Brakes

Model no. 6272



3" Swivel Wheel Attachments with Rear Glide Tips

Model no. 6273

Fixed-Wheel Attachments



3" Single Fixed-Wheel Attachments with Rear Glide Tips


Model no. 6270**



5" Single Fixed Wheel Attachments with Rear Glide Tips

Model no. 6271**

- Use 13" rear leg extensions for maximum adjustability



Tall 3" Single Fixed-Wheel Attachments with Tall Rear Legs and Glide Tips

Model no. 6278 with 13" rear leg extensions

- Adjusts to accommodate persons up to 6' 6" tall



Bariatric Wheel Kit

Model no. 6372 only to be used with Invacare bariatric walker 6441-A

- Dual 5" fixed wheels for added stability
- Large number of height adjustments

**Also available in retail packaging with a minimum purchase of four sets. Order as Model nos. 6270-4 and 6271-4.

Rear Brake Attachments

Brake attachments operate without interfering with regular walking motion.



10" Rear Wheel Brakes

Model no. 6264

- Use with 3" or 5" single fixed-wheel attachments

13" Rear Wheel Brakes

Model no. 6265 (not shown)

- Use with 3" or 5" swivel wheel attachments or tall 3" single fixed-wheel attachments

Leg Extensions



Tall Leg Extensions

Model no. 6275

- Adjusts to accommodate persons up to 6' 6" tall

13" Rear Leg Extensions with Glide Tips

Model no. 6266

- Use with swivel or fixed wheels



Glide Tips

Model no. 6280

- Packaged 24 per carton
- Durable construction for long wear

Walker Accessory Specifications*

*Walker accessories are not to be used with WalkLite™ Walker.

Model no.	Height adjustment			Limited warranty
	Adult walkers 6291-1, 6291-3F, 6291-5F, 6291-A, 6281-A, 6240-A, 6240-5F, 6291-HDA, 6441-A	Junior walkers 6291-JR, 6281-JR, 6291-JR3F, 6291-JR5F, 6240-JR, 6240-JR5F		
Fixed-wheel				
6270	3" single fixed-wheels w/rear glide tips	31" - 37"	26" - 32"	lifetime
6271	5" single fixed-wheels w/rear glide tips	33" - 39"	28" - 33"	lifetime
6277	tall 3" single fixed-wheels w/glide tips	31" - 37"	26" - 32"	lifetime
6278	tall 3" single fixed-wheels w/tall rear leg ext. w/glide tips	33" - 40"	31" - 35"	lifetime
6372*	double 5" fixed-wheels for 6441-A bariatric walker	33 ½" - 43 ½"	—	3 years
Swivel wheel				
6267	5" swivel wheels w/rear glide tips	37" - 40"	32" - 35"	lifetime
6268	5" swivel wheels w/13" cambered SureGlide™ brakes	37" - 40"	32" - 35"	lifetime
6269	3" swivel wheels w/13" rear leg extensions w/glide tips	34" - 40"	29" - 35"	lifetime
6272	3" swivel wheels w/13" cambered SureGlide brakes	34" - 40"	29" - 35"	lifetime
6273	3" swivel wheels w/rear glide tips	34" - 37"	29" - 32"	lifetime
Leg extensions				
6266	13" rear leg extensions w/glide tips	33" - 40"	29" - 35"	lifetime
6275	tall front and rear leg extensions	33" - 40"	28" - 35"	lifetime
Rear brake attachments				
6264	10" rear wheel brakes	33" - 37"	28" - 32"	lifetime
6265**	13" rear wheel brakes	33" - 40"	29" - 35"	lifetime

**For use with 3" or 5" swivel wheels or tall 3" fixed wheels.

The Invacare Walker Platform attachment is now even easier to attach and detach. Invacare makes a good product even better with a 2" clamp-on and tighten feature that does not require tools to assemble. The clamp-on feature also helps with ease of adjustment when determining correct height.



Model no. 6027



Platform Attachment

Model no. 6027
(Compatible with all walkers except model nos. 6291-HDA, 6441-A, 6252 or 6300-JRA, 6300-ATA)

- Easy tool-free installation
- For adult and junior walkers
- Thick, vinyl-covered pad with heavy-duty Velcro® strap
- 300 lb. weight capacity
- One per carton
- One year limited warranty



Assistance Attachments

All attachments mount easily to most Invacare® walkers.



Walker Basket

Model no. 6098
(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA, 6240 walkers)

- Mounts easily with snap hooks
- Rust-resistant, gray vinyl coating
- Basket can be installed on the inside or outside of the walker
- 5 lb. weight capacity
- One per carton
- One year limited warranty



Walker Tray

Model no. 6007
(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA, 6240, 6281 walkers)

- Provides convenient assistance for a variety of everyday activities
- Molded with details for cups and pen holders
- Folds out of the way when not in use
- 5 lb. weight capacity
- One per carton
- One year limited warranty



Walker Pouch

Model no. 6015-4
(Compatible with all walkers except model nos. 6291-HDA, 6252, 6300-JRA, 6300-ATA)

- Made of forest green/deep sea blue denier nylon
- Inside and outside pockets
- Easily attaches over front crossbar
- Lightweight and water-repellent
- 5 lb. weight capacity
- Four per carton
- One year limited warranty

NEW Knee Walker



Model no. 65950

Introducing the new Invacare® Knee Walker. The Knee Walker is a great alternative for those who are weary of the discomfort related to crutches. It features a comfortable knee pad for resting the wounded limb, and it is easily maneuverable with its five-inch front swivel wheels. Its easy-to-operate ergonomic hand brakes keep the consumer in control. In addition, the Knee Walker folds easily for storage or transport, and it comes with a basket for the consumer's convenience.



Specifications

	65950
Depth opened	29 1/2"
Width opened	18"
Adj. handle height min.-max.	30" - 38"
Client approx. height min.-max.	5' - 6' 3"
Wheels	5"
Seat width x depth	7" x 15 1/2"
Seat-to-floor height	17" - 20 1/2"
Product weight	16.5 lb.
Weight capacity	300 lb.
Limited warranty	
frame	3 years
other components	1 year
Units/carton	1

NEW Value-Line Rollators

Introducing the two newest rollators to Invacare's Value-Line. These rollators will meet all your needs with a focus on style and affordability. These rollators feature easy-to-operate ergonomic hand brakes and easy folding for storage and transport. They offer great maneuverability and stability, and are height adjustable for a custom fit. Basket is standard on both.



65350GR Junior Rollator Features

Model no. 65350GR

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up seat
- Lightweight aluminum frame
- Basket hooks easily on front seat brace
- Removable foam padded straight backrest enhances comfort and can be easily removed



65851R Adult Rollator Features

Model no. 65851R

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up seat
- Lightweight aluminum frame
- Basket hooks easily on front seat brace
- Removable foam-padded straight backrest enhances comfort and can be easily removed

Other Value-Line Rollators include:



65410 Features

Model no. 65410

- Easy-to-operate, ergonomic hand brakes
- Rollator folds easily with tote bag attached
- Indoor/outdoor 8" x 2" solid rubber wheels
- Basket standard, attaches easily to the frame
- Tray standard, attaches easily to the basket
- Tote bag standard, attaches easily to the frame

65450 Features

Model no. 65450

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Sturdy seat
- Rollator folds easily with basket attached
- Basket and tray standard

65650 / 65650R Features

Model no. 65650 Green, 65650R Red

- Easy-to-operate, ergonomic hand brakes lock the rear wheels for security
- Flip-up padded seat with zipper compartment for storage of valuables
- Lightweight aluminum frame
- Removable foam padded curved backrest enhances comfort and can be easily removed
- Rollator folds easily with basket attached
- Universal basket fits underneath seat or hooks easily on front seat brace

66500 Bariatric Rollator Features

Model no. 66500

- 500 lb. weight capacity
- Wide, deep frame
- Flip-up padded seat
- Removable padded backrest
- Folds for ease of storage and transport
- Large 8" wheels for improved mobility over rougher terrain
- Height adjustable loop handbrakes with foam handgrips
- Durable steel construction
- Basket included



Specifications

Baskets are standard on model nos. 65410, 65350GR, 65450, 65650, 65650R, 65851R and 66500. Tray is standard on model nos. 65410 and 65450.

	65410	65410 Tote bag	65350GR	65450	65650/65650R	65851R	66500
Depth opened	22.85"	14 3/16" depth	24 1/2"	26 1/2"	28 3/32"	27"	27"
Width opened	26"	14 3/16" width	26"	23 1/2"	24 1/4"	28 1/2"	31 1/4"
Length		10 5/8"					
Adj. handle height min.-max.	31 1/2" - 36 1/4"		29 1/2" - 34"	31 3/8" - 38 1/8"	32 1/16" - 37 5/16"	34" - 38"	34" - 39"
Client approx. height min.-max.	5' 3" - 6'		5' - 5' 6"	5' 3" - 6' 5"	5' 5" - 6' 2"	5' 6" - 6' 3"	5' 8" - 6' 5"
Distance between handles	14"		13 1/2"	14 3/4"	18"	16"	22 1/2"
Wheels	8" x 2"		6"	7 1/2"	7 1/2"	6"	8"
Seat width x depth	N/A		14" x 13"	15" x 6 3/8"	13 3/4" x 12 1/2"	14" x 13"	19 1/4" x 13 1/4"
Seat-to-floor height	N/A		19 1/2"	24 3/8"	23 1/4"	22"	22"
Product weight	17 lb.		13 lb.	25 lb.	19 lb.	13.4 lb.	26 lb.
Weight capacity	264 lb.	11 lb.	250 lb.	300 lb.	300 lb.	300 lb.	500 lb.
Limited warranty							
frame	3 years	1 year	3 years	3 years	3 years	3 years	3 years
other components	1 year		1 year	1 year	1 year	1 year	3 years
Units/carton	1	standard on 65410	1	1	1	1	1

The Invacare line of Rollite® Rollators is unlike any other rollator line in the industry. Invacare has taken the Rollite one step further, offering the sizes, wheel types, and colors that consumers want. The Rollite line offers adult and junior sizes in different colors. For improved mobility over rougher terrain, large 8" wheel models are available in adult and tall sizes (Midnight Blue color). Along with additional features, the Rollite line provides the comfort, durability and quality that consumers have come to expect from Invacare products.



Color: Electric Red
Model nos. 65100R Rollite Adult,
65100R-JR Rollite Junior (shown)



Color: Midnight Blue
Model nos. 65100 Rollite Adult (shown)
65100-JR Rollite Junior



Color: Midnight Blue
Model nos. 68100 Rollite Adult w/8" wheels,
68100-TA Rollite Tall Adult w/8" wheels (shown).

Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

Basket and Tray Accessories

- Basket (model no. 65110) For use on Rollite models 65100, 65100R and 68100-TA
- Tray (model no. 65408)
 - For food and beverage placement
 - For use on all Rollite models
- Basket (model no. 65109) For use on Rollite models 65100-JR, 65100R-JR and 68100
 - 10 lb. weight capacity
 - Position basket underneath seat to secure valuables



Folds into a compact unit for storage. Built-in seat handle allows for easy transport.

Specifications

	65100R Rollite Adult (Electric Red)	65100R-JR Rollite Junior (Electric Red)	65100 Rollite Adult (Midnight Blue)	65100-JR Rollite Junior (Midnight Blue)	68100 Rollite Adult w/8" wheels (Midnight Blue)	68100-TA Rollite Tall Adult w/8" wheels (Midnight Blue)
Depth opened at base	24"	24"	24"	24"	26"	26"
Depth folded	7"	7"	7"	7"	7"	7"
Width inside handgrips	18"	18"	18"	18"	18"	18"
Width inside base of legs	22"	22"	22"	22"	22"	22"
Width overall	27"	27"	27"	27"	27"	27"
Height adjustment min. - max.	32" - 36"	29" - 33"	32" - 36"	29" - 33"	32" - 36"	35" - 39"
*Approx. client height min. - max.	5' 4" - 6'	4' 11" - 5' 7"	5' 4" - 6'	4' 11" - 5' 7"	5' 4" - 6'	5' 11" - 6' 6"
Seat-to-floor height min. - max.	19" - 23"	16 ½" - 20 ½"	19" - 23"	16 ½" - 20 ½"	19 ½" - 23 ½"	22" - 26"
Seat width x depth	17" x 9 ½"	17" x 9 ½"	17" x 9 ½"	17" x 9 ½"	17" x 9 ½"	17" x 9 ½"
Wheels	5"	5"	5"	5"	8"	8"
Product weight	13 ½ lb.	13 ½ lb.	13 ½ lb.	13 ½ lb.	13 ½ lb.	16 ½ lb.
Weight capacity	300 lb.	300 lb.	300 lb.	300 lb.	300 lb.	300 lb.
Units/carton	1	1	1	1	1	1
Limited warranty						
frame	lifetime	lifetime	lifetime	lifetime	lifetime	lifetime
brakes	6 years	6 years	6 years	6 years	6 years	6 years
other components	1 year	1 year	1 year	1 year	1 year	1 year

*Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

Invacare offers a stylish and unique line of canes. With a wide assortment of handgrip styles, you can fit virtually any of your clients with ease. Our quad canes come in medium and large base varieties, have a 250-pound weight capacity, and allow for proper fit with right- or left-hand use.



Invacare® Standard Crook Canes

Invacare Grip
Model no. 8916 – Black



Invacare® Offset Canes

Invacare Grip with Strap
Model no. 40918-6 – Black
Model no. 8917 – Silver



1 Medium Base Quad Canes

Model no. 3920-2 - Invacare Grip

2 Large Base Quad Canes

Model no. 7827 - Invacare Grip

3 Invacare® Orthopedic Grip Canes

Model no. 8929 – Black (Right Hand)

- Our unique orthopedic grip is ergonomically shaped to spread pressure evenly across the palm

Invacare® Canes



Cane Tips

Model nos. 40850-4, 40851-4, 40852-4, 40853-4

- Each pair of cane tips is retail-packaged

Specifications

Model no.	Units/ carton	Fits cane diameter
40850-4	8	1"
40851-4	8	7/8"
40852-4	8	3/4"
40853-4	8	5/8"



The Invacare® Grip

- Our familiar triangular-shaped grip. The flat top helps reduce the stress a user's hand may experience through everyday cane use.



Cane Tip and Collar

- Locking adjustment collar secures cane at desired length. Tip features reinforced metal insert for longer wear.

Cane Specifications

Model no.	Style	Color	Grip type	Adjustment	*Client height (approx.)	Weight	Weight capacity	Units/ carton	Limited warranty
8916	crook	black	Invacare	30" - 39"	4' 6" - 6' 6"	15 oz.	250 lb.	6	lifetime
40918-6	offset	black	Invacare	29" - 38"	4' 4" - 6' 4"	15 oz.	250 lb.	6	lifetime
8917	offset	silver	Invacare	29" - 38"	4' 4" - 6' 4"	15 oz.	250 lb.	6	lifetime

*All client height ranges are approximate. Each individual should always consult with his/her physician or therapist to determine proper adjustment and usage.

Quad Cane Specifications

Model no.	Base	Color	Grip type	Adjustment	Footprint	Weight	Weight capacity	Units/ carton	Limited warranty
3920-2	medium	silver-vein finish	Invacare	29 1/4" - 38 1/4"	10 1/4" x 6"	1 lb. 13 oz.	250 lb.	2	lifetime
7827	large	silver-vein finish	Invacare	29 1/4" - 38 1/2"	13" x 7 3/4"	2 lb. 10 oz.	250 lb.	2	lifetime

Orthopedic Grip Cane Specifications

Our unique orthopedic grip is ergonomically shaped to spread pressure evenly across the palm.

Model no.	Handle style	Color	Adjustment	height* (approx.)	Weight	Weight capacity	Units/ carton	Limited warranty
8929	orthopedic (right)	black	31" - 40"	4' 6" - 6' 6"	15 oz.	250 lb.	6	3 years

*Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.