

KENDALL™ Neonatal Tape



KENDALL™ neonatal tape is designed with our special clear hydrogel formula to gently yet firmly secure catheters, I.V. lines, lead wires and other medical devices to the fragile, sensitive skin of infants. This tape may be repositioned and reapplied to minimize waste without causing skin irritation or damage. The clear hydrogel also allows the patient's skin condition to be monitored without tape removal. This product is completely latex and DEHP free.

Order Code	Description	Ship Case
MI00677	KENDALL™ Hydrogel Baby Tape Strips	500

Features and Benefits

Features	Benefits
Clear adhesive hydrogel	Repositionable and Re applicable

Frequently Asked Questions

Question	Answer
Is the KENDALL™ Neonatal tape latex and DEHP free?	Yes.
How is the tape packaged?	KENDALL™ Neonatal tape is packaged as individual strips.
What is the shelf life?	KENDALL™ Neonatal tape has a shelf life of 24 months.
What size is the tape?	KENDALL™ Neonatal tape come in 3/4" x 2" strips.