

KODAK DIRECTVIEW CLASSIC V SYSTEM

Kodak – a Canadian market leader in the computed radiology world has announced a CR unit for the Veterinary market. Using the same high quality components and technology, the **Kodak Classic V** brings value and image quality together in one very economical package.

The color coded user interface makes the operation easy using intuitive tools and functionality to provide enhanced workflow and quick image processing and review. In less than 60 seconds, the high quality image is on your screen ready for manipulation, review or electronic transmission.

This compact, fully DICOM CR package includes:

- 4 computed radiology cassettes – 2 x 14 x 17 and 2 x 24 x 30
- Directview software including grid suppression and low exposure optimization software
- Complete installation and warranty services
- Full in house training session
- Remote diagnostics



Shown with wall mount