

## **Konica Minolta ImagePilot LE Medical SERIES CR System** **The Complete Digital Integrated Radiography Solution**

The ImagePilot LE Solution is ideal for smaller clinics with modest imaging volume. It includes Konica's new generation CR with AutoPilot Image Processing which dramatically simplifies CR image acquisition to a simple push of a button – perfect, consistent images every time! This integrated system offers patient registration, image acquisition, clinical review and long term image storage all in one easy to use and maintain system. It's as easy as 1-2-3!

### **Features:**

- ◆ Patient registration
- ◆ CR Acquisition with AutoPilot image processing
- ◆ Image review with annotation and measurement tools
- ◆ Image manipulation tools including window level settings
- ◆ Export images – DICOM or JPG to disk or CD/DVD with DICOM viewer( PDI)
- ◆ Windows paper printing
- ◆ Compare with prior exams
- ◆ Advanced search for reviewing past exams
- ◆ Long term storage capability and ability to expand
- ◆ One DICOM Store license – i.e. remote PACS or off site storage
- ◆ One Client licence for distributed workflow - i.e. exam room review
- ◆ Simple Note Tool
- ◆ ImagePilot Worklist Toolkit – includes DICOM MML and Study mapping tool
- ◆ DVD back up tool
- ◆ Grid line suppression

### **Regius Model 110 CR Reader**

- ◆ Up to 60 plates per hour (14 x 17)
- ◆ User operated optical unit cleaner for cleaning CR optics and maintaining image quality

### **Image Pilot Imaging Station – Small Form Factor**

- ◆ Lenovo PC – desktop, Intel Core 2 Duo CPU, 2GB Ram, Window XP Pro
- ◆ 250GB HD – holds approximately 10,000 14 x 17 images
- ◆ Choose from 19" Lenovo, 21" or 24" Eizo Color Monitor

### **Cassettes included:**

- ◆ Two 14 x 17 Cassettes and Plates
- ◆ Two 10 x 12 Cassettes and Plates



**Not exactly as shown**