

STAINLESS STEEL EQUIPMENT STANDS

A mobile stand moves your equipment easily wherever it is needed.

COMMON FEATURES: STAINLESS STEEL EQUIPMENT STANDS

- ✓ Versatile mobile equipment stands come in six different variations for your best fit
- ✓ Models with drawers offer enclosed storage. Drawer is 15 1/2" W x 18"D x 4 1/2"H.
- ✓ Stainless steel is durable for carrying heavy equipment and easy to clean and sanitize
- ✓ Optional locks available for one and 2 drawer models
- ✓ Optional braking casters (2)



#356

#359

#474

FREIGHT CLASS
125

STANDARD PRODUCTS
SHIP IN **3** DAYS OR LESS

MODEL	WEIGHT CAPACITY	# OF DRAWERS	# OF SHELVES	SHELF CLEARANCE	DIMENSIONS			CASTER TYPE AND SIZE	SHIP WEIGHT (LBS.)	FedEx UPS
					W	D	H			
359	300	0	3	10 3/4"	16 3/4"	18 3/4"	30 1/8"	all swivel - 3 1/2" dia.	25	OK
469	500	0	3	10 3/4"	16 3/4"	18 3/4"	31"	all swivel - 4" dia.	33	OK
356	300	1	3	8 1/8"	16 3/4"	18 3/4"	30 1/8"	all swivel - 3 1/2" dia.	36	OK
466	500	1	3	8 1/8"	16 3/4"	18 3/4"	31"	all swivel - 4" dia.	50	OK
358	300	2	2	17 1/8"	16 3/4"	18 3/4"	35"	all swivel - 3 1/2" dia.	44	OK
474	500	2	3	10"	20 1/8"	36 3/8"	35"	all swivel - 4" dia.	76	NA

STAINLESS STEEL INSTRUMENT TABLES

These extremely durable tables will serve you well for years to come.

- 18-gauge stainless steel top (and shelf if provided) is easy to clean and sanitize
- 1" square 16-gauge stainless steel legs for superior strength
- 18-gauge reinforcements supply extra security for larger models
- All welded construction increases durability and makes surface smoother for cleaning
- 3" locking casters ensure table is stationary when needed
- Nesting tables save space

STANDARD PRODUCTS
SHIP IN **7** DAYS OR LESS

SPECIFICATIONS

- Top shelves turned down 1" and hemmed back 3/4" around bottom
- Distance between top shelf and bottom shelf is 16"
- Distance from H-Brace to floor is 9 1/8"



#8353



#8383

MODEL	DESCRIPTION	DIMENSIONS			SHIP WEIGHT (LBS.)
		W	L	H	
8353	Instrument Table w/ Shelf	16"	20"	34"	40
8354	Instrument Table w/ Shelf	16"	30"	34"	44
8355	Instrument Table w/ Shelf	18"	33"	34"	52
8356	Instrument Table w/ Shelf	20"	36"	34"	54
8357	Instrument Table w/ Shelf	24"	36"	34"	56
8358	Instrument Table w/ Shelf	20"	48"	34"	66
8359	Instrument Table w/ Shelf	24"	48"	34"	70
8398	Instrument Table w/ Shelf	24"	60"	34"	82
8399	Instrument Table w/ Shelf	24"	72"	34"	103
8386	H-Brace Instrument Table	16"	20"	34"	45
8387	H-Brace Instrument Table	16"	30"	34"	35
8388	H-Brace Instrument Table	18"	33"	34"	40
8390	H-Brace Instrument Table	20"	36"	34"	43
8391	H-Brace Instrument Table	24"	36"	34"	43
8392	H-Brace Instrument Table	20"	48"	34"	47
8393	H-Brace Instrument Table	24"	48"	34"	56
8394	H-Brace Instrument Table	24"	60"	34"	64
8395	H-Brace Instrument Table	24"	72"	34"	71
8380	Nesting Table	14"	28"	32"	38
8381	Nesting Table	16"	32"	34"	38
8382	Nesting Table	18"	36"	36"	45
8383	Nesting Table	20"	40"	38"	49
8384	Nesting Table	22"	44"	40"	52
8385	Nesting Table	24"	48"	42"	38
8389	Set - All 6				260