

## CURAD® Aloetouch latex



- Inside coated with pure organic aloe to moisturize skin
- Textured throughout
- 6.7 mil finger thickness; 6.3 mil palm thickness

## CURAD® Latex



- Textured throughout
- 6.7 mil finger thickness; 5.9 mil palm thickness

## MediGuard® Latex



- Textured throughout
- 5.9 mil finger thickness; 5.1 mil palm thickness

## MediGuard® Powdered Latex



- Powdered with medical-grade absorbable dusting powder (USP) for added comfort and easy donning
- Smooth — not textured
- 5.5 mil finger thickness; 4.3 mil palm thickness

## Eudermic MP™ Heavy-Duty Latex



- Extra thick latex glove for heavy-duty protection
- 12" extended cuff
- Textured fingertips
- 15.0 mil finger thickness; 13.0 mil palm thickness