



Toll Free: 866-529-3526 Fax Free: 888-433-5174

6248 Bittersweet Lane Nisswa, MN 56468

## MRI Step Stools

### Specifications:

16" x 12" Base  
8 1/2" Step Height  
Made of Heavy Duty Stainless Steel  
**500 lb Weight Capacity**

Model	Description	Price
SS-1000	▲ Without 41" Handrail	\$177.45 ea.
SS-1001	▲ With 41" Handrail	\$219.90 ea.

Pricing in this catalog is not current

## MRI Step Stool with Dual 41" Handrail

### Specifications:

30" x 12" Base  
8 1/2" Step Height  
Made of Heavy Duty Stainless Steel  
**500 lb Weight Capacity**

Model	Description	Price
SS-1005	▲ With Dual 41" Handrail	\$374.00 ea.

## MRI Double Step Stool With 41" Handrail

### Specifications:

36" L X 12" W X 8" H Top Step  
36" L X 24" W X 8" H Base  
Weight: 47 lbs  
**400 lb Weight Capacity**

Model	Description	Price
SS-1002	▲ With Dual 41" Handrail	\$715.00 ea.



More options and sizes of Step Stools available online at [www.MRIequip.com](http://www.MRIequip.com), including PVC options with casters and higher weight capacities.

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current



## MRI Stools

Stainless steel stools with solid non-slip material mounted into recess of stool seat. Seat diameter is 14". Comes with rubber tips or casters. Choose from a fixed height stackable stool or heights that adjust with threaded stem, from 15" to 21" or 21" to 27".



CH-8058



CH-8052



CH-8056



CH-8057

Model	Description	Price
CH-8058	▲ Fixed Height Stool w/ Rubber Tips, 18"	\$276.75 ea.
CH-8050	▲ Adjustable Stool w/ Rubber Tips, 15"-21"	\$405.00 ea.
CH-8054	▲ Adjustable Stool w/ 2" Casters, 15"-21"	\$433.00 ea.
CH-8051	▲ Adjustable Stool w/ Rubber Tips, 15"-21" & Back	\$545.00 ea.
CH-8055	▲ Adjustable Stool w/ 2" Casters, 15"-21" & Back	\$573.00 ea.
CH-8052	▲ Adjustable Stool w/ Rubber Tips, 21"-27"	\$410.70 ea.
CH-8056	▲ Adjustable Stool w/ 2" Casters, 21"-27"	\$459.26 ea.
CH-8053	▲ Adjustable Stool w/ Rubber Tips, 21"-27" & Back	\$555.00 ea.
CH-8057	▲ Adjustable Stool w/ 2" Casters, 21"-27" & Back	\$580.00 ea.

Pricing in this catalog is not current

## Storage and Transport Cases

Holes on each side of latch opening allow for padlocks or other securing devices. Top storage area for reaching smaller items without having to open the case. Overall Outside Box Dimensions: 15 1/4"W x 10 1/8"H x 7 5/8"D. Removable top tray. Tray Outside Dimensions: 13 5/8"W x 3"H x 7 1/8"D. Large carry handle. Special rubber sealing gasket keeps contents dry and clean.

Model	Description	Price
MT-1609	Orange	\$36.81 ea.



## Non-Magnetic Safety Padlock

Plastic shackle helps prevent electrocution and arc flashing, also prevents electrical current from shackle to key. All components are non-sparking and non-magnetic - plastic, copper, brass. 1-1/2" wide 1-3/4" tall Xenoy green or red body, plastic shackle w/1-1/2" clearance. Designed exclusively for Lockout/Tagout applications. Durable, lightweight, non-conductive Xenoy lock body. Compliance with OSHA "one employee, one lock, one key" directive. Key retaining-ensures that padlock is not left unlocked. Keyed alike 6-pin tumbler cylinder

Model	Description	Price
PL-1000	Green	\$18.95 ea.
PL-1001	Red	\$18.95 ea.





### Stools

4" Padded Seat

Burgundy Naugahyde Upholstery

Stainless Steel Construction

5 Leg Design for Maximum Stability

Padded Armrests (on selected models)

2" Dual Wheel Casters (on selected models)

Pricing in this catalog is not current



### Adjusts 15" to 21" with Rubber Tips

Model	Description	Price
CH-8030 ▲	Adjustable Stool w/ Rubber Tips	\$388.00 ea.
CH-8031 ▲	Adjustable Stool w/ Rubber Tips & Back	\$458.00 ea.
CH-8032 ▲	Adjustable Stool w/ Rubber Tips, Back & Arms	\$498.00 ea.

### Adjusts 21" to 27" with Rubber Tips

Model	Description	Price
CH-8033 ▲	Adjustable Stool w/ Rubber Tips	\$408.00 ea.
CH-8034 ▲	Adjustable Stool w/ Rubber Tips & Back	\$478.00 ea.
CH-8035 ▲	Adjustable Stool w/ Rubber Tips, Back & Arms	\$518.00 ea.



### Adjusts 16" to 22" with Casters

Model	Description	Price
CH-8036 ▲	Adjustable Stool w/ Casters	\$438.00 ea.
CH-8037 ▲	Adjustable Stool w/ Casters & Back	\$508.00 ea.
CH-8038 ▲	Adjustable Stool w/ Casters, Back & Arms	\$549.00 ea.

### Adjusts 22" to 28" with Casters

Model	Description	Price
CH-8039 ▲	Adjustable Stool w/ Casters	\$435.00 ea.
CH-8040 ▲	Adjustable Stool w/ Casters & Back	\$513.00 ea.
CH-8041 ▲	Adjustable Stool w/ Casters, Back & Arms	\$575.00 ea.



### Custom Stenciling Option

Customize your doctor stool. Add your facility name, department name or any information you would like.

Model	Description	Price
CH-8914	Custom Stenciling Option	\$75.00 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current