



Fixed Gurney

- Stainless Steel construction
- Drop down side rails
- 5" Total Lock Casters
- Grounding chain
- Storage tray
- Full surround rubber bumper
- Adjustable headrest
- 3" pad included

Specifications:

- Height 29 1/2"
- Height with Rails up 38 1/2"
- Height of Rails 9"
- Length 77"
- Rail Length 57 1/2"
- Width: 30"
- Storage Tray 22" x 31"
- Weight: 100 lbs.
- Back Adj. Range 0-75 Degrees
- Number of Positions: 11
- Weight Capacity: 600 lbs**

Warranty

- 2 years on frame
- 90 days on soft goods and moving parts



**Optional Oxygen Holder Available
See page 7**

Model	Description	Price
ST-2000	▲ Fixed Gurney	\$3,219.00 ea.

Pricing in this catalog is not current



MRI Stretcher/Gurney Accessories



Contoured Shifter

Model	Description	Price
PV-3000	■ Contoured Shifter	\$229.00 ea.

Specifications:

Width 22"
Length 72"
Thick 3/16"
See page 11 for more details

I.V. Pole

Model	Description	Price
ST-2006	▲ I.V. Pole	\$175.00 ea.

Fits the ST-2000 Only
Can be attached at any of the four corners
Multiple I.V. Poles can be attached to a single gurney

Side Rail Pads

Model	Description	Price
ST-2007	■ Side Rail Pads	\$224.00 pr.

Fits the ST-2000 Only

Specifications:

Width 8½"
Length 58"

Safety Straps

Model	Description	Price
ST-2008	■ Safety Straps	\$82.50 set

Fits the ST-2000 Only
Set of 3 Straps

Replacement Pad

Model	Description	Price
ST-2005	■ Replacement Pad	\$240.00 ea.

Specifications:

Height 3"
Width 24½"
Length 73½"

Oxygen Holder

Model	Description	Price
ST-2009	▲ Oxygen Holder	\$179.50 ea.

Fits the ST-2000 Only





Ferno MRI Gurney 400 lb Capacity

Model 30NM The "Standard" and the lightest lift-in cot in the industry. Ferno's Model 30 Multi-Level, Lift-in Cot served as the industry standard for many years and is still used by many ambulance companies and medical facilities. The Ferno Model 30 cot adjusts to six height levels, enabling patient transfer from most nursing home beds or hospital treatment tables. To operate the height adjustment mechanism, the EMTs must assume the patients weight prior to raising or lowering the cot-bed frame. Because the Model 30 must be lifted to be loaded into the vehicle, it can be used with ambulances that have various floor heights. The Model 30 cot will work well in station wagons, vans or modular ambulances. The telescoping head-end main frame section allows the Model 30 cot to be shortened for maneuverability around tight corners and spaces. The Model 30NM is the Non-Magnetic version for medical treatment areas where the cot must be non-magnetized.



Contoured Style Pad

Specifications:

Max Height: 32"
Min Height: 12"
Length: 78"
Width: 23"
Weight: 56 lbs
Number of Positions: 6
Casters: All Locking
Weight Capacity: 400 lbs

Warranty

2 Years on frame Pricing in this catalog
90 Days on pad and straps is not current

Model	Description	Price
ST-1000 ▲	With Contoured Pad	\$3,999.00 ea.
ST-1001 ▲	With Flat Style Pad	\$3,999.00 ea.

Ferno MRI Gurney 500 lb Capacity

Model 35ANM The lightest cot available with a 500 lb weight capacity. The Ferno Model 35ANM Mobile Transporter Multi-Level Roll-In Cot combines the benefits of a multi-level cot with the easy loading feature of a roll-in style cot. Designed for a variety of ambulance floor heights, the Mobile Transporter will work well with station wagons, vans or modular ambulances. EMTs no longer have to lift from the side to load the cot. One EMT is positioned at the foot of the cot, while the other EMT simply raises the undercarriage during loading. Eight height positions enable easy patient transfer from any couch, bed or treatment table. The Mobile Transporter should be rolled in a level position. The two inclined positions are designed as loading positions for either vans or modular ambulances. The Model 35ANM is the Non-Magnetic version for medical treatment areas where the cot must be non-magnetized.



Contoured Style Pad

Specifications:

Max Height: 40"
Min Height: 9"
Max Length: 63" -79"
Width: 24"
Weight: 74 lbs
Number of Positions: 8
Casters: All Locking
Weight Capacity: 500 lbs

Warranty

2 Years on frame
90 Days on pad and straps

Model	Description	Price
ST-1015 ▲	With Contoured Pad	\$3,950.00 ea.
ST-1017 ▲	With Flat Style Pad	\$3,950.00 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current



Toll Free: 866-529-3526 Fax Free: 888-433-5174
6248 Bittersweet Lane Nisswa, MN 56468

MRlequip
Non-Magnetic Equipment .com

Ferno Gurney Parts & Accessories



I.V. Pole
For Ferno Gurneys Only

Warranty
1 Year

Contoured Shifter

Specifications:
Width 17 3/4"
Length 72"
Thick 3/16"
See page 11 for details



Model	Description	Price
ST-1004	▲ IV Pole	\$354.60 ea.

Model	Description	Price
PV-3001	■ Contoured Shifter	\$229.00 ea.

Replacement Safety Belts

Set of 3



Replacement Casters



Replacement Pads



Model	Description	Price
ST-1009	■ Safety Belts	\$294.00 ea.
ST-1007	▲ Locking Caster	\$212.00 ea.
ST-1008	▲ Non-Locking Caster 400lb Capacity Gurney	\$159.00 ea.
ST-1020	▲ Non-Locking Caster 500lb Capacity Gurney	\$274.00 ea.

Model	Description	Price
ST-1005	▲ Contoured Pad	\$246.30 ea.
ST-1006	▲ Flat Style Pad	\$222.50 ea.

Replacement Stretcher Pads

Replacement stretcher pads covered with Sure-Check Fusion III material. Available with square corners or angle cut corners. Stretcher pads are **Non-Returnable**. Custom sizes available



Square Corner Pad



Cut Corners Pad

Model	Description	Price
ST-4000	■ 24"W x 75"L x 2" T	\$199.70 ea.
ST-4001	■ 24"W x 75"L x 3" T	\$234.00 ea.
ST-4002	■ 24"W x 75"L x 4" T	\$244.00 ea.
ST-4003	■ 26"W x 75"L x 2" T	\$199.70 ea.
ST-4004	■ 26"W x 75"L x 3" T	\$234.00 ea.
ST-4005	■ 26"W x 75"L x 4" T	\$244.00 ea.
ST-4006	■ 28"W x 75"L x 2" T	\$199.70 ea.
ST-4007	■ 28"W x 75"L x 3" T	\$234.00 ea.
ST-4008	■ 28"W x 75"L x 4" T	\$244.00 ea.

Model	Description	Price
ST-4009	■ 24"W x 75"L x 2" T	\$199.70 ea.
ST-4010	■ 24"W x 75"L x 3" T	\$234.00 ea.
ST-4011	■ 24"W x 75"L x 4" T	\$244.00 ea.
ST-4012	■ 26"W x 75"L x 2" T	\$199.70 ea.
ST-4013	■ 26"W x 75"L x 3" T	\$234.00 ea.
ST-4014	■ 26"W x 75"L x 4" T	\$244.00 ea.
ST-4015	■ 28"W x 75"L x 2" T	\$199.70 ea.
ST-4016	■ 28"W x 75"L x 3" T	\$234.00 ea.
ST-4017	■ 28"W x 75"L x 4" T	\$244.00 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current