



Toll Free: 866-529-3526 Fax Free: 888-433-5174
6248 Bittersweet Lane Nisswa, MN 56468

MRlequip
Non-Magnetic Equipment .com

MRI Folding Walkers

Pricing in this catalog is not current

18" Wide, Depth 17"
250 lb weight capacity

Model	Description	Price
WA-1000 ▲	Adult - 32" to 36"	\$195.00 ea.
WA-1001 ▲	Junior - 28" to 32"	\$195.00 ea.

MRI Heavy Duty / Extra Wide Folding Walker

24" Wide, Depth 21"
500 lb weight capacity

Model	Description	Price
WA-1002 ▲	Extra Wide - 32" to 38½"	\$227.50 ea.

5" Adjustable walker Wheels

Converts folding walker into wheeled walker.
Adjustable for 8 height adjustments

350 lb weight capacity

Model	Description	Price
WA-1010 ▲	5" Walker Wheels	\$49.75 pr.

MRI Transfer / Gait Belt

2" x 48"

Model	Description	Price
GB-1000 ■	2" x 48" Transfer/Gait Belt	\$18.00 ea.

20" Wide Transport Wheelchair

Removable Full Length Padded Arms, Swingaway, Detachable Footrests or Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 20"
Overall Width: 25"
Overall Depth: 27" w/o footrest
Seat Depth: 16"
Overall Height: 35.75"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 12.50" x 2.50"
Wheel Locks: pull to lock
Weight: 45lbs

Weight Capacity: 250 lbs

Warranty

5 years on frame and crossbraces

6 months on soft goods and moving parts

Model	Description	Price
WC-1027 ▲	Chair With Footrests	\$1,250.00 ea.
WC-1028 ▲	Chair With Legrests	\$1,375.00 ea.

See page 5 for wheelchair parts and accessories

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current



NEW



Back View

Chair With Footrests



Chair With Legrests



Chair With Footrests



Back View

Chair With Legrests

24" Wide Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 24"
Overall Width: 32"
Overall Depth: 33" w/o footrest
Seat Depth: 18"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 46 lbs w/o footrest

Weight Capacity: 350 lbs

Pricing in this catalog

Warranty

is not current

5 years on frame and crossbraces

6 months on soft goods and moving parts

Model	Description	Price
WC-1024	▲ Chair With Footrests	\$2,148.00 ea.
WC-1025	▲ Chair With Legrests	\$2,348.00 ea.

See page 5 for wheelchair parts and accessories



Chair With Footrests



Back View

Chair With Legrests

26" Wide Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 26"
Overall Width: 35"
Overall Depth: 36" w/o footrest
Seat Depth: 20"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 61 lbs w/o footrest

Weight Capacity: 550 lbs

Warranty

5 years on frame and crossbraces

6 months on soft goods and moving parts

Model	Description	Price
WC-1003	▲ Chair With Footrests	\$2,595.00 ea.
WC-1026	▲ Chair With Legrests	\$2,795.00 ea.

See page 5 for wheelchair parts and accessories

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current



20" Wide Wheelchair, With Full Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 20"
Overall Width: 28"
Overall Depth: 33" w/o footrest
Seat Depth: 16"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 46 lbs w/o footrest

Weight Capacity: 250 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts

Model	Description	Price
WC-1004	▲ Chair With Footrests	\$1,575.00 ea.
WC-1022	▲ Chair With Legrests	\$1,670.00 ea.

Chair With
Footrests



Back View

Chair With
Legrests



See page 5 for wheelchair parts and accessories

Pricing in this catalog is not current

22" Wide Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 22"
Overall Width: 30"
Overall Depth: 33" w/o footrest
Seat Depth: 18"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 46 lbs w/o footrest

Weight Capacity: 350 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts

Model	Description	Price
WC-1002	▲ Chair With Footrests	\$1,998.00 ea.
WC-1023	▲ Chair With Legrests	\$2,198.00 ea.

Chair With
Footrests



Back View

Chair With
Legrests



See page 5 for wheelchair parts and accessories

Pricing in this catalog is not current

■ MR Safe ▲ MR Conditional ● Not MR Safe



18" Wide Wheelchair, With Fixed Footrest

Removable Full Length Padded Arms
Fixed Footrests
Color: Burgundy

Specifications

Seat Width: 18"
Overall Width: 26"
Overall Depth: 44" w/footrest
Seat Depth: 16"
Overall Height: 35"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 44 lbs w/footrest

Weight Capacity: 250 lbs Pricing in this catalog is not current

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Back View

Model	Description	Price
WC-1000	▲ Chair With Footrests	\$1,398.00 ea.

See page 5 for wheelchair parts and accessories



Chair With Footrests

18" Wide Wheelchair, With Flip-Up Arms

Flip-Up Desk Length Padded Arms
Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests
Color: Burgundy

Specifications

Seat Width: 18"
Overall Width: 26"
Overall Depth: 31" w/o footrest
Seat Depth: 16"
Overall Height: 35"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 40 lbs w/o footrest

Weight Capacity: 250 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Back View

Model	Description	Price
WC-1001	▲ Chair With Footrests	\$1,570.00 ea.
WC-1021	▲ Chair With Legrests	\$1,670.00 ea.

See page 5 for wheelchair parts and accessories



■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current

Chair With Legrests



MRI Wheelchair Accessories



H Custom Stenciling, add your company name and/or other info on back of wheelchair

Fig	Model	Description	Price
A	RP-1053	▲ Anti-Tippers	\$44.00 pr.
B	RP-1085	▲ Elevating Legrest for 18" wide flip back arms, 20" wide, 22" wide & 24" wide wheelchair	\$225.00 pr.
B	RP-1086	▲ Elevating Legrest for 26" wide wheelchair	\$300.00 pr.
C	WC-1102	▲ Oxygen Holder	\$90.45 ea.
D	WC-1100	▲ I.V. Pole	\$101.85 ea.
E	WC-1101	▲ I.V. & Oxygen Holder	\$163.00 ea.
F	RP-1054	▲ Anti-Theft Bar for 18" Wheelchair	\$101.85 ea.
F	RP-1077	▲ Anti-Theft Bar for 20" Wheelchair	\$138.50 ea.
F	RP-1055	▲ Anti-Theft Bar for 22" Wheelchair	\$138.50 ea.
F	RP-1078	▲ Anti-Theft Bar for 24" Wheelchair	\$138.50 ea.
G	RP-1114	▲ Anti-Theft Clamps for 26" Wheelchair, set of 4	\$74.50 set
H	WC-1110	Custom Stenciling (Company Name, etc.)	\$175.00 ea.

Pricing in this catalog is not current

MRI Wheelchair Parts - MR Conditional

Model	Description	Cost	Model	Description	Cost
RP-1000	Footrest for 18" wide flip back arm and 20" wide chair, left	\$69.25 ea.	RP-1018	Bearing for rear wheel - 7/16" Axle	\$8.95 ea.
RP-1001	Footrest for 18" wide flip back arm and 20" wide chair, right	\$69.25 ea.	RP-1019	Bearing for rear wheel - 5/8" Axle	\$8.95 ea.
RP-1002	Footrest for 22" & 24" wide chair, Left	\$119.10 ea.	RP-1025	Brake for all MRI chairs, left	\$26.75 ea.
RP-1003	Footrest for 22" & 24" wide chair, right	\$119.10 ea.	RP-1026	Brake for all MRI chairs, right	\$26.75 ea.
RP-1004	Footrest for 18" wide fixed chair, left	\$63.50 ea.	RP-1028	Seat uph, for 18" wide chair, Burgundy	\$64.50 ea.
RP-1005	Footrest for 18" wide fixed chair, right	\$63.50 ea.	RP-1029	Seat uph, for 22" wide chair, Burgundy	\$85.00 ea.
RP-1006	Front roller bumper for footrest	\$9.75 ea.	RP-1030	Back uph, for 18" std chair, Burgundy	\$75.00 ea.
RP-1007	Heel loop and hardware for footrest	\$13.85 ea.	RP-1031	Back uph, for 18" flip back arm, Burgundy	\$77.40 ea.
RP-1009	Front fork assy. fork, wheel, bearings, Hardware	\$65.45 ea.	RP-1032	Back uph, for 22" wide chair, Burgundy	\$77.40 ea.
RP-1010	8" Wheel assy. with bearings	\$18.25 ea.	RP-1033	Seat uph, for 18" wide chair, Dark Blue	\$64.50 ea.
RP-1011	Bearing for front caster wheel 5/16"	\$7.95 ea.	RP-1034	Back uph, for 18" std chair, Dark Blue	\$75.00 ea.
RP-1012	8" front fork with nut	\$41.75 ea.	RP-1035	Padded armrest, full length, Burgundy	\$12.40 ea.
RP-1013	Fork stem bearing	\$7.95 ea.	RP-1036	Padded armrest, full length, Dk Blue	\$12.40 ea.
RP-1016	24" complete rear wheel, with 7/16" bearings	\$155.00 ea.	RP-1037	Plastic armrest, Black for flip back arm chair	\$7.95 ea.
RP-1017	24" complete rear wheel, with 5/8" bearings	\$155.00 ea.	RP-1038	Padded armrest, desk length, Burgundy	\$12.40 ea.
			RP-1046	Button lock for removable arms 22" wide	\$16.25 ea.
			RP-1115	Locking kits for foot/legrest	\$11.45 pr.
			RP-1117	Vice Tire Setter for Snap on Tires	\$27.50 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe

Pricing in this catalog is not current