



- Quality X-Ray Grids
- Radiology Accessories
- Cassette Repair
- Barium Swallow Chair

- Grid Protection Devices
- Mammography Chairs
- Ultrasound Table
- Ultrasound Accessories



DELUXE HYDRAULIC MAMMOGRAPHY CHAIR

MXE also provides two hydraulic positioning chairs designed with the patient and the technologist in mind. To assist with stressed patients, the Deluxe Reclining Mammo Chair back reclines to a maximum of 170 degrees. The five star base with locking castors provides additional stability.

Click on
Thumbnails to
See Detailed
Pictures!



Criteria	Deluxe Reclining	Deluxe Non-reclining
Seat Height	23" to 32"	23" to 32"
Mechanism	Hydraulic lift	Hydraulic lift
Arm Rests	Removable	Removable
Head Rest	Adjustable removable	N/A
Foot Rest	Removable	Removable
Chair Back	Reclines to 170 degrees	N/A
Chair Base	Five star, porcelain covered cast iron	Five star, porcelain covered cast iron
Castors	Locking, five	Locking, five
Foot Control	Raises/lowers patients, locks seat rotation	Raises/lowers patients, locks seat rotation
Color	Light gray	Light gray
Finish	High density powder coated	High density powder coated
Seat Cover	Fire retardant, stain resistant vinyl	Fire retardant, stain resistant vinyl
Shipping Weight	Approximately 95 lbs.	Approximately 95 lbs.
Weight Limit	300 lbs. maximum	300 lbs. maximum