



MXE's Ultrasound Table is the perfect table for ultrasound in ob/gyn and vascular procedures. The MXE Ultrasound Table features three adjustable sections for patient comfort: sitting, partial reclining, and supine. For added comfort, the patient's feet can be elevated. The MXE Ultrasound Table has a motorized hydraulic lift that is operated by a simple foot pedal providing effortless height adjustment. The Table adjusts smoothly from 27" to 37" for patient positioning and proper scanning height.

The MXE Ultrasound Table's sturdy, all steel frame allows up to 375 lbs. of patient weight and a 3" foam pad assures patient comfort even during lengthy procedures. For easy cleaning, the Table's foam pad and removable arms are covered with washable vinyl. Four locking, swivel castors provide maximum mobility plus stability and security.



Features:

- **Adjustable to sitting position with simple hand operated locking lever**
- **Sturdy all steel frame**
- **Four locking 3" castors**
- **Paper roll holder**
- **3" foam pad**
- **Durable vinyl cover**
- **Color: white only**
- **Safety strap**
- **Weight limit up to 375 lbs.**
- **Motorized hydraulic lift with foot pedal**
- **Padded arm rests, adjustable and removable**
- **Optional stirrups and/or IV Pole**



Dimensions:

- 28" wide x 6' long
- 33" wide at arm rests
- 26" to 37" adjustable height (w/castors)

Shipping Weight: 138 lbs.

