

mediport

The Maxant **MediPort 3000** is an *All-In-One*, integrated **Monitor** and **CPU** workstation in a standard x-ray viewer footprint. Designed for space and workflow efficiency.

All **MediPort 3000** workstations are Medical Grade and safe to be installed in any patient area care or surgical environment. It meets the stringent UL 60601 Hospital Standards and has 501k Clearance from the FDA as a Medical Device.

Each **MediPort 3000** is equipped with a high performance CPU featuring the following specifications:

- Intel Core 2 Duo Processor
- Windows 7 Professional
- Up to 4 GB of RAM
- 80 GB Hard Drive
- CD/DVD ROM
- 10/100/1000 LAN Controller
- Keyboard
- Mouse
- Wall Mounting Bracket

MediPort 3000 monitor options include:

- 1.3 MP Color
- 2.0 MP Color
- 3.0 MP Grayscale



mediport
A solution that always fits

CAN-med Healthcare

Unit 1, 130 Matheson Boulevard East, Mississauga, Ontario, Canada L4Z 1Y6
Tel: (905) 502 1166 • Toll Free: (866) 277 3577 • Fax: (905) 502 1167