

# mediport

by MAXANT

Introducing the Newest Generation  
in Medical Display Workstations  
*MediPort 3000*



The solution that always fits.

MediPort delivers convenient access to medical images and data wherever you need it—even sensitive patient care areas. Ideal for PACS, EMR and virtually all other medical software applications, the new MediPort 3000 series is a comprehensive line of all-in-one PC workstations offering medical-grade technology, seamless compatibility and flexible configuration. Now more than ever, MediPort is *the solution that always fits*.

## Why “Medical Grade” Matters

Along with many benefits, digital imaging also presents new risks for many medical facilities. For instance, off-the-shelf PCs are safe for general business use, but they’re tested to standards that are much lower than those of products used in medical environments.

Studies have shown that patients with even minor health issues are many times more susceptible to harm from electrical shock. That’s why products certified to medical-grade electrical safety standards are so important: they virtually eliminate the risk to patients, doctors and other medical personnel.

MediPort, one of the first medical-grade digital imaging workstations, complies with even the most stringent requirements for medical and electrical safety:

- 510k FDA clearance
- UL 60601-01 medical-grade safety requirements
- National Fire Protection Association’s NFPA 99 electrical safety standards
- National Electrical Code

Designed specifically with patient safety in mind, MediPort 3000 is safe for *every* patient care environment within a healthcare facility—including surgery, ED and all other critical care areas.

## Smarter Use of Your Space

With single-, double- and triple-head workstations, MediPort 3000 can be configured for any space where physicians, nurses, technologists and practitioners need convenient access to images and patient data.

It’s easy to install, and it takes no more space than a wall-mounted film viewer.

Simply put, MediPort can be incorporated easily into virtually any layout, any workflow, for any need.



*Optional large display monitor\**

# Great Versatility in a Small Package

When MAXANT set out to develop the MediPort 3000 series, we listened closely to our customers' needs. The result is an expanded line of medical display workstations that are ideal for a wider-than-ever variety of medical environments and applications.

Why choose the 3000 series for your facility?

- High-performance, all-in-one PC and monitor
  - Reduces renovation costs
  - Clean, unobtrusive look
  - Hassle-free installation in your existing space
- Thin-profile, 3.8-inch depth
- Diagnostic or clinical review configurations
- Compliance with DICOM standards
- Medical-grade, safe for patient environments
- Can be used on a mobile cart\*
- Can be used with a large display monitor\*



## Inside the Technology

MediPort 3000 models feature a high-performance, integrated computer and monitor. Each system meets medical-grade specifications for the OR, ED, ICU and other patient care environments. For greater display flexibility, MediPort can be configured with monitors ranging from 1.3 MP color to 3.0 MP gray.

### Computer Specifications and Upgrades

- Intel Core 2 Duo or Dual-Core Processor
- Intel chipset and motherboard
- 1GB RAM (up to 4GB\*)
- 10/100/1000 LAN controller (Ethernet)
- 80GB SATA hard drive (up to 500GB\*)
- CD-R/DVD ROM (DVD-RW\*)
- Windows XP Pro 32 bit

### Display Options

- Choose 1, 2 or 3 monitors
- Monitor 1: 1.3 MP color; or 2 MP color or gray
- Monitors 2 and 3: 1.3 MP color; or 2 MP color or gray; or 3 MP gray

### Standards and Certifications

- 510k FDA clearance
- 60601-1
- NFPA 99
- CAN/CSA C22.2 No. 601.1

(\*indicates optional feature)

DISCLAIMER: Specifications may change. Contact Maxant to verify specifications.



Mobile cart\*


## Choose Your MediPort:



MediPort 3100



MediPort 3200



MediPort 3300

### An ideal solution for:

- ICU, ED and OR
- Patient rooms
- Imaging centers
- Radiology
- Clinical work areas
- Consultation suites
- Physician offices
- Wherever access to images is needed



### MediPort 1100

Shown here with two features available only on the MediPort 1000 series:

- Recess mount (projects 2.75" from wall)
- Optional MAXANT TR401 film viewer

Discover more about the MediPort 1000 series at [www.maxant.com](http://www.maxant.com).

## About MAXANT

MAXANT is a leading manufacturer of high-quality diagnostic imaging devices and accessories. Founded in 1976, we continue to build our business on understanding healthcare professionals' evolving needs and meeting them with innovative, practical products. Visit us today at [www.maxant.com](http://www.maxant.com).

# MAXANT

MAXANT Technologies, Inc.  
7540 Caldwell Avenue  
Niles, IL 60714  
USA  
Tel: 800.307.4190  
Fax: 847.588.1920  
Web: [www.maxant.com](http://www.maxant.com)

**MAXANT**



**mediport**

