

*Customize compression levels
and adjust range of motion*

ORT28100



Ankle Supports



A Neoprene Ankle Support with Open Heel is low profile and easy to wear

- Made of a four-way stretch, 1/8" (0.3 cm) neoprene that gives uniform compression
- Provides warmth and compression
- Open heel and toe construction for comfort
- Fits right or left ankle
- Measurements based upon relative shoe sizing

Description	Item Number	Qty
Small, Shoe size: Men's 5-6, Women's 5-6½	ORT26200S	ea.
Medium, Shoe size: Men's 6-8½, Women's 6½-9	ORT26200M	ea.
Large, Shoe size: Men's 8½-10, Women's 9-11	ORT26200L	ea.
X-Large, Shoe size: Men's 10-12, Women's 11-13	ORT26200XL	ea.
2X-Large, Shoe size: Men's 12-14, Women's 13-14	ORT262002XL	ea.



B Elastic Ankle Support helps control edema by providing even compression

- Made with latex-free elastic that provides support to weak or injured ankles
- Open heel and toe construction for comfort
- Low profile design allows for use inside of a shoe
- Fits right or left ankle
- Measurements based upon ankle circumference

Description	Item Number	Qty
7-8½" (17.8-21.6cm), Small	ORT26100S	ea.
8½-10" (21.6-25.4 cm), Medium	ORT26100M	ea.
10-11½" (25.4-29.2cm), Large	ORT26100L	ea.
11½-13" (29.2-33 cm), X-Large	ORT26100XL	ea.
13-14½" (33-36.8 cm), 2X-Large	ORT261002XL	ea.



C Neoprene Ankle Support with Zipper makes it easy to apply to an injured ankle

- Made of a four-way stretch, 1/8" (0.3 cm) neoprene for uniform compression
- Zipper closure increases ease of application
- Unique soft-seam construction is smooth to reduce irritation caused by conventional methods of seaming
- Fits right or left ankle
- Measurements based upon ankle circumference

Description	Item Number	Qty
7-8" (17.8-20.3 cm), X-Small	ORT26210XS	ea.
8-9" (20.3-22.9 cm), Small	ORT26210S	ea.
9-10" (22.9-25.4 cm), Medium	ORT26210M	ea.
10-11" (25.4-27.9 cm), Large	ORT26210L	ea.
11-12" (20.3-22.9 cm), X-Large	ORT26210XL	ea.
Over 13" (33 cm), 2X-Large	ORT262102XL	ea.



D Customize the amount of compression with Medline's Stirrup Ankle Splint with Air/Foam Bladders

- Combination air and foam bladders can be inflated and deflated for desired level of compression
- Foot strap adjusts to accommodate many widths
- Hook-and-loop straps for secure application
- Fits right or left ankle

Description	Item Number	Qty
10" (25.4 cm) H, Universal	ORT27200	ea.
9½" (24.1 cm) H, Trainer	ORT27200T	ea.
8¾" (22.2 cm) H, Small	ORT27200S	ea.

F Active Ankle® Multiphase Support has a hinge to allow plantar/dorsiflexion

- Anatomically hinged U-Shape design relieves pressure from the ankle joint
- Provides superior inversion/eversion protection
- Anatomical bilateral hinge allows for plantarflexion/dorsiflexion
- Fits right or left ankle
- Measurements based upon relative shoe sizing

Description	Item Number	Qty
X-Small, Shoe Size: Women's 5-6	ORT27400XS	ea.
Small, Shoe Size: Men's 5-8½, Women's 6½-9½	ORT27400S	ea.
Medium, Shoe Size: Men's 9-12, Women's 10-13	ORT27400M	ea.
Large, Shoe Size: Men's 12½-16, Women's 13½-16	ORT27400XL	ea.



E Air/Gel Ankle Stirrup conforms to the ankle while limiting ankle rotation

- Can be used for cold therapy to help treat edema
- Air & gel bladders are separated with the gel bladder on the bottom and the air bladder on the top
- Fits right or left ankle

Description	Item Number	Qty
10" (25.4 cm) H, Universal	ORT27220	ea.
9½" (24.1 cm) H, Trainer	ORT27220T	ea.
8¾" (22.2 cm) H, Small	ORT27220S	ea.



Ankle Supports

A Deluxe Foam Ankle Stirrup provides stability to the ankle by minimizing ankle rotation

- Made of 2 polypropylene outer shells attached with a soft adjustable foot strap
- Memory foam pads contour to the ankle while providing uniform compression and a customized fit
- Fits right or left ankle



Description	Item Number	Qty
10" (25.4 cm) H, Universal, Black	ORT27110	ea.
9½" (24.1 cm) H, Trainer, Black	ORT27110T	ea.
8¾" (22.2 cm) H, Small, Black	ORT27110S	ea.

B



B Canvas Ankle Lace-Up Support has removable stays for just the right amount of support

- Exterior is a high-quality canvas outer with a soft flannel lining for patient comfort
- Removable and moldable medial, lateral, and tongue aluminum stays
- Lace closure allows for proper fit around swollen ankle
- Anatomically designed to fit the contours of the ankle and heel
- Fits right or left ankle
- Measurements based upon ankle circumference

Description	Item Number	Qty
8-9" (20.3-22.9 cm), Small	ORT26300S	ea.
9-11" (22.9-27.9 cm), Medium	ORT26300M	ea.
11-13" (27.9-33 cm), Large	ORT26300L	ea.
Over 13" (33cm), X-Large	ORT26300XL	ea.

C



C Stirrup Ankle Splint with Memory Foam conforms easily to the ankle for a comfortable fit

- Memory foam contours to the swelling of the injured ankle
- Removable fabric covered foam liner can be hand washed
- Foot strap adjusts to accommodate many widths
- Hook-and-loop straps for secure application
- Fits right or left ankle

Description	Item Number	Qty
10" (25.4 cm) Height, Universal	ORT27100	ea.



D Medline's Short Leg Walker is designed to help maintain a normal gait

- Low profile design minimizes the difference in leg height to help maintain a normal walking gait and balance
- Wide distance between the uprights accommodates extra dressings and swelling
- Three leg straps and two foot straps are adjustable for desired size and compression
- Nonskid rocker bottom sole for safety
- Fits left or right leg

Description	Item Number	Qty
X-Small, Shoe size: Women's 4-5½	ORT28100XS	ea.
Small, Shoe size: Men's 5-7½, Women's 6-8½	ORT28100S	ea.
Medium, Shoe size: Men's 8-11, Women's 9-12	ORT28100M	ea.
Large, Shoe size: Men's 12-14, Women's 13	ORT28100L	ea.
X-Large, Shoe size: Men's 14+	ORT28100XL	ea.



E Ankle Walker Offers the same great features as the Short Leg Walker in a shorter height

- Low profile design minimizes the difference in leg height to help maintain a normal walking gait and balance
- Wide distance between the uprights accommodates extra dressings and swelling
- One leg strap and two foot straps are adjustable for desired size and compression
- Nonskid rocker bottom sole for safety
- Fits left or right leg

Description	Item Number	Qty
X-Small, Shoe size: Women's 4-5½	ORT28200XS	ea.
Small, Shoe size: Men's 5-7½, Women's 6-8½	ORT28200S	ea.
Medium, Shoe size: Men's 8-11, Women's 9-12	ORT28200M	ea.
Large, Shoe size: Men's 12-14 Women's 13	ORT28200L	ea.
X-Large, Shoe size: Men's 14+	ORT28200XL	ea.

* Measurements based upon relative shoe sizing



A Medline's Deluxe Short Leg Walker is wider on the sides, conforming to your anatomy to ensure a better fit

- Contoured sides allow the walker frame to conform to the patient's anatomy
- Closed heel design helps protect the heel and prevent migration
- Nonskid molded rocker bottom sole for safety
- Low-profile design minimizes the difference in leg height to help maintain a normal walking gait
- Fits left and right leg

Description	Item Number	Qty
X-Small, Shoe size: Men's 2-4, Women's 3½-5½	ORT28110XS	ea.
Small, Shoe size: Men's 4½-7, Women's 6-8	ORT28110S	ea.
Medium, Shoe size: Men's 7½-10½, Women's 8½-11½	ORT28110M	ea.
Large, Shoe size: Men's 11-12½, Women's 11½-13½	ORT28110L	ea.
X-Large, Shoe size: Men's 12½+, Women's 13½+	ORT28110XL	ea.



B Deluxe Pneumatic Short Leg and Ankle Walkers allow you to pump air into the walker for desired compression

- Pneumatic system increases stability while decreasing pain and swelling
- Air is easily adjusted to achieve desired compression and fit
- Low profile design minimizes the difference in leg height to help maintain a normal walking gait and balance
- Nonskid rocker bottom sole for safety
- Accommodates swelling patterns that occur throughout the rehabilitation process
- Fits left or right leg

Description	Item Number	Qty
Short Leg Walker		
X-Small, Shoe size: Men's 2-4, Women's 3½-5½	ORT28610XS	ea.
Small, Shoe size: Men's 4½-7, Women's 6-8	ORT28610S	ea.
Medium, Shoe size: Men's 7½-10½, Women's 8½-11½	ORT28610M	ea.
Large, Shoe size: Men's 11-12½, Women's 11½-13½	ORT28610L	ea.
X-Large, Shoe size: Men's 12½+, Women's 13½+	ORT28610XL	ea.
Ankle Walker		
X-Small, Shoe size: Men's 2-4, Women's 3½-5½	ORT28620XS	ea.
Small, Shoe size: Men's 4½-7, Women's 6-8	ORT28620S	ea.
Medium, Shoe size: Men's 7½-10½, Women's 8½-11½	ORT28620M	ea.
Large, Shoe size: Men's 11-12½, Women's 11½-13½	ORT28620L	ea.
X-Large, Shoe size: Men's 12½+, Women's 13½+	ORT28620XL	ea.

* Measurements based upon relative shoe sizing



C



C Articulating (ROM) Short Leg and Ankle Walkers allow for range of motion adjustment through recovery

- Range of motion is adjustable from fixed to 45 degrees
- Closed heel design helps protect the heel and prevent migration
- Low profile design minimizes the difference in leg height to help maintain a normal walking gait and balance
- Nonskid rocker bottom sole for safety
- Fits left or right leg

Description	Item Number	Qty
Short Leg Walker		
Small, Shoe size: Men's 4½-6, Women's 5½-7	ORT28800S	ea.
Medium, Shoe size: Men's 6½-10, Women's 7½-11	ORT28800M	ea.
Large, Shoe size: Men's 10½-12, Women's 11½-13	ORT28800L	ea.
Ankle Walker		
Small, Shoe size: Men's 4½-6, Women's 5½-7	ORT28900S	ea.
Medium, Shoe size: Men's 6½-10, Women's 7½-11	ORT28900M	ea.
Large, Shoe size: Men's 10½-12, Women's 11½-13	ORT28900L	ea.

D



D Medline's Diabetic Walker was designed to help with the treatment of plantar ulcers

- Low profile design to allow for comfort and increased stability
- Pressure relief insole and removable pegs help off-load plantar pressures and shear reduction
- Fits left or right leg

Description	Item Number	Qty
Small, Shoe Size: Men's 4½-6, Women's 5½-7	ORT28120S	ea.
Medium, Shoe Size: Men's 6½-10, Women's 7½-11	ORT28120M	ea.
Large, Shoe Size: Men's 10½-12, Women's 11½-13	ORT28120L	ea.

* Measurements based upon relative shoe sizing

Night Splints



A

- A** Dorsal Night Splint designed to hold the foot in a neutral position during sleep for added comfort
- Low profile, less bulk/weight and less heat retention
 - Helps alleviate morning pain and discomfort associated with plantar fasciitis

Description	Item Number	Qty
Small, Shoe Size: Men's 5-10, Women's 5-9½	ORT28400S	ea.
Large, Shoe Size: Men's 10½-16, Women's 10-16	ORT28400L	ea.



B

- B** Adjustable Night Splint provides relief from the discomfort of plantar fasciitis
- Lightweight and flexible one-piece shell
 - Adjustable flexion straps help provide gentle stretching of the plantar fascia and Achilles tendon
 - Adjusts from 10° to 90° dorsi-flexion
 - Includes a cool foam liner and padded calf and ankle straps for added comfort

Description	Item Number	Qty
Small, Shoe Size: Men's up to 7, Women's up to 7½	ORT28500S	ea.
Medium, Shoe Size: Men's 7½-9, Women's 8-9½	ORT28500M	ea.
Large, Shoe Size: Men's 9½-11½, Women's 10-12½	ORT28500L	ea.

* Measurements based upon relative shoe sizing

Cast protection and foot support for any weather

ORT30300

Cast Boots



A Black Canvas Cast Sandal offers a lightweight way to protect your cast

- Made of a quality canvas duck upper for durability
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Open toe and open heel for comfort and ventilation
- Heel and toe pitch promotes a normal gait
- XXS size comes in pediatric print

Description	Item Number	Qty
5½" (14 cm) L, 3" (7.6 cm) W, 2X-Small Shoe size: Child's 5½-7½	ORT29100XXS	ea.
7¼" (18.4 cm) L, 3½" (8.9 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT29100XS	ea.
8¾" (22.2 cm) L, 4½" (11.4 cm) W, Small Shoe size: Women's 5-7	ORT29100S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium Shoe size: Men's 9-10, Women's 10-12	ORT29100M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large Shoe size: Men's 12-15½	ORT29100L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT29100XL	ea.



B Vinyl Closed Toe Cast Boot protects cast from cold, damp weather

- Closed toe and heel is ideal for colder climates
- Vinyl upper resists water and easily wipes clean after exposure to the elements
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Heel and toe pitch promotes a normal gait

Description	Item Number	Qty
5½" (14 cm) L, 3" (7.6 cm) W, 2X-Small Shoe size: Child's 5½-7½	ORT29210XXS	ea.
7¼" (18.4 cm) L, 3⅝" (9.2 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT29210XS	ea.
8¾" (22.2 cm) L, 4¾" (11.1 cm) W, Small Shoe size: Women's 5-7	ORT29210S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium Shoe size: Men's 9-10, Women's 10-12	ORT29210M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large Shoe size: Men's 12-15½	ORT29210L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT29210XL	ea.

C Vinyl Open Toe Cast Boot protects cast while providing some ventilation

- Closed heel design helps keep the cast dry in damp weather
- Open toe for comfort and breathability
- Vinyl upper resists water and easily wipes clean after exposure to the elements
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Heel and toe pitch promotes a normal gait

Description	Item Number	Qty
4½" (11.4 cm) L, 3½" (8.9 cm) W, 2X-Small, Shoe size: Child's 5½-7½	ORT29200XXS	ea.
7¼" (18.4 cm) L, 3⅝" (9.2 cm) W, X-Small, Shoe size: Women's 3-4½, Child's 10½-12½	ORT29200XS	ea.
8¾" (22.2 cm) L, 4¾" (11.1 cm) W, Small, Shoe size: Women's 5-7	ORT29200S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium, Shoe size: Men's 9-10, Women's 10-12	ORT29200M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large, Shoe size: Men's 12-15½	ORT29200L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large, Shoe size: Men's 16+	ORT29200XL	ea.



* Measurements based upon product size dimensions and relative shoe sizing



D Canvas Molded Rocker Cast Sandal offers a secure fit without the bulk

- Rocker sole promotes a more normal gait with heel and toe pitch
- Two canvas straps secure with hook-and-loop over the top of the foot
- Hook-and-loop heel strap for a secure fit



E Rainbow Cast Sandal is perfect for warm climates

- Multicolored elastic with hook-and-loop closures secure around the cast providing a high level of ventilation
- Thick rocker sole is made of a lightweight EVA
- Rocker sole promotes a more normal gait by incorporating both heel and toe pitch

Description	Item Number	Qty
7¼" (18.4 cm) L, 3" (7.62 cm) W, X-Small Shoe size: Child's 10½-12½	ORT290010XS	ea.
9½" (24.1 cm) L, 4¼" (10.8 cm) W, Small Shoe size: Men's 5-7½, Women's 6-8½	ORT290010S	ea.
11¼" (28.6 cm) L, 5" (12.7 cm) W, Medium Shoe size: Men's 10½-12½, Women's 11½-13	ORT290010M	ea.
12¾" (32.4 cm) L, 5½" (14 cm) W, Large Shoe size: Men's 15-16	ORT290010L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT290010XL	ea.

Description	Item Number	Qty
7¼" (18.4 cm) L, 3½" (8.9 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT290040XS	ea.
9½" (24.1 cm) L, 4½" (10.5 cm) W, Small Shoe size: Men's 5-7½, Women's 5-8½	ORT290040S	ea.
10½" (26.7 cm) L, 4¾" (11.8 cm) W, Medium Shoe size: Men's 8-10, Women's 9-11	ORT290040M	ea.
11" (27.9 cm) L, 5" (12.7 cm) W, Large Shoe size: Men's 10½-12½, Women's 11½-13½	ORT290040L	ea.
12½" (31.8 cm) L, 5¾" (13.7 cm) W, X-Large Shoe size: Men's 13-16	ORT290040XL	ea.

* Measurements based upon product size dimensions and relative shoe sizing

Post-Op Shoes



- A** Vinyl Hook-and-Loop Post-Op Shoe is easy to apply
- Rigid sole prevents excessive foot motion after surgery
 - Heel and toe pitch encourages a natural gait
 - Vinyl upper resists water and wipes clean
 - Skid resistant bottom provides safe traction
 - Two hook-and-loop straps make it quick and easy to apply

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30200MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30200MM	ea.
Large: 13-14	ORT30200ML	ea.
X-Large: 14½-15½	ORT30200MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30200WXS	ea.
Small: Women's 6-7½	ORT30200WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30200WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30200WL	ea.



- C** Semi-Rigid Post-Op Shoe is our most popular style
- Cushioned rocker sole with heel and toe pitch for natural gait
 - Durable nylon upper with a soft inner lining
 - Molded heel collar secures foot placement and reduces heel slipping
 - Loop-lock closure for an easy, secure fit
 - Removable tongue allows for optimum placement and comfort

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30300MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30300MM	ea.
Large: 13-14	ORT30300ML	ea.
X-Large: 14½-15½	ORT30300MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30300WXS	ea.
Small: 6-7½	ORT30300WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30300WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30300WL	ea.



- B** Vinyl Lace-Up Post-Op Shoe offers protection and security
- Lace-up design provides security and a customizable size
 - Rigid sole prevents excessive foot motion after surgery
 - Heel and toe pitch encourages a natural gait
 - Vinyl upper resists water and wipes clean
 - Skid-resistant bottom provides safe traction

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30100MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30100MM	ea.
Large: 13-14	ORT30100ML	ea.
X-Large: 14½-15½	ORT30100MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30100WXS	ea.
Small: 6-7½	ORT30100WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30100WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30100WL	ea.



- D** Medline's Rigid Post-Op Shoe offers protection and support during recovery
- Rigid sole reduces unwanted flexion during recovery from broken bones
 - Nylon upper with soft inner lining for patient comfort
 - Loop-lock closure for easy, secure fit

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30400MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30400MM	ea.
Large: 13-14	ORT30400ML	ea.
X-Large: 14½-15½	ORT30400MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30400WXS	ea.
Small: 6-7½	ORT30400WS	ea.
Medium: 8-9 (Men's 7-9)	ORT30400WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30400WL	ea.

* Measurements based upon relative shoe sizing