

ORT31300

*Meeting your traction needs  
from head to toe*

# Traction Products



A

**A** Bucks Traction Boot is universally sized to fit most patients

- Removable aluminum medial and lateral stays
- Pile foam outer with convoluted foam inner lining
- Four 3" (7.6 cm) latex-free elastic straps with hook-and-loop closure prevent migration
- 6" (15.2 cm) spreader bar attached
- Measurements based upon calf circumference

Description	Item Number	Qty
Up to 20" (50.8 cm) Circ., Universal	ORT31100	ea.



B

**B** Over-the-Door Cervical Traction Kit

- 10 feet of traction cord (3.04 m)
- 12" (30.5 cm) spreader bar
- Two heavy-duty pulleys
- Door bracket
- 20 lb (9.1 kg) water bag
- Deluxe head halter

Description	Item Number	Qty
One Size	ORT31400	ea.



C

Deluxe Head Halter is made with high quality materials for patient comfort

- C**
- Constructed of tri-laminate tricot, foam and cotton duck material for security and stretch-resistance
  - Hook-and-loop closures for easy application
  - Vertical pulls have metal D-rings for attachment to traction cord

Description	Item Number	Qty
One Size Fits Most	ORT31300	ea.



D

**D** Medline's Traction Cord is made with strong and durable nylon

Description	Item Number	Qty
100' (30.48 m) Roll	ORT31500	ea.