



*Basic to deluxe – many styles
to fit your needs and budget*

ORT11300

Arm Slings



A Deep Pocket Arm Sling maximizes patient comfort and ease of use

- Large pocket made of durable poly/cotton easily accommodates bandages and casts
- Wide shoulder strap with shoulder pad reduces neck and shoulder strain
- Thumb loop and adjustable shoulder strap inhibits sling migration
- Double D-ring with hook-and-loop closure make application quick and easy

Description	Item Number	Qty
10 ³ / ₄ " L x 6 ¹ / ₂ " W (27.3 x 16.5 cm), X-Small	ORT11300XS	ea.
13 ³ / ₄ " L x 7 ¹ / ₂ " W (34.9 x 19.1 cm), Small	ORT11300S	ea.
16 ¹ / ₂ " L x 8 ¹ / ₂ " W (41.9 x 21.6 cm), Medium	ORT11300M	ea.
20" L x 9 ³ / ₄ " W (50.8 x 24.8 cm), Large	ORT11300L	ea.
23 ¹ / ₄ " L x 9 ³ / ₄ " W (59.1 x 24.8 cm), X-Large	ORT11300XL	ea.
26 ¹ / ₂ " L x 9 ³ / ₄ " W (67.3 x 24.8 cm), 2X-Large	ORT113002XL	ea.

Also available in Dispenser Boxes

Description	Item Number	Qty
13 ³ / ₄ " L x 7 ¹ / ₂ " W (34.9 x 19.1 cm), Small	ORT11300SCS	12/box, 2 bx/cs
16 ¹ / ₂ " L x 8 ¹ / ₂ " W (41.9 x 21.6 cm), Medium	ORT11300MCS	12/box, 2 bx/cs
20" L x 9 ³ / ₄ " W (50.8 x 24.8 cm), Large	ORT11300LCS	12/box, 2 bx/cs
23 ¹ / ₄ " L x 9 ³ / ₄ " W (59.1 x 24.8 cm), X-Large	ORT11300X LCS	12/box, 2 bx/cs



B Keep your patients cool and comfortable with Medline's Cool Mesh Arm Sling

- Offers the same generous design and features as our Deep Pocket Arm Sling (above)
- Cool polyester mesh increases ventilation for patient comfort

Description	Item Number	Qty
7 ³ / ₄ " L x 5 ¹ / ₂ " W (19.7 x 14 cm), X-Small	ORT11200XS	ea.
10 ³ / ₄ " L x 6 ¹ / ₂ " W (27.3 x 16.5 cm), Small	ORT11200S	ea.
13 ³ / ₄ " L x 7 ¹ / ₂ " W (34.9 x 19.1 cm), Medium	ORT11200M	ea.
16 ¹ / ₂ " L x 8 ¹ / ₂ " W (41.9 x 21.6 cm), Large	ORT11200L	ea.
20" L x 9 ³ / ₄ " W (50.8 x 24.8 cm), X-Large	ORT11200XL	ea.

Also available in Dispenser Boxes

Description	Item Number	Qty
10 ³ / ₄ " L x 6 ¹ / ₂ " W (27.3 x 16.5 cm), Small	ORT11200SCS	12/box, 2 bx/cs
13 ³ / ₄ " L x 7 ¹ / ₂ " W (34.9 x 19.1 cm), Medium	ORT11200MCS	12/box, 2 bx/cs
16 ¹ / ₂ " L x 8 ¹ / ₂ " W (41.9 x 21.6 cm), Large	ORT11200LCS	12/box, 2 bx/cs
20" L x 9 ³ / ₄ " W (50.8 x 24.8 cm), X-Large	ORT11200X LCS	12/box, 2 bx/cs



C Universal Arm Sling adjusts to accommodate many sizes

- Double D-ring with hook-and-loop closure makes it easy to apply
- Made of a breathable polyester knit mesh
- Wide shoulder strap adjusts in length for appropriate fit
- Pocket can be folded in and secured with hook-and-loop tabs for size adjustment
- Includes shoulder pad

Description	Item Number	Qty
19 ¹ / ₈ " L x 7 ³ / ₄ " W (49.1 x 19.7 cm), Universal	ORT11010	ea.

Also available in Dispenser Boxes

Description	Item Number	Qty
19 ¹ / ₈ " L x 7 ³ / ₄ " W (49.1 x 19.7 cm), Universal	ORT11010CS	12/box, 2 bx/cs

* Measurements based upon the length between the end of the patient's elbow and the center of their shortest finger with an extended hand

D Standard Arm Sling makes your patients and your budget happy

- Constructed of cool and comfortable poly/cotton
- Strap adjusts with a metal slide buckle for proper fit

Description	Item Number	Qty
10 ² / ₅ " L x 5" W (27.1 x 12.7 cm), Child	ORT11100XS	ea.
13 ¹ / ₂ " L x 5 ¹ / ₂ " W (34.3 x 14 cm), Small	ORT11100S	ea.
16" L x 6 ¹ / ₄ " W (40.6 x 15.9 cm), Medium	ORT11100M	ea.
18 ¹ / ₂ " L x 6 ³ / ₄ " W (46.6 x 17.2 cm), Large	ORT11100L	ea.
20 ¹ / ₈ " L x 7 ¹ / ₂ " W (51.6 x 19.1 cm), X-Large	ORT11100XL	ea.

Also available in Dispenser Boxes

Description	Item Number	Qty
10 ² / ₅ " L x 5" W (27.1 x 12.7 cm), X-Small	ORT11100XSCS	12/box, 2 bx/cs
13 ¹ / ₂ " L x 5 ¹ / ₂ " W (34.3 x 14 cm), Small	ORT11100SCS	12/box, 2 bx/cs
16" L x 6 ¹ / ₄ " W (40.6 x 15.9 cm), Medium	ORT11100MCS	12/box, 2 bx/cs
18 ¹ / ₂ " L x 6 ³ / ₄ " W (47 x 17 cm), Large	ORT11100LCS	12/box, 2 bx/cs
20 ¹ / ₈ " L x 7 ¹ / ₂ " W (51.6 x 19 cm), X-Large	ORT11100XLCS	12/box, 2 bx/cs



E Universal Cradle Style Arm Sling allows for easy application

- Made of durable poly/cotton material
- Open elbow design
- Straps can be adjusted to accommodate many sizes

Description	Item Number	Qty
One Size Fits Most, Universal	MDT823210	ea.



E Medline's Universal Triangular Arm Sling is our most affordable arm sling alternative

- Made of durable cotton
- Easy to wash

Description	Item Number	Qty
36" L x 18" D (91.4 x 46 cm)	MDT824085	12/cs



* Measurements based upon the length between the end of the patient's elbow and the center of their shortest finger with an extended hand

Arm Slings



A Make your young patients smile with Medline's Pediatric Print Arm Sling

- Poly/cotton blended fabric with fun print
- 1"(2.5 cm) wide adjustable shoulder strap with foam pad ensure proper fit while reducing neck and shoulder strain
- Hook-and-loop closure for quick and easy application
- Thumb loop inhibits sling migration

Description	Item Number	Qty
9" L x 5" W (22.9 x 12.7 cm), Infant	ORT11400XXS	ea.
11" L x 6" W (27.9 x 15.2 cm), Toddler	ORT11400XS	ea.
13" L x 7" W (33 x 17.8 cm), Child	ORT11400S	ea.

Also available in Dispenser Boxes

Description	Item Number	Qty
9" L x 5" W (22.9 x 12.7 cm), Infant	ORT11400XXSCS	12/box, 2 bx/cs
11" L x 6" W (27.9 x 15.2 cm), Toddler	ORT11400XSCS	12/box, 2 bx/cs
13" L x 7" W (33 x 17.8 cm), Child	ORT11400SCS	12/box, 2 bx/cs



B Help reduce edema of the hand and arm with Medline's Arm Elevator

- Use as a sling or as an elevator
- Two 12" (30.5 cm) pockets allow for convenient application of ice bags or instant cold packs
- Made with durable cotton material
- Hook-and-loop closure for securing the arm pocket
- Metal slide buckle adjusts shoulder strap for use as a sling
- Vertical strap with D-ring attaches easily to traction devices or IV stands for elevation

Description	Item Number	Qty
16" L x 10" W (40.6 x 25.4 cm), Small	ORT11600S	ea.
17½" L x 11" W (44.5 x 27.9cm), Large	ORT11600L	ea.



C Hemiplegic Sling offers shoulder support for hemiplegic patients

- Simple, one-strap design makes it easy to apply
- Metal slide buckle adjusts strap length for proper fit
- Elbow and wrist cuffs support the weight of the arm

Description	Item Number	Qty
One Size Fits Most, Universal	ORT11700	ea.

* Measurements based upon the length between the end of the patient's elbow and the center of their shortest finger with an extended hand

D The low profile of Medline's Elastic Shoulder Immobilizer makes it cool and comfortable

- Contoured design fits both male and female patients
- Elastic chest strap stretches to best fit the patient
- Humeral and wrist cuffs are made of soft foam backed with a durable poly/cotton material
- Hook-and-loop closure on the cuffs adjust for proper fit
- Designed for right or left arm

Description	Item Number	Qty
23-27" (58.4-68.6 cm), X-Small	ORT16100XS	ea.
27-31" (68.6-78.7 cm), Small	ORT16100S	ea.
31-35" (78.7-88.9 cm), Medium	ORT16100M	ea.
35-39" (88.9-99.1 cm), Large	ORT16100L	ea.
39-43" (99.1-109.2 cm), X-Large	ORT16100XL	ea.
43-47" (109.2-119.4 cm), 2X-Large	ORT161002XL	ea.
47-51" (119.4-129.5 cm), 3X-Large	ORT161003XL	ea.



E Universal Cut-Away Shoulder Immobilizer provides a customizable fit

- Made of soft foam that can be cut to customize the fit
- Contoured design to fit men and women
- Cuffs and shoulder straps attach with hook-and-loop so they can be positioned where they are most effective
- Extra wide shoulder strap adjusts with metal slide buckle

Description	Item Number	Qty
Up to 53" (134.62 cm), Universal	ORT16000	ea.



F Sling and Swathe Immobilizer provides comfort and security during rehabilitation

- Made with soft pile foam that is comfortable against the skin
- Two over the shoulder straps criss-cross and attach with contact closures
- Swathe is made with a contact closure for accurate fit
- May be used as an arm sling if swathe is not desired
- Fits right or left arm

Description	Item Number	Qty
30-46" (76.2-116.8 cm), Small/Medium	ORT16020SM	ea.
32-52" (81.3-132.1 cm), Large/X-Large	ORT16020LXL	ea.



G Our Universal Sling and Swathe Immobilizer can accommodate a wide variety of patients

- Offers many of the same features as our sized Sling and Swathe
- Metal slide buckles secure the shoulder straps and attach to the sling with contact closures

Description	Item Number	Qty
Up to 53" (134.62 cm), Universal	ORT16010	ea.



* Measurements based upon chest size circumference

Shoulder Immobilizers & Elbow Supports



A Sling Style Shoulder Immobilizer provides the comfort and support of a sling

- Waist and shoulder straps attach using a metal clip buckle
- Made of durable cotton
- Metal slide buckle on body strap allows for easy size adjustment
- Fits right or left arm
- Measurements based upon distance between the end of patient's elbow and center of their shortest finger with an extended hand

Description	Item Number	Qty
12" (30.5 cm) L x 4½" (11.4 cm) W, X-Small	ORT16200XS	ea.
14" (35.6 cm) L x 5½" (14 cm) W, Small	ORT16200S	ea.
16" (40.6 cm) L x 6½" (16.5 cm) W, Medium	ORT16200M	ea.
18" (45.7 cm) L x 7¾" (17.8 cm) W, Large	ORT16200L	ea.
18" (45.7 cm) L x 7¾" (17.8 cm) W, X-Large	ORT16200XL*	ea.

* The X-Large has a 5" (12.7 cm) longer shoulder strap



B Shoulder Immobilizer with Abduction Pillow provides post-surgical support

- Adjustable shoulder strap includes neck pad to reduce shoulder and neck strain
- Waist strap attaches to the pillow using a plastic delrin buckle
- Removable pillow provides 25 degrees of abduction
- Fits right or left arm
- Measurements based upon distance between the end of patient's elbow and center of their shortest finger with an extended hand

Description	Item Number	Qty
14" (35.6 cm) L x 6" (15.2 cm) W, Small	ORT16300S	ea.
16" (40.6 cm) L x 7" (17.8 cm) W, Medium	ORT16300M	ea.
18" (45.7 cm) L x 8" (20.3 cm) W, Large	ORT16300L	ea.
20" (50.8 cm) L x 9" (22.9 cm) W, X-Large	ORT16300XL	ea.



C Medline's Tennis Elbow Strap help reduce pain caused by tennis elbow

Standard (pictured above, left)

- 2" W x 29" L (5.1 x 73.7 cm) compression strap
- Loop lock closure

Deluxe (not pictured)

- 2" W x 32.5" L (5.1 x 82.6 cm), neoprene cuff with sewn-in pressure pad
- Hook tab on cuff facilitates one-hand application
- 15" (38.1 cm) attached compression strap with loop lock closure
- Measurements based upon circumference of the forearm

D Neoprene Elbow Support provides warmth, support and compression

- Made of ⅛" (0.32 cm) black neoprene fabric
- Anatomically contoured to reduce bunching
- Measurements based upon circumference below the elbow

Description	Item Number	Qty
7-9" (17.8-22.9 cm) Circ., Small	ORT17200S	ea.
9-11" (22.9-27.9 cm) Circ., Medium	ORT17200M	ea.
11-13" (27.9-33 cm) Circ., Large	ORT17200L	ea.
13-15" (33-38.1 cm) Circ., X-Large	ORT17200XL	ea.
15-17" (38.1-43.1 cm) Circ., 2X-Large	ORT172002XL	ea.

Description	Item Number	Qty
7-14" (17.8 x 35.6 cm), Standard	ORT17100	ea.
8½"-14½" (21.6 X 36.8 cm), Deluxe	ORT17110	ea.



Specialized features provide maximum support and comfort

ORT18210

Wrist & Forearm Splints



- A** 10" Universal Wrist & Forearm Splint is easily adjusted for comfort and support
- Made with thick plush foam for optimum comfort
 - Top stay is moveable and attaches with a hook backing to allow for size adjustment
 - Sewn-in aluminum palmar stay supports the wrist in a natural cock-up position



- B** Medline's 10" Vinyl Wrist & Forearm Splint brings durability and comfort to a new level
- Flannel-lined, perforated vinyl resists wear and tear while providing necessary ventilation
 - Circumferential wrist strap and two forearm straps for compression and support
 - Splint does not extend beyond the palmar crease, allowing full finger range-of-motion
 - Slip-on design for ease of application
 - Sewn-in dorsal & palmar stays for immobilization

Description	Item Number	Qty
Right		
10" (25.4 cm) L, 4-11" (10.16 -27.9 cm) Circ., Universal	ORT18000R	ea.
Left		
10" (25.4 cm) L, 4-11" (10.16 -27.9 cm) Circ., Universal	ORT18000L	ea.

Description	Item Number	Qty
Right		
10" (25.4 cm) L, 4-5" (10.16-12.7 cm) Circ., X-Small	ORT18100RXS	ea.
10" (25.4 cm) L, 5-6" (12.7-15.2 cm) Circ., Small	ORT18100RS	ea.
10" (25.4 cm) L, 6-7" (15.2-17.8 cm) Circ., Medium	ORT18100RM	ea.
10" (25.4 cm) L, 7-8" (17.8-20.3 cm) Circ., Large	ORT18100RL	ea.
10" (25.4 cm) L, 8-9" (20.3-22.9 cm) Circ., X-Large	ORT18100RXL	ea.
Left		
10" (25.4 cm) L, 4-5" (10.16 -12.7 cm) Circ., X-Small	ORT18100LXS	ea.
10" (25.4 cm) L, 5-6" (12.7-15.2 cm) Circ., Small	ORT18100LS	ea.
10" (25.4 cm) L, 6-7" (15.2-17.8 cm) Circ., Medium	ORT18100LM	ea.
10" (25.4 cm) L, 7-8" (17.8-20.3 cm) Circ., Large	ORT18100LL	ea.
10" (25.4 cm) L, 8-9" (20.3-22.9 cm) Circ., X-Large	ORT18100LXL	ea.

C Deluxe 10" Wrist & Forearm Splint looks sharp and works hard

- Lined, faux suede material reduces migration
- Circumferential wrist strap and two forearm straps for compression and support
- Splint does not extend beyond the palmar crease, allowing full finger range-of-motion
- Slip-on design for ease of application
- Sewn-in dorsal and palmar stays for immobilization



Description	Item Number	Qty
Right		
10" (25.4 cm) L, 5-6½" (12.7-16.5 cm) Circ., X-Small	ORT18110RXS	ea.
10" (25.4 cm) L, 6-7" (15.2-17.8 cm) Circ., Small	ORT18110RS	ea.
10" (25.4 cm) L, 6½-8" (16.5-20.3 cm) Circ., Medium	ORT18110RM	ea.
10" (25.4 cm) L, 7-9" (17.8-22.9 cm) Circ., Large	ORT18110RL	ea.
10" (25.4 cm) L, 8½-10" (21.6-25.4 cm) Circ., X-Large	ORT18110RXL	ea.
Left		
10" (25.4 cm) L, 5-6½" (12.7-16.5 cm) Circ., X-Small	ORT18110LXS	ea.
10" (25.4 cm) L, 6-7" (15.2-17.8 cm) Circ., Small	ORT18110LS	ea.
10" (25.4 cm) L, 6½-8" (16.5-20.3 cm) Circ., Medium	ORT18110LM	ea.
10" (25.4 cm) L, 7-9" (17.8-22.9 cm) Circ., Large	ORT18110LL	ea.
10" (25.4 cm) L, 8½-10" (21.6-25.4 cm) Circ., X-Large	ORT18110LXL	ea.

D Immobilize injured wrists and thumbs with our 8" Wrist & Forearm Splint with Abducted Thumb

- Nylon outer material is lined with a plush foam for exceptional comfort
- Dual stays support both the palmar surface of the wrist and the extensor surface of the thumb
- Thumb support adjusts with hook-and-loop closure
- Three forearm straps secure with loop lock closure for a high level of compression



Description	Item Number	Qty
Right		
8" (20.3 cm) L, 5-6½" (12.7-16.5 cm) Circ., X-Small	ORT18210RXS	ea.
8" (20.3 cm) L, 6-7" (15.2-17.8 cm) Circ., Small	ORT18210RS	ea.
8" (20.3 cm) L, 6½-8" (16.5-20.3 cm) Circ., Medium	ORT18210RM	ea.
8" (20.3 cm) L, 7-9" (17.8-22.9 cm) Circ., Large	ORT18210RL	ea.
8" (20.3 cm) L, 8½-10" (21.6-25.4 cm) Circ., X-Large	ORT18210RXL	ea.
Left		
8" (20.3 cm) L, 5-6½" (12.7-16.5 cm) Circ., X-Small	ORT18210LXS	ea.
8" (20.3 cm) L, 6-7" (15.2-17.8 cm) Circ., Small	ORT18210LS	ea.
8" (20.3 cm) L, 6½-8" (16.5-20.3 cm) Circ., Medium	ORT18210LM	ea.
8" (20.3 cm) L, 7-9" (17.8-22.9 cm) Circ., Large	ORT18210LL	ea.
8" (20.3 cm) L, 8½-10" (21.6-25.4 cm) Circ., X-Large	ORT18210LXL	ea.

* Measurements based upon the circumference of patient's wrist

Wrist & Forearm Splints



E

E 7" Lace-Up Wrist Splint is remarkably easy to self-apply

- Durable cotton canvas is lined with soft flannel
- Laces pull tight and secure with hook-and-loop tabs
- Aluminum palmar stay provides added support

Description	Item Number	Qty
Right		
7" (17.8 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19300RXS	ea.
7" (17.8 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19300RS	ea.
7" (17.8 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19300RM	ea.
7" (17.8 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19300RL	ea.
7" (17.8 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19300RXL	ea.
Left		
7" (17.8 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19300LXS	ea.
7" (17.8 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19300LS	ea.
7" (17.8 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19300LM	ea.
7" (17.8 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19300LL	ea.
7" (17.8 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19300LXL	ea.

G Reduce the number of products you keep on hand with our 6" Universal Wrist Splint

- Made with thick plush foam for optimum comfort
- Top stay is moveable and attaches with a hook backing to allow for size adjustment
- Sewn-in aluminum palmar stay supports the wrist in a natural cock-up position



G

Description	Item Number	Qty
Right		
6" (15.2 cm) L, 4 1/2-10 1/2" (11.4-26.7 cm) Circ., Universal	ORT19000R	ea.
Left		
6" (15.2 cm) L, 4 1/2-10 1/2" (11.4-26.7 cm) Circ., Universal	ORT19000L	ea.

F Lightweight and secure is what Medline's 8" Wrist Splint is all about

- Nylon lined foam material is soft and breathable
- Slip-on design with three loop lock closures to achieve the desired level of compression
- Removable and adjustable aluminum palmar stay



F

Description	Item Number	Qty
Right		
8" (20.3 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19400RXS	ea.
8" (20.3 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19400RS	ea.
8" (20.3 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19400RM	ea.
8" (20.3 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19400RL	ea.
8" (20.3 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19400RXL	ea.
Left		
8" (20.3 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19400LXS	ea.
8" (20.3 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19400LS	ea.
8" (20.3 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19400LM	ea.
8" (20.3 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19400LL	ea.
8" (20.3 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19400LXL	ea.

H 7" Elastic Wrist Splint provides compressive support

- Latex-free cotton elastic is soft and durable
- Contact closure for ease of application
- Malleable palmar stay for added support and proper fit



H

Description	Item Number	Qty
Right		
7" (17.8 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19100RXS	ea.
7" (17.8 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19100RS	ea.
7" (17.8 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19100RM	ea.
7" (17.8 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19100RL	ea.
7" (17.8 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19100RXL	ea.
Left		
7" (17.8 cm) L, 4 1/2-5 1/2" (11.4-14 cm) Circ., X-Small	ORT19100LXS	ea.
7" (17.8 cm) L, 5 1/2-6 1/2" (14-16.5 cm) Circ., Small	ORT19100LS	ea.
7" (17.8 cm) L, 6 1/2-7 1/2" (16.5-19.1 cm) Circ., Medium	ORT19100LM	ea.
7" (17.8 cm) L, 7 1/2-8 1/2" (19.1-21.6 cm) Circ., Large	ORT19100LL	ea.
7" (17.8 cm) L, 8 1/2-9 1/2" (21.6-24.1 cm) Circ., X-Large	ORT19100LXL	ea.

* Measurements based upon the circumference of patient's wrist