


**In the
selection
of syringes,
some
choices
are
clear**

medrad®

...Medrad makes them clearer and therefore safer.

Ever since *Medrad* introduced the first disposable angiographic syringes, we've taken special pride in making them better. As needs have changed over the years, we've made *Medrad Syringes* easier to use, more convenient, and clearer, all of which makes them **safer. The key to *Medrad's* superior performance and reliability: *our* syringes are designed and calibrated as part of a **Total Fluid Delivery System**, meeting the highest standards of *Medrad* injectors. Because *Medrad* believes some choices should be clearer than others.**

CE 0123



With angiographic syringes, greater clarity often means greater safety. *Medrad's* Syringes are designed to make it easier for you to detect air bubbles. You can see the contrast inside the syringe, and a clear nut at the tip provides the best air bubble visualization available.

Quality is also clearly evident in every *Medrad* Syringe. Ribbed, reinforced necks provide extra strength. Operating plungers are designed to perform without sticking, alleviating uneven delivery. And the quick fill tubing shortens your preparation time, while giving you the most efficient setup possible.

Convenience is important in contrast delivery...it directly contributes to the quality of services you provide. That's why *Medrad* gives you a variety of syringe sizes to meet your patients' individual needs. We also design and calibrate *our* syringes so they are easier to connect to *Medrad* injectors and deliver contrast smoothly according to your exact requirements.

- It's more than another syringe—it's designed and calibrated to work with the *Medrad* injector as an integral part of our ***Total Fluid Delivery System.***
- *Medrad's* history of consistent quality ensures that our syringes provide superior performance with each procedure.
- Because procedures and dosages differ, *Medrad* Syringes offer you flexibility with multiple sizes.
- *Medrad's* complete range of syringes fit all major injectors, so you can choose the proper syringe to meet the exact needs of *the procedure* and *the patient.*


At *Medrad*, we believe there can be no compromise in quality when designing a contrast delivery system. It's why *Medrad* is committed to *Quality for Life.*®

To find out more about *Medrad* Syringes and the benefit of utilizing our *Total Fluid Delivery Systems*, contact your local *Medrad* Representative at 1-412-767-2400. Or visit our website: www.medrad.com

How to Order MEDRAD Quality.

The Qwik-Fit Syringe® Disposable provides:

- Quick disconnect, greater visualization, and reduced injector maintenance
- Reduced exposure to blood-borne pathogens
- One handed connection of low-pressure tubing
- Efficient contrast delivery with shortneck design



Medrad 200ml Qwik-Fit Syringe® Disposables CT Custom Kit



Catalog Number

MEDRAD®

The only source for Medrad Engineering, Medrad Quality, Medrad Service & Medrad Support.



Medrad is dedicated to continually improving the quality of all our products and services such that our customers' satisfaction, loyalty, and respect are unsurpassed. It is our policy to clearly understand and agree upon the valid requirements of the work we perform for our customers, both internal and external, and to pursue 100% conformance to those requirements...on time, every time.

Medrad reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Medrad representative for the most current information.

Medrad, Inc. is EN-ISO 9001/EN 46001 Certified
Medrad Europe B.V. is EN-ISO 9002/EN 46002 Certified

Medrad, Quality for Life, Q and design, Qwik-fit Syringe, Spectris, Mark V ProVis, Medrad Vistron CT and Medrad Pulsar are federally registered trademarks and Mark IV, Mark V, Mark V Plus, MCT, MCT Plus, EnVision CT and Solaris are trademarks and Medrad is a service mark of Medrad, Inc. U.S.A.

Medrad, Inc.

One Medrad Drive
Indianola, PA 15051-0780 U.S.A.
Telephone: 1 (412) 767-2400

Customer Service/Orders

1-800-MEDRAD-1
(1-800-633-7231)

Customer Service FAX

1 (412)767-4120

International FAX

1 (412) 767-4128

Medrad Europe B.V.

P.O. Box 205 • 6190 AE Beek
The Netherlands
Telephone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

www.medrad.com

E-mail address for more information: info@medrad.com

Injector	Vol.	Style	Syringe Alone	Syringe with Quick Fill Tube	Syringe with Filling Device and Low Pressure Tubing
Mark V ProVis® & Mark V Plus™**	60 ml	FasTurn		60-FT-Q	
	150 ml	FasTurn	150-FT	150-FT-Q	
	200 ml	FasTurn	200-FT	200-FT-Q	
Mark V™*	150 ml	FasTurn	150-FT	150-FT-Q	
	200 ml	FasTurn	200-FT	200-FT-Q	
Mark IV™ – All Series	130 ml	FasTurn	DSK 130	DSK 130-Q	
	130 ml	Push-N-Turn		130 AML-P-Q	
Medrad CT 200/300 Series Injector*	200 ml	FasTurn	200-FT	200-FT-Q	CTP-100
100 Series	130 ml	FasTurn	DSK 130	DSK 130-Q	
	130 ml	Push-N-Turn		130 AML-P-Q	
Medrad MCT™/MCT Plus™ Front Load Series	200 ml	Luer	200-FLS	200-FLS-Q	CTP-200-FLS
Medrad EnVision CT™ & Medrad Vistron CT®	125 ml	Luer		125-FLS-Q	CTP-125-FLS
	200 ml	Luer	200-FLS	200-FLS-Q	CTP-200-FLS
Medrad Spectris® MR	65 ml	Luer			SQK 65s (note 1) SQK 65vs (note 2)
Medrad Spectris® Solaris™ MR	65 ml 115 ml	Luer Luer			SSQK 65/115vs
Medrad Pulsar® Ultrasound Injector	30 ml	Luer			UDK 50 (note 3)
Model 6000 and 6000 CT Liebel-Flarsheim	150 ml	FasTurn		150-FT-A-Q	CTP-150-A
	260 ml	FasTurn		260-FT-Q	CTP-260
Angiomat Model 3000 and Viamonte-Hobbs 2000**	100 ml	FasTurn	100-FT-A	100-FT-A-Q	

Notes: 1. Kit contains 2 syringes, spikes for loading contrast media and saline, and 96" Y tubing **without** back check valve.
2. Kit contains 2 syringes, spikes for loading contrast media and saline, and 96" Y tubing **with integrated back check valve.**

3. Medrad Pulsar Ultrasound injector is available in Europe, Japan and Canada. Not available in the U.S.A.

* Medrad Syringes 150-FT, 200-FT are patented products.
** Angiomat is a trademark of the Liebel-Flarsheim Corporation.