


MILTEX® HI-LIGHT

SURGICAL INSTRUMENT CASSETTES

MILTEX® INSTRUMENT STERILIZATION & STORAGE SYSTEM

FEATURING THE "TUCK-AWAY™" LID



FEATURES AND BENEFITS

1. Flow-through design allows for effective rinsing and cleaning.
2. Medical grade silicone for secure instrument retention (comes standard).
3. Full length hinge provides stable and secure attachment to cassette.
4. Easy push button release provides a positive lock.
5. Pre-stamped tabs for custom coding of cassettes.
6. Double-hinged lid folds for easy storage.
7. Self-lubricating plastic hinge pin will not corrode from repeated sterilization.
8. Passivated stainless steel provides corrosion protection and aesthetic finish.

UNIVERSAL CASSETTE WITHOUT MAT

REF	Size
4-6835	6-1/2 x 6-1/4 x 1-3/16
4-083000	8 x 3-13/16 x 1-3/16
4-084000	8 x 4-7/8 x 1-3/16
4-085000	8 x 6 x 1-3/16
4-089100	8 x 9 x 1-3/16
4-008122	8 x 11 x 1-3/16
4-008115	8 x 11 x 2
4-009126	9-1/8 x 13 x 1-3/16
4-6815	15 x 8-1/4 x 1-3/16

SLIMLINE SERIES WITHOUT MAT

REF	Size
3-072007	7 x 2-5/8 x 11/16*
3-082008	8 x 2-5/8 x 11/16*
STD5L27	7 x 2-5/8 x 1*
3-072014	7 x 4-5/8 x 11/16*

*Single-hinge design

See reverse side for additional information

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. Authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.




+H834SURGCASSHI0T

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

LATEX-FREE SILICONE MATS



FEATURES AND BENEFITS

Cushion and protect your delicate instruments with these Latex-free mats. Customized to fit perfectly into your Miltex® Surgical Cassettes.

1. The raised mat (or feet) assures:
 - Adhesion of the mat to the tray surface is prevented
 - Allows for complete circulation under the mat
2. Holes positioned on mat for proper drainage
3. Soft, cushion fingers protect your instrument investment
4. Customized to fit perfectly into the Miltex® Cassettes



SILICONE MAT ONLY

Item#	Size	Description
4-6835MAT	6-1/2 x 6-1/4	for Cassette 4-6835
4-083000MAT	8 x 3-13/16	for Cassette 4-083000
4-084000MAT	8 x 4-7/8	for Cassette 4-084000
4-085000MAT	8 x 6	for Cassette 4-085000
4-089100MAT	8 x 9	for Cassette 4-089100
4-008122MAT	8 x 11	for Cassette 4-008122 & 4-008115
4-009126MAT	9-1/8 x 13	for Cassette 4-009126
4-6815MAT	15 x 8-1/4	for Cassette 4-6815

See reverse side for additional information

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. Authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.