

OLYMPIC SMART SCALE®

*A pediatric scale that safely
weighs all your patients...
from infants to toddlers*

Accurate & Easy-to-use

**Unique Design
Ensures Flexibility**

Safe & Efficient



PEDIATRIC SMART SCALE®

Model 50



Accurate Measurements

- Computes weight to the nearest ounce
- Automatically zeros itself when turned on and signals after three seconds when zeroing is complete
- Zero (tare) button allows extra weight, such as blankets or toys, to be cancelled out to accurately weigh child only

Easy to Use

- Fully Automatic — After placing child on the scale, weight is displayed continuously
- Averaging Button — Allows active children to be weighed within seconds
- Memory Recall Button — Enables display of last averaged weight, even if the scale has been turned off

Flexible & Efficient

- Weigh children within seconds, even when they are moving
- Using fold-down ends, patients can be weighed seated or lying down – from infant to 3 years old
- Operates on rechargeable batteries or connected to an electrical outlet while batteries are recharging
- Available as either a table-top model or as a roll-around model on a sturdy cart

Safe

- Features rounded corners & comfortable shape
- Fold-up sides ensure patient placed in scale is securely enclosed
- Stable weighted base provides anti-slip & anti-slide properties

natus
pediatrics

Natus Medical Incorporated
1501 Industrial Road • San Carlos, CA 94070 USA
1-800-303-0306 • 1-650-802-0400
www.natus.com

ROLL-AROUND PEDIATRIC SMART SCALE

Model 51



- Patient is protected with scale securely anchored to the sturdy roll-around base
- Padded crash-bar handles at both ends protect the scale when moved around
- Lower height enables easier placement of patients on the scale, while reducing potential back strain for the clinician
- Four fully-swiveling, ball-bearing casters enable easy maneuvering
- Supplies can be stored in large enclosed cabinet with transparent doors

Capacity: Infant scales safely hold up to 25,000 grams (55 lbs)

Battery Operation: Internal Ni-Cad batteries included. Automatically recharges when connected to an electrical outlet.

Electrical Requirements: 120V~, 60 Hz, 6-ft power cord with hospital-grade connector.

Electrical Safety: CSA C/US certified to C22.2 601.1 M90 and UL 60601-1.

Calibration: Recalibration procedure can be performed by hospital technicians. (Note: All scales require periodic calibration checks using standard weights.) NIST-traceable calibration weight available; see Ordering information.

Ordering Information

Olympic Smart Scale Model 50 – Table Top

Complete with internal rechargeable batteries.

CAT. NO. 56350\$3695.00

Olympic Smart Scale Model 51 – Roll-Around

Complete with internal rechargeable batteries.

CAT. NO. 56351\$4495.00

Calibration Weight

12.5 kg ±1.0 gm. NIST-traceable. Two nickel-plated steel weights.

CAT. NO. 56357\$325.00

For information or ordering
Call 1-800-303-0306
(toll-free in US/Canada)