

OLYMPIC SMART SCALES®

*The ultimate standard
for accurate weight*



SMART SCALE.
Weigh active babies
quickly & accurately



WARM SCALE.
Gently warm babies
during weighing



OLYMPIC SMART SCALES®

Olympic Smart Scales & Warm Scales combine remarkable accuracy with simple operation. To weigh a baby, simply place the child in the tray, and the scale will automatically compute weight to the nearest gram – even if he/she is moving.

With the Smart Scale, changing numbers do not have to be tracked — the final weight is automatically locked on the display. One can continue caring for the baby and chart the weight when convenient. Even if the infant is removed from the scale, weight remains locked on the display.

SMART SCALE & WARM SCALE

Accurate Measurements

- Automatically computes weight to the nearest gram
- Automatically zeros itself before each weighing for accurate reproducible weight measurements
- Zero (tare) button allows extra weight, such as blankets or toys, to be cancelled out to accurately weigh only the baby

Easy To Use

- Fully Automatic — After placing baby on the scale, weight is automatically computed & locked on the display
- Automatic Weight Averaging — When the scale detects motion, it automatically averages & computes within seconds, and locks the final weight on the display
- Memory Recall button — Allows display of the last averaged weight, even if the scale has been turned off
- Easy Lift Tray (Smart Scale only) — Features two hand ports that enable easier lifting of a baby in & out
- Smooth contours, molded shape & edgeless tray enable easy cleaning

Efficient

- Weigh babies within seconds, even when they are moving
- Weighs on rechargeable batteries or connected to an electrical outlet while batteries are charging
- Automatic shut down feature conserves battery life
- Available as either a table-top style or as a roll-around model on a sturdy cart



WARM SCALE




Warm & Safe

- Warmed air enters bottom of tray, rising all around baby to help protect against heat loss & eliminate need for other warming equipment
- The air warming is isolated from the tray for safety, while air temperature is thermostatically controlled
- No temperature adjustments required — Air temperature is pre-set & maintained automatically
- Warm Scale functions even when baby is on a blanket or liner, by keeping both warm
- Liquids, such as urine or disinfectants spilled in the tray, will not interfere with the warming function

Safe, secure & efficient weighing solutions



Intelligent Status Lights

-  Automatically zeroing.
-  Baby moving; scale computing weight.
-  Final, accurate weight is now displayed



ROLL-AROUND SMART SCALE / WARM SCALE

- Patient is protected with scale securely anchored to the sturdy roll-around base
- Padded crash-bar handles at both ends protects the scale when moved around
- Lower height makes lifting infants in & out of the scale safer for the baby, while reducing potential back strain for the clinician
- Four fully-swiveling, ball-bearing casters enable easy maneuvering
- Supplies can be stored in large enclosed cabinet with transparent doors

OLYMPIC SMART SCALES®

Additional Information

Capacity: Infant scales safely hold up to 15,000 grams (33 lbs)

Battery Operation: Internal Ni-Cad batteries included. Automatically recharges when connected to an electrical outlet.

Front panel “charging” signal: Scale can be used while batteries are charging.

Calibration: Easy calibration procedure can be performed by local technician. Standard calibration weight available; see ordering information.

Servicing: Field service, if required, can be performed by hospital technicians. Toll-free hotline to Olympic Medical.

Electrical Safety: CSA C/US certified to CSA 60601.1 and UL 60601-1.

Electrical Requirements: 120V~, 60 Hz, 0.2A, six-foot power cord with hospital-grade connector.

Ordering Information

Olympic Smart Scale Model 60 – Table-Top
Complete with nursery tray & rechargeable batteries.
CATALOG NO. 56360

Olympic Smart Scale Model 65 – Roll-Around
Complete with nursery tray & rechargeable batteries.
CATALOG NO. 56365

Olympic Warm Scale Model 23 – Table-Top
Complete with newborn safety warming tray & rechargeable batteries.
CATALOG NO. 56323

Olympic Warm Scale Model 28 – Roll-Around
Complete with newborn safety warming tray & rechargeable batteries.
CATALOG NO. 56328

Calibration Weight
7500 ±1 gm. Nickel plated steel.
CATALOG NO. 56315

Also Available:

OLYMPIC PEDIATRIC SMART SCALE MODELS 50 & 51

Table Top & Roll-Around Options

- Unique fold-down ends enable weighing of patients seated or lying down
- Weighs both infants & toddlers up to 25,000 grams (55 lbs)

Ordering Information

Olympic Smart Scale Model 50 – Table-Top
Complete with internal rechargeable batteries.
CAT. NO. 56350

Olympic Smart Scale Model 51 – Roll-Around
Complete with internal rechargeable batteries.
CAT. NO. 56351



*Olympic Smart Scale
Model 50*

*Olympic Smart Scale
Model 51*

natus
pediatrics

Natus Medical Incorporated
1501 Industrial Road • San Carlos, CA 94070 USA
1-800-303-0306 • 1-650-802-0400
www.natus.com

For information or ordering
Call 1-800-303-0306
(toll-free in US/Canada)