



11-675D Square Face Wall Aneroid

Durable and reliable, Omron's Wall Aneroids offer a LIFETIME calibration warranty and a long-life nylon cuff in a sturdy wall-mount to maintain the performance requirements of a busy clinical practice.

Features

- No-stop, 300 mmHg gauge and a deluxe air-release valve
- Latex (11-675D) and Latex FREE (11-675DNL) versions available
- Adult Cuff included; Infant, Pediatric, Large Adult, and Thigh size cuffs sold separately
- Includes tool for on-site calibration
- LIFETIME calibration warranty +/- 3 mmHg



11-685NL Floor Aneroid

Durable and reliable, Omron's Floor Aneroids offer a LIFETIME calibration warranty and a long-life nylon cuff in a sturdy wall-mount to maintain the performance requirements of a busy clinical practice.

Features

- Large, easy-to-read 6 faceplate
- High-impact housing is lightweight, yet durable
- Lifetime calibration warranty ± 3 mmHg
- Deluxe air release valve for optimal performance
- Latex and latex free versions
- Includes adjustable mobile floor stand