



HEM-907XL

Clinically validated for accuracy, the Omron IntelliSense® Professional Digital Blood Pressure Monitor is designed for clinical use in a professional setting including a physician's office, hospital, or other point of care environment.

Features

- Automatic cuff inflation and deflation eliminates the need to predetermine inflation level setting
- Average mode measures up to three readings and averages the total, for a more accurate measurement
- Greater choice of arm cuff circumferences with four cuff sizes included -- small, medium, large and extra-large
- Reliable power supply through included AC adapter or rechargeable battery pack
- Five-Year warranty