



G10VP

Longer Containers

Counter-top

How It Works



G14KA

Longer Containers

Counter-top

Abdominal Probes

Specifications:

Model	G10VP	G14KA
Number of probes soaked	1	2
Face Velocity FPM	45	45
Electrical (115V 50/60Hz)	.5A	.5A
Space & Weight Requirements		
Width	13"	18.5"
Depth	10"	10"
Height	59"	59"
Weight (lbs)	24	28

UL and CE approved. 230V models available. Counter top units have different space requirements