

NEW!

GUS™ Wall Mounted Soak Stations for Vaginal and Rectal Ultrasound Probes



Model G10VP for a single probe. One container is for soaking, the other is for the initial rinse.

Additional Benefits

- 12" custom molded containers fit all brands of vaginal and rectal probes
- Each container holds less than one quart, providing significant savings
- Completely self contained unit, no ductwork or plumbing required

For more information, please call
PCI medical at **800-862-3394**
www.pcimedical.com

Protect Your Staff

- Safely eliminates fumes from Glutaraldehyde and OPA
- Probes are suspended in containers preventing damage to the delicate tips
- Fumes are drawn away from the operator and into a patented carbon filter where they are neutralized
- Locking safety door allows you to install the system in a patient area
- Meets **OSHA** and **JCAHO** requirements



Model G14KA allows for two probes to soak at one time, middle container is for the initial rinse

PCI
MEDICAL