

Parker Labs Electromedical Products

SPECTRA 360® ELECTRODE GEL



The only salt-free and chloride-free electrically conductive gel, recommended for all electromedical procedures, except defibrillation. Salt-free characteristics make it particularly suitable for long-term applications. Spectra 360 differs significantly from all other electrically conductive media...it works by wetting the skin, thereby reducing skin resistance.

Features

- Non-gritty STAY-WET® formula allows for prolonged use without re-application
- Product of choice for ECG, TENS, monitoring and pediatrics
- Non-irritating, hypoallergenic, bacteriostatic
- Salt-free, no sodium ion transfer
- Product of choice for conductive rubber/carbon electrodes
- Not recommended for defibrillation

Product

Package Size

12-02

60 g (2 oz.) tube, 12 per box

12-08

250 g (8.5 oz.) tube, 12 per box



SIGNAGEL® ELECTRODE GEL

Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.

Features

- Highly conductive
- Hypoallergenic, bacteriostatic
- See-thru tube with stand-up cap lets you know how much SIGNAGEL is left
- Water soluble, non-staining, non-gritty

Product

Package Size

15-60

60 g (2 oz.) tube, 12 per box

15-25

250 g (8.5 oz.) tube, 12 per box



SIGNAPAD® PREMOISTENED ELECTRODE PADS

The ultimate in pre-moistened ECG electrode pads for single or multi-channel use. Dispensed individually or several-at-a-time from our sturdy multipad "pop-up" container.

Features

- Highly conductive, non-irritating
- Saves preparation and clean up time
- Easy to use, faster ECG's
- Eliminates residue build-up, electrodes stay clean
- Bacteriostatic
- Not recommended for defibrillation

Product # Package Size

16-40	44 x 50 mm pad, 400 pads per dispenser, 2 dispensers per box, 10 boxes per case
-------	---



SIGNACREME® ELECTRODE CREAM

Highly conductive cream electrolyte for use in electromedical procedures where a cosmetic quality cream is preferred.

Features

- Cosmetic quality, pleasing to patients
- Non-irritating, hypoallergenic, bacteriostatic
- Conveniently packaged in refillable see-thru dispenser with flip-top cap
- Not recommended for defibrillation

Product # Package Size

17-05	142 g (5 oz.) dispenser, 12 per box
17-20	Half U.S. Gallon(1.9 liter) with 2 dispensers and 1 dispenser pump



SIGNASPRAY® ELECTRODE SOLUTION AND SKIN PREP*

Highly conductive, economical spray electrolyte and non-gritty skin prep.

Features

- The ideal electrolyte for use with muscle stimulation equipment
- Superior non-gritty skin prep
- Bacteriostatic, hypoallergenic, non-staining
- No residue build-up
- Easy to apply, easy to clean up
- Clear bottle lets you know how much SIGNASPRAY is left
- Not recommended for defibrillation

Product # Package Size

18-25	250 ml (8.5 fl. oz.) clear spray dispenser bottle, 12 per box
18-28	3.8 liter (1 U.S. gallon) with dispenser, 4 per case (Dispenser pump not included)
18-04	SIGNASPRAYPAC® contains: 4 SIGNASPRAY gallons, 2 refillable spray dispensers and 1 dispenser pump



TENSIVE® CONDUCTIVE ADHESIVE GEL

Safe, non-flammable, odorless, conductive adhesive gel. Recommended for TENS, EMG, EEG and other electromedical procedures.

Features

- Eliminates tape and tape irritation
- Conductive immediately, no need to wait
- Non-flammable, no solvent odor
- Best adhesive gel available
- Hypoallergenic, bacteriostatic, non-irritating
- Water soluble, easily removed with water
- Not recommended for defibrillation

Product # Package Size

22-60	50g tube, 12 per box
-------	----------------------



REDUX® ELECTROLYTE GEL

Non-abrasive, water-based gel for routine ECG and defibrillation.

Features

- Highly conductive
- Water soluble, non-staining, non-gritty

Product #	Package Size
-----------	--------------

65-04

114g (4 oz.) tube, 10 per box



REDUX® ELECTROLYTE CREME

For routine ECG.

Features

- Cosmetic quality, pleasing to patients
- Non-irritating
- Highly conductive electrode cream
- Not recommended for defibrillation

Product #	Package Size
-----------	--------------

66-04

114g (4 oz.) tube, 10 per box



REDUX® ELECTROLYTE PASTE

Water-based, abrasive gel for short term monitoring, ECG and defibrillation.

Features

- Highly conductive
- Water soluble

Product #	Package Size
67-05	142g (5 oz.) tube, 10 per box