

Up

About Us

Specifications

Features

Our Price List

Replacement Parts

Our Warranty

Privacy Policy

Welcome To PediaGraph.com

The PediaGraph Immobilizer (PGI) has received high marks after being field tested at LeBonheur Children's Medical Center. Many of the features were developed in response to welcome suggestions from the professionals who actually use pediatric x-ray positioners on a daily basis. As infants don't understand the radiographic process, the fast positioning and easy - adjusting features of the PGI save time and tempers and get the job done properly.



- ◆ The PediaGraph Immobilizer incorporates ergonomic engineering features for operational ease and patient comfort.
- ◆ The PGI imposes no radiation risk to medical personnel.
- ◆ The positive locking turntable helps to eliminate motion in x-ray films.
- ◆ The turntable rotates on non-stick glide pads for easy positioning.

[Up](#)[About Us](#)[Specifications](#)[Features](#)[Our Price List](#)[Replacement Parts](#)[Our Warranty](#)[Privacy Policy](#)

LIMITED WARRANTY

Plastics Unlimited, Inc. (PUI) warrants this Pedia Graph Immobilizer against any defects in material or workmanship as follows.

1. For a period of (1) year from the date of purchase, if this Immobilizer or any parts thereof are determined to be defective, PUI will repair or replace the defective part (s), at its option, free of charge.
2. For this warranty to be effective, you must give PUI notice of the defective part within one year of the date of purchase and promptly return the defective part to PUI postage prepaid.
3. Acrylic parts are not warranted against breakage.

This warranty does not cover cosmetic damage or damage due to acts of God, accident, misuse, negligence, or modification of the product. The product warranty does not cover damage due to improper operation or maintenance, abnormal working conditions, or transfers to someone other than the original consumer purchaser.

EXCEPT AS EXPRESSLY SET FORTH HEREIN, PUI MAKE NO WARRANTIES EITHER EXPRESS OR IMPLIED.

PLASTICS UNLIMITED, INC.
(931) 722-7181

jeff@PediaGraph.com

or write to:

Plastics Unlimited, Inc.
1242 Savannah Hwy.
Waynesboro, Tn. 38485