

T-Spray II

Disinfectant Cleanser for Ultrasound and Mammography



T-SPRAY II Disinfectant Cleanser for Ultrasound and Mammography is the essential medical accessory. T-SPRAY II cleans and disinfects all ultrasound probes, transducers, mammography compression plates and cassettes after each use and between patients. T-SPRAY II has a neutral pH and contains no alcohol, phenol, glutaraldehyde or ketones. It is safe for equipment and is an effective bactericide and virucide killing HIV-1. The active ingredient is conveniently user activated by adding tap water to the fill line indicated on the bottle.

Unit Sizes:

- 60 ml Tube
- 250 ml Bottle