

NICOLET POCKET-DOP 3™

.....> **Accurate, durable, compact**



Nicolet
V A S C U L A R
a division of VIASYS Healthcare

.....> Nicolet Pocket-Dop 3™

Interchangeable, Convenient Probes

The Pocket-Dop 3 is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker



amplifies fetal heart sounds to reassure the expectant mother. The Pocket-Dop 3 features interchangeable probes: the 3 MHz probe is for early fetal heartbeat detection and the 2 MHz probe for use in later pregnancy or with larger patients. The small "foot-

print" of these probes enables them to be easily angled down behind the pubic bone for optimum detection of early fetal heartbeats. When not in use, the probe is attached to the side of the unit by a unique magnetic holder. Obstetrical exams are quick, easy and reliable with the Pocket-Dop 3.

Accuracy of Autocorrelation

The Pocket-Dop 3 features autocorrelation which is a pattern recognition technique that quickly acquires and accurately tracks the complex signal of the fetal heart rate. This reduces the time wasted waiting for the system to acquire

the best signal to read the fetal heart rate. The system averages the heart rate on command and then shows it in a high-contrast display panel. After three minutes of non-use, the unit turns itself off.



Smart Recharging System and Wall Mount

The Pocket-Dop 3 comes with a convenient recharging stand to be placed on a counter or mounted to a wall. This smart recharging system won't overcharge the Doppler which ensures maximum battery life. The Pocket-Dop 3 is always handy and ready for use.



Reliability and Customer Satisfaction

Nicolet Vascular is a leading manufacturer of hand-held Doppler systems. The Pocket-Dop 3 follows in the Nicolet tradition of technological innovation and superior customer service. The Pocket-Dop 3 is manufactured in our ISO 9001 certified facility and is covered by a 5-year parts and 1-year labor warranty. Like all of our products, it is backed by our Total Customer Satisfaction pledge: *At Nicolet Vascular, we're not satisfied until you are.*



Technical Specifications

Dimensions: 12.7 x 8.38 x 2.5 cm
5 x 3.3 x 1 inches

Weight: 0.285 Kg.
10 oz.

Power Source: 110 volt recharger, 9 volt rechargeable battery
220 volt recharger available

Catalog Numbers

2 MHz OB probe:	J200
3 MHz OB probe:	J300
2 & 3 MHz probes:	J230

CE
0344

ISO
9001
CERTIFIED

Nicolet
V A S C U L A R
a division of VIASYS Healthcare

Nicolet Vascular • 5225 Verona Road, Building 2 • Madison, WI 53711
1-800-525-2519 • Tel: 608-441-2266 • Fax: 608-441-2232
Websites: www.NicoletVascular.com • www.NicoletBiomedical.com
www.VIASYSHealthcare.com • Email: info@NicoletVascular.com

VIASYS™
HEALTHCARE
Excellence For Life
NEUROCARE GROUP