

Aneroid

Sphygmomanometers



Deluxe Aneroid Sphygmomanometer

Offering the most features, the Pro Advantage Deluxe Aneroid Sphygmomanometer includes a precision crafted manometer, cuff size marking system with extended graduation marks, filter protected bulb, dip molded inflation bladder for heavy use, chrome-plated brass air release valve and a zippered carry case. Lifetime calibration warranty, 2 year system warranty. Contains latex.

- P548200** Aneroid, Adult, Black, Latex
- P548210** Aneroid, Child, Black, Latex
- P548215** Aneroid, Large Adult, Black, Latex

Premium Pocket Aneroid Sphygmomanometer

The Pro Advantage Premium Pocket Aneroid Sphygmomanometer is a professional inflation system featuring a manometer, a cuff size marking system, chrome-plated brass air release valve, dip molded inflation bladder and bulb for heavy use and a zippered carry case. Lifetime calibration warranty, 2 year system warranty. Available in latex and latex-free.

- P548310** Aneroid, Adult, Black, Latex
- P548320** Aneroid, Child, Black, Latex
- P548350** Aneroid, Large Adult, Black, Latex
- P548340** Aneroid, Adult, Black, Latex-Free
- P548330** Aneroid, Child, Black, Latex-Free
- P548360** Aneroid, Large Adult, Black, Latex-Free



Standard Pocket Aneroid Sphygmomanometer

The Pro Advantage Standard Pocket Aneroid Sphygmomanometer is an economical inflation system featuring a manometer, chrome-plated deflation valve, nylon cuff, inflation bladder and bulb and a zippered carry case. 10 year calibration warranty, 1 year system warranty. Contains latex.

- P548440** Aneroid, Adult, Navy Blue, Latex
- P548445** Aneroid, Child, Navy Blue, Latex
- P548450** Aneroid, Large Adult, Navy Blue, Latex

Aneroid

Sphygmomanometers



Pocket Gauge



Cuff & Bladder, Latex



Cuff & Bladder, Latex-Free



Inflation System, Latex



Inflation System, Latex-Free

Sphygmomanometer Accessories

Accessories for Pro Advantage Sphygmomanometers. Available choices are a pocket gauge, cuffs and bladders, and inflation systems. Gauge features a black enamel finish. Cuff and bladder features a nylon cuff with marking system and a dip molded bladder. Inflation system features a nylon cuff with marking system, inflation bladder, filter protected bulb and valve. Two year warranty; Pocket Gauge has a lifetime calibration warranty. Available in latex and latex-free.

P548900 Pocket Gauge, Black
P549200 Cuff & Bladder, Adult, 2-Tube, Black, Latex
P549220 Cuff & Bladder, Adult, 2-Tube, Black, Latex-Free
P549240 Cuff & Bladder, Large Adult, 2-Tube, Black, Latex-Free

P549500 Inflation System, Adult, Black, Latex
P549505 Inflation System, Child, Black, Latex
P549520 Inflation System, Large Adult, Black, Latex
P549515 Inflation System, Adult, Black, Latex-Free
P549510 Inflation System, Child, Black, Latex-Free
P549525 Inflation System, Large Adult, Black, Latex-Free

Distributed By:

Pro
ADVANTAGE[®]
by **NDC**